

MAGNAMED

Operations Manual

Electronic Pulmonary Ventilator

FlexiMag Max 700

Fleximag Max 500

Fleximag Max 300

This operating manual includes the 3 models of pulmonary ventilators of FlexiMag Max's family: Fleximag Max 700, Fleximag Max 500 and Fleximag Max 300 developed and manufactured by Magnamed Tecnologia Medica S / A.

Review of this operation manual: 02

CE 2460

ANVISA Registration nº 80659160005

All rights reserved by:

Magnamed Tecnologia Médica S/A



Rua Santa Mônica, 801, 831

06715-865 – Parque Industrial San José -Cotia – SP Brazil

Tel/Fax: +55 11 4615-8500

E-mail: magnamed@magnamed.com.br

Website: www.magnamed.com.br

CNPJ: 01.298.443/0002-54

State registration: 149.579.528.111

Technical Responsible: Toru Miyagi Kinjo

Inscription CREA-SP: 5061555031

Legal Responsible: Wataru Ueda

Table of Contents

1 SAFETY NOTES	14
1.1 DEFINITIONS	14
1.2 WARNING.....	15
1.3 ATTENTION	17
1.4 OBSERVATION	17
2 FEATURES	18
2.1 INTENDED USE.....	18
2.2 PRINCIPLE OF OPERATION	18
2.3 GENERAL FEATURES	20
2.4 TECHNICAL CHARACTERISTICS	21
2.5 SAFETY FEATURES	22
3 UNPACKING THE PRODUCT	23
3.1 INITIAL VERIFICATIONS	23
3.2 LIST OF COMPONENTS.....	23
3.3 OPTIONAL PARTS AND ACCESSORIES	25
4 IDENTIFICATION OF COMPONENTS	27
4.1 FRONT VIEW	27
4.1.4 <i>Front panel</i>	28
4.2 BACK VIEW.....	29
4.3 VIEW FROM THE BOTTOM OF THE DISPLAY	30
5 PREPARATION FOR USE	31
5.1 ASSEMBLY	31
5.2 ELECTRICAL NETWORK CONNECTION.....	41
5.3 VERIFICATION BEFORE USE	42
6 INSTRUCTIONS FOR USE	44
6.1 INITIAL SEQUENCE	44

6.2 RESTRICTED MENU	47
6.2.1 <i>Restricted menu</i>	47
6.2.1.1 Patient	47
6.2.1.2 Autotest.....	47
6.2.1.3 Calibration	47
6.2.1.4 Status.....	48
6.2.1.5 Battery.....	48
6.2.1.6 System	48
6.2.1.7 Advanced.....	48
6.2.1.8 Ventilation.....	48
6.3 MAIN SCREEN	49
6.3.1 <i>Battery status area</i>	49
6.3.2 <i>Patient information area</i>	49
6.3.2.1 Volume x weight patient definition.....	49
6.3.3 <i>Historical events area</i>	50
6.3.4 <i>Alarm setting area</i>	50
6.3.5 <i>Mute alarm button</i>	50
6.3.6 <i>Manual cycle button</i>	50
6.3.7 <i>Inspiratory hold button</i>	50
6.3.8 <i>Expiratory hold button</i>	51
6.3.9 <i>O₂+ button (O₂ flush)</i>	51
6.3.10 <i>Nebulizer / TGI button</i>	51
6.3.11 <i>Graphic freeze button</i>	51
6.3.12 <i>Menu access area</i>	51
6.3.12.1 Layout.....	51
2 waveforms	51
Big numbers	51

Monitor	52
3 waveforms	52
Lung	52
2 loops	52
6.3.12.2 Resources	52
6.3.12.2.1 Noninvasive Ventilation (NIV)	53
6.3.12.2.2 External Auxiliary pressure	54
6.3.12.2.3 Tube compensation	54
6.3.12.3 Maneuver	54
6.3.12.3.1 P0.1	55
6.3.12.3.2 Slow Vital Capacity	55
6.3.12.3.3 P/V Flex	55
6.3.12.3.4 Pi Max	56
6.3.12.3.5 Gas Trapping Volume	56
6.3.12.4 Trend	56
6.3.12.5 Calibration	56
6.3.12.5.1 O ₂ cell (galvanic cell only)	57
6.3.12.6 Advanced	57
6.3.13 Lock Button	57
6.3.14 Stand by button	57
6.3.15 Bargraph or monitored parameters area	57
6.3.16 Graphics area	57
6.3.17 Ventilation Mode Setting Area	58
6.3.17.1 Ventilatory modes available	58
6.3.17.2 Adjusting ventilatory modes	59
6.4 CALIBRATIONS	60
6.4.1 Distal Flow Sensor	60
6.4.2 Integrated valve	60

6.4.3 <i>O₂ cell (galvanic cell only)</i>	60
6.5 AUTOMATIC COMPENSATION OF THE RESPIRATORY CIRCUIT	61
7 TROUBLESHOOTING	62
8 CLEANING, DISINFECTION AND STERILIZATION	64
8.1 EXTERNAL PARTS.....	64
8.2 COMPONENTS	64
8.2.1 <i>PROXIMAL flow sensor, silicone line and integrated valve</i>	64
8.2.1.1 Wash	64
8.2.1.2 Rinse	65
8.2.1.3 Drying	65
COMPOSITION OF THE ENZYMATIC DETERGENT.....	65
8.3 DISINFECTION.....	65
8.3.1 <i>External Parts</i>	65
8.3.2 <i>Exhalation valve, proximal flow sensor and silicone line</i>	66
8.3.3 <i>Distal flow sensor (Envitec SpiroQuant A+)</i>	66
8.3.4 <i>IRMA CO₂ sensor</i>	66
8.3.5 <i>Oximetry sensor (oximeter)</i>	67
8.4 STERILIZATION	68
8.4.1 <i>Steam Sterilization</i>	68
8.5 IMPORTANT ADVICES	68
8.6 PROCESSING METHODS.....	70
9 PREVENTIVE MAINTENANCE	71
9.1 VERIFICATIONS	71
9.2 PREVENTIVE MAINTENANCE SCHEDULE	72
9.3 CONSUMABLE ITEMS.....	72
9.4 INTERNAL BATTERIES.....	72
9.5 WATER COLLECTORS WITH COALESCING FILTER.....	73

9.6 MESH FILTER	74
9.7 O ₂ CONCENTRATION CELL.....	75
10 TECHNICAL SPECIFICATIONS.....	77
10.1 CLASSIFICATION	77
10.1.1 <i>Class I Equipment</i>	77
10.1.2 <i>Protection class of applied parts</i>	77
10.1.3 <i>CE – Class IIB</i>	77
10.1.4 <i>ANVISA – Class III</i>	77
10.1.5 <i>FDA – Class II</i>	77
10.2 APPLICABLE STANDARDS.....	77
10.3 PHYSICAL AND ENVIRONMENTAL SPECIFICATIONS	78
10.4 ELECTRICAL SPECIFICATIONS	79
10.4.1 <i>Electrical network</i>	79
10.4.2 <i>Batteries</i>	80
10.4.3 <i>External DC Power Supply</i>	80
10.4.4 <i>Connectors</i>	81
10.4.4.1 <i>Protocol used for data communication with external devices</i>	82
10.4.5 <i>Gas Inlet Connector</i>	83
10.5 PNEUMATIC SPECIFICATIONS	84
10.5.1 <i>Pneumatic layout</i>	84
10.6 INTERNAL FLOW TRANSDUCER SPECIFICATIONS	85
10.7 DISTAL FLOW SENSOR SPECIFICATIONS.....	86
10.8 PROXIMAL FLOW SENSOR SPECIFICATIONS	86
10.9 PRESSURE SENSOR SPECIFICATIONS.....	86
10.10 GALVANIC O ₂ CELL SPECIFICATIONS.....	87
10.11 PARAMAGNETIC O ₂ CELL SPECIFICATION	88
10.12 VENTILATION MODES SPECIFICATIONS	89

10.12.1 VCV.....	89
10.12.2 PCV.....	91
10.12.3 PLV	93
10.12.4 PRVC.....	95
10.12.5 V-SIMV	97
10.12.6 P-SIMV.....	100
10.12.7 CPAP/PS	103
10.12.8 DualPAP	106
10.12.9 APRV.....	109
10.12.10 MMV	111
10.12.11 VS	112
10.12.12 Nasal CPAP.....	113
10.12.13 NIV	113
10.12.14 VG.....	114
10.12.15 O ₂ THERAPY.....	115
10.13 OXYGEN ADJUSTMENT RESPONSE TIME	116
10.14 ACCURACY OF SETTINGS.....	116
10.15 ADJUSTABLE PARAMETER SPECIFICATIONS	117
10.16 MONITORED PARAMETER SPECIFICATIONS	120
10.17 SAFETY SYSTEM SPECIFICATIONS	122
10.18 ALARM SYSTEM SPECIFICATION	123
10.18.1 Adjustable alarm specifications	124
10.18.2 Ventilator alarm messages	125
10.18.3 Ventilator Alert Messages.....	128
10.18.4 IRMA CO ₂ sensor alarm messages.....	129
10.18.5 Oximeter Alarm Messages	131
10.18.6 Alarm Testing.....	132

10.18.6.1 FiO ₂ alarm	132
10.18.6.2 Airway pressure alarm	133
10.18.6.3 Volume alarm	133
10.18.6.4 No AC power alarm	133
10.18.6.5 Low battery alarm	133
10.18.6.6 Gas supply failure alarm	133
10.18.6.7 PEEP alarm	133
10.18.6.8 Disconnection alarm.....	133
10.18.6.9 Obstruction alarm	134
10.18.6.10 Minute volume alarm	134
10.18.6.11 Respiratory rate alarm	134
10.18.6.12 Apnea alarm	134
10.18.6.13 EtCO ₂ alarm	134
10.18.6.14 ICO ₂ alarm	135
10.18.6.15 Heart rate alarm	135
10.18.6.16 SpO ₂ alarm.....	135
<i>10.18.7 Battery test</i>	<i>135</i>
10.19 PERFORMANCE SPECIFICATIONS	135
10.20 RESPIRATORY CIRCUIT SPECIFICATIONS	136
10.21 MAINTENANCE AND CALIBRATION SPECIFICATIONS	137
10.22 IRMA CO ₂ SENSOR SPECIFICATIONS	137
10.23 OXIMETER SPECIFICATIONS.....	140
10.24 ELECTROMAGNETIC COMPATIBILITY	141
<i>10.24.1 Guidelines and manufacturer's declaration - Electromagnetic emissions</i>	<i>142</i>
<i>10.24.2 Guidelines and manufacturer's declaration - Electromagnetic immunity.....</i>	<i>142</i>
<i>10.24.3 Electrical safety</i>	<i>146</i>
10.25 AUDIBLE ACOUSTIC ENERGY	146

10.26 MASK TO NON-INVASIVE VENTILATION	146
10.27 RESPIRATORY CIRCUIT	147
10.28 HME FILTER	147
10.29 HEPA FILTER	148
11 IRMA CO₂ SENSOR (OPTIONAL)	149
11.1 INTENDED USE.....	149
11.2 INSTRUCTIONS FOR USE	149
11.2.1 Reset procedure	150
11.2.2 LED status information.....	151
11.3 PREVENTIVE MAINTENANCE.....	151
11.4 IMPORTANT NOTICES.....	151
12 PULSE OXIMETER (OPTIONAL).....	153
12.1 INTENDED USE.....	153
12.2 PRINCIPLE OF OPERATION	153
12.3 IMPORTANT NOTICES.....	155
13 TECHNICAL SERVICE	157
14 DISCARD.....	158
15 SYMBOLOGY	159
15.1 SYMBOLS USED IN THE EQUIPMENT	159
15.2 SYMBOLS USED IN PACKAGING AND LABELING	160
16 ABBREVIATIONS AND TERMS USED.....	162
17 BIOCOMPATIBILITY STATEMENT	164
18 WARRANTY	165
19 TRAINING	166
ANNEX I – DIFFERENCES BETWEEN MODELS	167

Table index

Table 1 - Initial Verifications	23
Table 2 - Components that accompany FlexiMag Max	23
Table 3 - OPTIONAL parts and accessories.....	25
Table 4 - Verification before use	42
Table 5 - Ventilatory modes	58
Table 6 – Sterilizable accessories	68
Table 7 – Processing methods	70
Table 8 – Preventive maintenance schedule	72
Table 9 – Consumable items replacement	72
Table 10 - Physical and environmental specifications	78
Table 11 – Electrical network.....	79
Table 12 - Batteries.....	80
Table 13 – External DC power supply	80
Table 14 - Connectors.....	81
Table 15 - Gas Inlet Connector.....	83
Table 16 - Internal Flow Transducer - General Specifications.....	85
Table 17 - Distal Flow Sensor - General Specifications	86
Table 18 - Proximal Flow Sensor - General Specifications ADULT	86
Table 19 - Pressure sensor - General specifications.....	86
Table 20 – Galvanic O ₂ cell specifications.....	87
Table 21 – Paramagnetic O ₂ cell specification.....	88
Table 22 – Oxygen adjustment response time	116
Table 23 – Accuracy of settings.....	116
Table 24 – Uncertainty of settings	117
Table 25 - Adjustable parameters.....	117
Table 26 - Monitored ventilatory parameters	120

Table 27 - Alarm Priority	123
Table 28 - Alarms Characteristics.....	123
Table 29 - Adjustable alarms	124
Table 30 - High priority alarms.....	125
Table 31 - Medium Priority Alarms.....	126
Table 32 - Low priority alarms.....	127
Table 33 - Alert messages	128
Table 34 - High priority alarms.....	129
Table 35 - Medium Priority Alarms.....	130
Table 36 - Low priority alarms.....	130
Table 37 - High priority alarms.....	131
Table 38 - Medium Priority Alarms.....	131
Table 39 - Low priority alarms.....	132
Table 40 - Performance specifications.....	135
Table 41 - Expiratory and inspiratory branch resistance specifications.....	136
Table 42 - Respiratory circuit compliance specification.....	136
Table 43 - Operation Specifications.....	137
Table 44 - Maintenance and calibration specifications	137
Table 45 - IRMA CO ₂ sensor - General specifications.....	137
Table 46 - IRMA CO ₂ Sensor - Outputs.....	138
Table 47 - IRMA CO ₂ sensor - Gas analyzer.....	138
Table 48 - IRMA CO ₂ sensor – Accuracy I	139
Table 49 - IRMA CO ₂ sensor - Accuracy II	139
Table 50 - IRMA CO ₂ sensor - Effects of gas and vapor interference	139
Table 51 – Quantitative effects of humidity and condensation	140
Table 52 – Performance specifications.....	140
Table 53 – Environmental specifications	140

Table 54 – Audible acoustic energy specification	146
Table 55 – Mask to NIV ventilation specification	146
Table 56 – Respiratory circuit specification	147
Table 57 – HME filter specification	147
Table 58 – HEPA filter specification.....	148
Table 59 - IRMA CO ₂ Status LEDs	151
Table 60 - Symbols used in the equipment.....	159
Table 61 - Symbols used in packaging and labeling	160
Table 62 - Abbreviations and terms used	162

1 Safety Notes

1.1 Definitions



WARNING

- *Informs the user of the possibility of injury, death or other serious adverse reaction associated with the use or misuse of the equipment.*
-



ATTENTION

- *Informs the user of the possibility of equipment failure associated with the use or misuse, such as: malfunction, damage to the equipment itself, damage to third party assets and indirect injury to the patient.*
-




Observations

- *Important information to observe.*
-

1.2 Warning



Warning

- Whenever there is a  symbol, consult this manual for more information.
- This manual shall be read in its entirety carefully to ensure correct and safe use of the equipment and to provide maximum safety and the best resources to the patient. Observe all "WARNING" and "ATTENTION" contained in this manual and on the labeling of the equipment.
- The equipment shall be used only for the purpose specified in Intended Use (Chapter 2.1) and with appropriate monitoring.
- The equipment shall be operated by trained and qualified personnel, who shall maintain vigilance during use, including in volume limited ventilation, so that in case of malfunction or failure, help is available immediately.
- **Explosion Hazard:** This equipment is not approved for use with flammable anesthetic agents.
- Equipment may be adversely affected and may be subject to interference from certain transmitting equipment, such as cellular telephones, walkie-talkie, cordless telephones, pagers, high frequency surgical equipment (diathermy), and shortwave therapy, which may interrupt the operation of the equipment. Do not use this transmission equipment near the ventilator.
- This equipment shall not be used during a nuclear magnetic resonance (MTR, NMR, NMI), as it may cause interference, causing adverse effects to the patient.
- Before first use or after use on each patient, clean the equipment as indicated in chapter 8.
- When powering up the equipment, perform basic scanning and adjustment procedures. Follow the instructions in chapter 5.
- ALARMS and WARNINGS shall be promptly serviced in order to maintain the integrity of the equipment operation and patient safety.
- Do not use electrically conductive antistatic hoses or tubes in the respiratory system of the ventilator.
- Always check that the equipment is correctly set before using it.
- After starting the ventilation, using monitoring verify if ventilatory parameters are adequate.
- Use only MAGNAMED accessories listed in this manual, which have been tested and approved for use in conjunction with this equipment. Otherwise, functionality of the equipment may be compromised.
- During prolonged use of the equipment with patients with excess secretion and / or using a heated humidifier, the condition of the flow sensor shall be checked frequently,

attempting to drain accumulated fluids in the respiratory circuit when necessary.

- *Keep the equipment connected to an electrical source even when it is turned off, so that the internal batteries are always charged.*
- *After using the equipment or after an extended period of storage, fully recharge the batteries.*
- *If LOW BATTERY alarm occurs, IMMEDIATELY connect the ventilator to an electrical Network. If this is not possible, provide another suitable means of ventilatory support and then DISCONNECT patient equipment.*
- *Recharge the batteries before the next use, otherwise, any electrical power failure may interrupt the operation of the ventilator.*
- *Always check that there are no obstructions, as it is extremely important for a correct monitoring of ventilatory parameters that the respiratory circuit is completely unobstructed.*
- *Never obstruct the pressure connections. The pressures measured at these points are used by the patient monitoring system.*
- *After use, reusable respiratory circuit components shall go through a high-level disinfection process before their next use.*
- *All parts of the equipment that come in contact with fluids from the patient SHALL be discarded, disinfected, sterilized or disposed of as potentially infectious hospital waste.*
- *All parts of the FlexiMag Max that are in the gas passage way, including accessories and applied parts, are made of non-toxic, latex-free and phthalate-free material and do not cause irritation or allergy to the patient.*
- *Common use non-exclusive FlexiMag Max accessories, such as masks, circuits, nebulizers, humidifiers, HME filters, among others, shall be registered with local government.*
- *Do not use the equipment if a problem cannot be solved.*
- *Have a manual ventilation device available for the case of complete discharge of the battery, lack of gases for the operation of the device or even for general failure of the ventilator. The absence of manual ventilation in these cases may result in death of the patient.*
- *Always use officially approved oxygen cylinders and pressure reducing valves that meet local government requirements.*
- *For proper ventilation, when adjusting the ventilator, take into consideration the dead spaces of the respiratory circuit, especially for low tidal volumes.*
- *The ventilator shall not be covered or positioned in such a way that operation or performance is affected.*
- *When components of the respiratory circuit or other components or subassemblies are added to the respiratory system of the ventilator, the pressure gradient across the respiratory system, measured relative to the patient's port of attachment, may increase, adversely affecting the performance of the ventilator.*

-
- *Nebulization or humidification can increase the resistance of respiratory system filters. The operator shall often monitor increased resistance and filter blockage of the respiratory system.*
 - *The ventilator cannot be used in a hyperbaric chamber.*
 - *Only use air and oxygen in the gas inlet of the equipment. Do not use nitric oxide, helium or other gases.*
 - *The expiratory branch may become contaminated with bodily fluids or expired gases during use of the equipment under normal conditions and under a single fault.*
 - *HME filter, HEPA filter and airway adapter are single use. The reuse of these accessories may cause cross contamination.*
 - *The distal flow sensor is for single patient use and shall be discarded when changing patients. Re-use for another patient may result in cross-contamination.*
 - *FlexiMag Max does not generate sub atmospheric pressure during the expiratory phase.*
-

1.3 Attention



ATTENTION

- *This ICU ventilator does not emit electromagnetic waves that may interfere with the operation of equipment in its vicinity.*
 - *Perform all maintenance according to the operating time specified in this manual.*
 - *All service, modification or maintenance on the ventilator can only be carried out by qualified technician, trained and duly authorized by MAGNAMED.*
 - *Only use parts, cables, sensors, filters and breathing circuits specified by MAGNAMED and compatible with the system. For purchase, please inform the codes presented in the corresponding chapter.*
-

1.4 Observation



Observations

- *The technical specifications of MAGNAMED products are subject to change without notice.*
-

2 Features

2.1 Intended Use

FlexiMag Max's ventilator family has been developed to provide invasive and non-invasive ventilatory support with complete ventilatory monitoring for neonatal, pediatric and adult patients with impaired respiratory function in postoperative intensive care units in post-anesthetic recovery rooms (RPA), or for in-hospital transport.

2.2 Principle of Operation

FlexiMag Max is a family of a micro processed electronic lung ventilators whose principle of operation is based on the integration between specific modules:

- Pneumatic module (manifold)
- Electronic Control Module
- Electronic interface module

At the entrance of the pneumatic module, two valves regulate the pressure of the gases coming from the hospital network or cylinders, trying to assure a suitable range to the equipment.

Together with these regulating valves, there are microswitches or limit switches that constantly monitor the gas pressure, so that the insufficiency or absence of pressure of one or both gases is immediately indicated by a priority alarm.

Subsequently, electronically controlled air¹ and O₂ proportional valves regulate the flow of gases accurately to ensure concentration and volume appropriate to each situation.

For the models Fleximag Max 300 and Fleximag Max 700, the ventilator relies on the FlowAir system², consisting of a turbine to electronically control the flow of air.

After their respective flows are adjusted, the gases are mixed to make the measurement of the O₂ concentration and the measurement of the resulting flow.

The O₂ concentration is obtained through a galvanic cell or optionally, a paramagnetic cell by indirect contact with the gas of the patient, through the passage of gas in the sensor.

¹ Only for Fleximag Max 500 and Fleximag Max 700 models

² Only for Fleximag Max 300 and Fleximag Max 700 models

The resulting flow measurement is performed by a high precision internal flow sensor, which guarantees an adequate reading, without the need for recalibrations and has mass flow technology, whose reading is independent of temperature or pressure.

The flow expired by the patient is measured by an external sensor, of the thermal type or hot wire anemometer, connected to the expiratory valve. Optionally, this flow can be read through a proximal flow sensor, connected to the patient's "Y" output, whose measurement is based on the pressure differential between two points.

The system pressures are taken through points in the pneumatic module, which are connected to the transducers in the electronic control module. The maximum ventilator pressure is guaranteed by the control system, which monitors the pressure in the system through transducers by setting the PR Limit or PR Control. The pressure can be relieved by expiratory valve opening and overpressure.

All these flow and pressure measurements are converted into digital signals by the electronic control module and serve to feed the control algorithm uninterruptedly, ensuring a gradual and safe adjustment of the ventilatory process.

The pneumatic module also includes safety valves, such as the overpressure valve and the anti-asphyxiation valve.

The input and output of information are processed by the electronic interface module. The information entered by the operator, via touch screen or via button, is translated, interpreted and sent to the electronic control module by serial communication, through secure protocols. Through this information the ventilator establishes the appropriate parameters to act in each different situation.

As soon as it receives information, the control module also sends it to the interface module. All measured or calculated data are also sent, via serial, to the interface module. It is up to this module to treat and display this information to the operator in a friendly and intuitive way.

All risk situations that require operator intervention are analyzed by the control module and sent to the interface module that then issues, according to the degree of risk, the necessary alarms or alerts.

Fleximag Max forms an electromedical system with the humidification system.

Observations

- *Parts and pieces indicated in this manual may be used in the patient's environment.*

2.3 General features

- Integration of the entire system into a single low volume and weight electronic pneumatic module;
- Electronic flow control system for the patient, using digital technology;
- Intelligent and integrated alarm system with all safety systems required by national and international standards;
- Power backup system;
- Possibility to use the parameter adjustment performed before the equipment was switched off;
- Volatile storage (while the equipment is not turned off), of all parameters, for each ventilatory mode, allowing a quick return to previously used settings (e.g. weaning);
- Fully adjustable equipment via both the touch screen and the turn-confirm button;
- Adjustable resistance compensation for both endotracheal and tracheostomy tubes;
- Volume compensation / compliance according to the respiratory circuit;
- Automatic altitude compensation;
- Automatic barometric pressure compensation.
- Leak Compensation;
- Trend charts, with memorization of the last 240 hours of ventilation;
- Resources for evaluating the patient's respiratory mechanics (P0.1, Vital Capacity, PV Inflection Points, Maximum Pi, Trapped Volume);
- Inspiratory and expiratory hold with variable time, determined by the operator;
- Freezing and saving up to 3 simultaneous loops (PxV and VxF);
- Complete monitoring;
- Additional monitor with data related to the respiratory mechanics of the patient;
- Breath adjustment option in controlled volume ventilatory modes;
- Alarm mute time adjustment option;
- Time adjustment and O₂ concentration for aspiration Option;
- Monitor with complete ventilation graphics:
 - Pressure x Time Curve
 - Flow x Time Curve
 - Volume x Time Curve
 - Pressure x Volume Loop
 - Volume x Stream Loop
 - Pressure x Flow loop
- Instant pressure bar graph with numerical indicator of peak, plateau or instantaneous pressure;

 **ATTENTION**

- *This equipment shall be operated only by qualified personnel and properly trained for its use.*
-

2.4 Technical Characteristics

- 15 "color liquid crystal display (LCD) with 1024 x 768-pixel resolution, touch sensitive;
- Screen configurable for night mode;
- Control Boards with:
 - Presentation of the data in the display;
 - Processor IMX6 SOM-IMX6;
- Hotkeys for:
 - STAND BY
 - SILENCE OF ALARMS (2 min)
 - O₂ + (oxygen flush)
 - MANUAL (manual trigger)
 - INSP HOLD (inspiratory hold)
 - EXP HOLD (expiratory hold)
 - FREEZE (freeze graphics)
 - LOCK (keyboard lock)
 - NEB / TGI (Nebulizer or TGI)
 - MODE (ventilatory modes)
 - MENU (displays more options)
 - EVENTS (displays active and historical alarms)
 - PATIENT (adjusts patient data)
 - BATTERY (Displays advanced battery information)
- Pressure reading in respiratory circuit;
- Regulated pressure reading;
- Network pressure reading;
- Auxiliary pressure reading;
- Barometric pressure reading;
- Reading of the O₂ concentration in the gas mixture administered;
- Speaker for alarms and alerts;
- High brightness LED for prompt identification of alarms with 360-degree visibility;
- GREEN LED indicator for power supply connection, being switched on when the equipment is connected to a power supply and off when the equipment is being powered by the battery only;

-
- Single distal flow sensor for all patient types;
 - Proximal flow sensor for each type of patient;
 - Internal flow sensor for all patient types;
 - External input 100-240 VAC - 50 - 60 Hz;
 - On-off key;
 - O₂ cell galvanic or optionally paramagnetic (non-consumable);
 - Nebulizer or TGI;
 - Volume, pressure and concentration compensation when used with nebulizer;
 - Stand with locking 5" wheels, heated humidifier support;
 - Oximetry sensor with dedicated input;
 - Capnography sensor with dedicated input;
 - Nursing Call Out;
 - FlowAir system consisting of high-performance high-flow turbine and anti-noise system¹;
 - High flow therapy (flowmeter mode);
 - Software Updates via USB;
 - Engineering mode;
 - Maintenance mode for technical assistance;
 - Low O₂ pressure inlet;
 - HL7 Protocol;
 - Standard Ethernet Connector that enables HL7 protocol communication;
 - RS-232C serial interface;
 - Adult, pediatric and neonatal respiratory circuit;

2.5 Safety features

- Anti-asphyxia valve for protection against gas supply failure.
- 100hPa relief valve, in compliance with the basic standard of ventilators, avoiding possible overpressure in the respiratory circuit.

¹ Only for Fleximag Max 300 and Fleximag Max 700

3 Unpacking the Product

3.1 Initial Verifications

Observations

- *If the package is damaged, DO NOT OPEN and report immediately to the carrier responsible and to Magnamed.*
-

Table 1 - Initial Verifications

Step	Procedure
1	Check that the packaging is intact, making sure there are no dents, punctures or any other form of damage.
2	Open the package carefully, following the indications on the box.
3	Check the contents of the package.

3.2 List of Components

The following items are an integral part of the equipment and are for the exclusive use of the equipment:

Table 2 - Components that accompany FlexiMag Max

Item	Code	Description	Qty.	UMI
1	1106630	FLEXIMAG MAX 700 - LUNG VENTILATOR NEONATE TO ADULT	01	PC
	1107270	FLEXIMAG MAX 500 - LUNG VENTILATOR NEONATE TO ADULT	01	PC

Item	Code	Description	Qty.	UMI
	1107240	FLEXIMAG MAX 300 - LUNG VENTILATOR NEONATE TO ADULT	01	PC
2	1703938	KIT 5 SENSORS SPIROQUANT ENVITEC CE 0123	01	PC
3	1705043	KIT FLOW SENSORS ADU INF NEO AUTOCLAVABLE 1,6M	1	PC
4	1707302	ARTICULATED ARM TO SUPPORT RESPIRATORY CIRCUITS	1	PC
5 ⁽¹⁾	1707451 ⁽¹⁾	RESPIRATORY CIRCUIT ADULT WITH WATER TRAP STRAIGHT Y ⁽¹⁾	1	PC
6	3206170	DIAPHRAGM FOR INTEGRATED VALVE	1	PC
7	3806167	INTEGRATED VALVE SET	1	PC
8	3002739	M5X12 SCREW PHS INOX	4	PC
9	9003365	3/16 PHILLIPS KEY	1	PC
10	2804669	3 WAY MOUNTED AC NETWORK CABLE 3M - NBR 14136	1	PC

¹ Not applicable for the European Community

Item	Code	Description	Qty.	UMI
11	1106630-NE-22-RR	OPERATION MANUAL	1	PC
12	1706378	MESH FILTER FOR FLEXIMAG MAX	1	PC
13	7007022	QUICK GUIDE - FLEXIMAG MAX	1	PC
14 ⁽¹⁾⁽²⁾	3903114 ^{(1) (2)}	COMPRESSED AIR DISS X2 3M HOSE ⁽¹⁾⁽²⁾	1	PC
15 ⁽²⁾	3902647 ⁽²⁾	O2 DISS X2 HOSE ⁽²⁾	1	PC

3.3 Optional Parts and Accessories



ATTENTION

- *Always use original parts and accessories to ensure the safety and efficacy of the equipment.*

Table 3 - OPTIONAL parts and accessories

Item	Code	Description	Qty.	UMI
1 ⁽²⁾	1707453 ⁽²⁾	RESPIRATORY CIRCUIT NEONATAL WITH WATER TRAP Y 90 ⁽²⁾	1	PC
2 ⁽²⁾	1707452 ⁽²⁾	RESPIRATORY CIRCUIT PEDIATRIC WITH WATER TRAP Y 90 ⁽²⁾	1	PC
3 ⁽²⁾	1704601 ⁽²⁾	RESPIRATORY CIRCUIT - ADULT 1.6M STRAIGHT Y AUTOCLAVABLE ⁽²⁾	1	PC

¹ Only for Fleximag Max 500 and Fleximag Max 700 models


² Not applicable for the European Community

Item	Code	Description	Qty.	UMI
4	1704603	RESPIRATORY CIRCUIT - PEDIATRIC 1.6M Y 90 - AUTOCLAVABLE ⁽¹⁾	1	PC
1	1704396	CO ₂ MAINSTREAM SENSOR WITH INTERCONNECTION CABLE AND 5 ADULT AIRWAY ADAPTER CE 0413	1	PC
2	1704409	PULSE OXIMETER SPO ₂ MASIMO - ADAPTER CABLE AND PED-ADU SENSOR CE 0123	1	PC
3	1704410	PULSE OXIMETER SPO ₂ MASIMO - ADAPTER CABLE AND NEONATAL SENSOR CE 0123	1	PC
4	3201100	ADULT FLOW SENSOR - AUTOCLAVABLE	1	PC
5	3201099	INFANT FLOW SENSOR - AUTOCLAVABLE	1	PC
6	3201098	NEONATAL FLOW SENSOR - AUTOCLAVABLE	1	PC
7	3802058	SILICONE LINE 1,6 M WITH UNIVERSAL CONNECTOR	1	PC
8	1704395	CO ₂ SENSOR ADULT AIRWAY ADAPTER CE 0413	1	PC
9	1704394	CO ₂ SENSOR NEONATAL AIRWAY ADAPTER CE 0413	1	PC

¹ Not applicable for the European Community


4 Identification of Components

4.1 Front view



- 4.1.1 Alarm indicator light
- 4.1.2 Color and touch screen LCD monitor
- 4.1.3 Turns and confirms button with led indicator of power supply
- 4.1.4 Front panel
- 4.1.5 Carrying handle
- 4.1.6 Stand
- 4.1.7 Wheels with brakes
- 4.1.8 Articulated arm

4.1.4 Front panel



4.1.4.1 Integrated valve with inspiratory and expiratory branch connector

4.1.4.2 Nozzles for connecting proximal flow sensor lines


4.1.4.3 Distal sensor

4.1.4.4 Connector for external sensors (capnograph and oximeter)

4.1.4.5 External auxiliary pressure reading nozzle

4.1.4.6 Connector for nebulizer or TGI

4.2 Back view¹



4.2.1 Low pressure O₂ inlet^{4,3}

4.2.2 Nursing call

4.2.3 Identification label

4.2.4 Air outlet for cooling

4.2.5 Compressed air inlet²

4.2.6 Push on/off power button

4.2.7 Input for external power supply

4.2.8 Input for electrical energy

4.2.9 Fuse Holders

4.2.10 Water collectors with coalescing filter for gas under high pressure


4.2.11 Mesh filter

4.2.12 Oxygen gas input

¹ Back view of Fleximag Max 700 model

² Only for Fleximag Max 500 and Fleximag Max 700 models

4.3 View from the bottom of the display



4.3.1 Standard HDMI connector


4.3.2 Standard USB Connectors



4.3.3 Standard RJ-45 Ethernet Connector



4.3.4. Standard RS-232 connector




5 Preparation for Use


5.1 Assembly


OK	Item	Assembly Sequence	Image
<input type="checkbox"/>	1	Use the bolt and wrench that came with the equipment to screw the wheel base to the ventilator module.	---
<input type="checkbox"/>	2	<p>Position the diaphragm in the expiratory branch on the integrated valve in the position shown in the figure.</p> <p>Connect the integrated valve on the front panel.</p> <p>Push the valve lock lever.</p>	



OK	Item	Assembly Sequence	Image
		<p>Check that the lever is fixed.</p>	
<input type="checkbox"/>	<p>3</p>	<p>Attach the articulated arm to the transport handle and rotate the handle to secure.</p>	



OK	Item	Assembly Sequence	Image
<input type="checkbox"/>	4	<p>Connect the patient circuit in the inspiratory and expiratory branches of the integrated valve and position on the articulated arm according to the image.</p>	
<input type="checkbox"/>	5	<p>Connect the flow sensor.</p> <p>A. Envitec Distal Sensor: Connect the sensor to the integrated valve.</p>	




OK	Item	Assembly Sequence	Image
		<p>B. Proximal sensor: Connect the proximal sensor line as shown in the figure to the right. Connect the flow sensor to the patient's respiratory circuit after the Y. Connect the other end of the proximal sensor line to the flow sensor in the position shown in the figure on the right.</p> <div data-bbox="363 1055 858 1115" style="background-color: black; color: white; padding: 5px;">  ATTENTION </div> <ul style="list-style-type: none"> <li data-bbox="411 1151 839 1462">• <i>The flow sensor connectors shall face up to avoid condensation and accumulation of secretion at the measurement points.</i> 	 



OK	Item	Assembly Sequence	Image
<input type="checkbox"/>	6	<p>When using a humidifier, connect the inspiratory branch to the humidifier as shown in the figure.</p> <p>If a heated humidifier is used, it shall comply with ISO 80601-2-74:2017</p>	

OK	Item	Assembly Sequence	Image
<input type="checkbox"/>	7	<p>When using the IRMA CO₂ sensor, connect the airway adapter to the IRMA CO₂ sensor.</p> <p>Connect the IRMA CO₂ sensor just after the proximal flow sensor, or directly to the Y connector.</p> <p>Connect the cable to one of the external sensor connectors on the front panel.</p>	

OK	Item	Assembly Sequence	Image
<input type="checkbox"/>	8	<p>If the respiratory circuit is to be used with the IRMA CO₂ sensor and the HME (Heat and Moisture Exchange) or HEPA filter (protection against viruses and bacteria), mount it according to the image.</p> <p>ATTENTION</p> <ul style="list-style-type: none"> • Only use filters specified by MAGNAMED. 	
<input type="checkbox"/>	9	<p>If respiratory circuit is to be used for NON-INVASIVE VENTILATION (NIV or NIV - Noninvasive Ventilation) with mask, in addition to IRMA CO₂ sensor and filter, follow the sequence according to the image.</p> <p>ATTENTION</p> <ul style="list-style-type: none"> • Use only MASKS specified by MAGNAMED. • Use MASK suitable for patient type. 	

OK	Item	Assembly Sequence	Image
<input type="checkbox"/>	10	<p>If respiratory circuit is to be used for NON-INVASIVE VENTILATION (NIV or Non-Invasive Ventilation) using the mask and without the filter, mount it according to the image.</p>	
<input type="checkbox"/>	11	<p>If the respiratory circuit is to be used for NON-INVASIVE VENTILATION (NIV or NIV - Non-Invasive Ventilation) without the CO₂ sensor. Mount it according to the image.</p>	

OK	Item	Assembly Sequence	Image
<input type="checkbox"/>	12	<p>If respiratory circuit is to be used for NON-INVASIVE VENTILATION (NIV) with mask and HME filter, mount it according to the image.</p>	
<input type="checkbox"/>	13	<p>If using a HEPA filter, connect the filter between the expiratory branch of the integrated valve and the expiratory branch of the patient's circuit.</p>	
<input type="checkbox"/>	14	<p>If you want to use the oximetry sensor, connect the sensor cable to the external sensor connector on the front panel.</p> <p>Position the oximetry sensor on the patient's finger.</p>	

OK	Item	Assembly Sequence	Image
<input type="checkbox"/>	15	<p>Connect the power cord to the machine and to the power supply.</p>	
<input type="checkbox"/>	16	<p>Connect the compressed air hoses and O₂ as indicated, according the ventilator model.</p> <p>ATTENTION</p> <ul style="list-style-type: none"> • <i>Input pressures exceeding the specified limit can damage the equipment.</i> • <i>The hoses shall comply with ISO 5359:2008.</i> • <i>The gas network connected to the equipment shall complies with ABNT NBR 12188:</i> 	

OK	Item	Assembly Sequence	Image
		<i>2012.</i>	

5.2 Electrical Network Connection

The equipment shall be connected to a grounded electrical supply network complying with ABNT NBR 13534: 2008 - Instalações elétricas de baixa tensão - Requisitos específicos para instalação em estabelecimentos assistenciais de saúde.

The internal batteries of the equipment shall always be charged and ready for use in any power failure or for use in external operations. To do this, keep your power source connected to the electrical network to charge the batteries, even if the equipment is off. Switching between the power supply and internal battery does not affect the performance of the equipment.

After prolonged use of the equipment using only the internal battery power, it is necessary to fully recharge the equipment, preparing the equipment for the next use. During the recharge of the internal batteries, the performance of the equipment is not affected.

If the equipment remains unplugged from the electrical network for a period of more than one month, the batteries must be fully recharged.

When switching from electrical network to battery, the performance of the equipment is not affected. The "no AC power" low priority alarm is triggered, and the battery icon is displayed. If the equipment remains out of the electrical network, when the battery charge is less than 50%, the "low battery" medium priority alarm will be displayed and at least 5 minutes before the battery runs out the high priority alarm "low battery".



WARNING

- *If LOW BATTERY alarm occurs, IMMEDIATELY connect the ventilator to the electrical network. If this is not possible, provide another suitable means of ventilatory support and DISCONNECT the equipment from the patient.*








Observation

- *After the electrical network has been interrupted and restored while the internal battery is charged with*

the equipment in operation, the performance of the equipment will not be affected, and the accuracies will be maintained.

- *After a period of long interruptions to the electrical network, connect the equipment to the electrical network, switch on the equipment and wait for 30 min. Perform calibrations and autotest.*

When switching from battery to electrical network, the performance of the equipment is not affected. Alarms related to the battery are ceased and the battery icon is displayed.



Icon	Description	Alarm
	Equipment charged and connected to electrical network	No alarm
	Equipment operating on battery, with load above 50%	Low priority "No AC power"
	Equipment operating on battery, with load between 50% and 20%	Medium priority "Low battery"
	Equipment operating on battery with load below 20%	High priority "Low battery"
	Equipment connected to electrical network and battery charging	No alarm

5.3 Verification Before Use

The purpose of this verification is to ensure greater safety through some simple and quick procedures that must be performed before each use or at least at the beginning of each work period.

Table 4 - Verification before use

Item	Procedure
1	Verify that the machine is turned off.
2	Perform a visual inspection of the equipment and its components, seeking to assess the integrity of the equipment.
3	Verify that all components are properly connected or inserted.
4	Verify that the integrated valve is installed correctly and is firmly seated. Also note the conditions of the diaphragm.
5	Verify the tight connection of the distal flow sensor to the expiratory valve.

Item	Procedure
6	Verify that the breathing circuit is firmly connected and is suitable for the patient.
7	Verify the tight connection of oxygen and compressed air hoses.
8	<p>Verify that the inlet pressure is within the specified range.</p> <p> ATTENTION</p> <ul style="list-style-type: none"> • <i>Input pressures exceeding the specified limit can damage the equipment.</i> • <i>For inlet pressures less than 250 kPa, the maximum flow shall be 120 L/min.</i>
9	<p>Check the power cord connection, if applicable.</p> <p>The ventilator can be used on battery for up to 210 minutes continuously, under normal ventilation conditions.</p> <p> WARNING</p> <ul style="list-style-type: none"> • <i>If LOW BATTERY alarm occurs, IMMEDIATELY connect the ventilator to the electrical network. If this is not possible, provide another suitable means of ventilatory support and then DISCONNECT the equipment from the patient.</i>
10	If all items are OK, then the equipment is ready for use.

 WARNING

- *Perform all verification procedures before each use.*
- *If any problem is found, correct it BEFORE USING THE DEVICE.*
- *If the problem cannot be resolved immediately, contact qualified technical assistance.*

6 Instructions for use

6.1 Initial Sequence

Turn on the ventilator via the push on/off button on the back of the equipment.

You will be presented with the initial screen that includes the patient options and the available services, depending on the model of ventilator purchased.

From the home screen, select the type of patient by tapping on the corresponding icon.

Once the patient type is identified, the flow sensor and humidifier adjustment screen will be displayed, and the operator will have the option to perform the autotest or initiate ventilation.

Select the type of humidifier or heat exchanger and the chosen flow sensor (distal or proximal).



ATTENTION

- *The choice of humidifier or heat exchanger is important for proper calculation of volumes according to the temperature and humidity conditions presented (STPD or BTPS).*
- *If you have opted for the proximal sensor, make sure it matches the type of patient selected (neonatal, pediatric, or adult).*
- *Only one type of external, distal OR proximal flow sensor will be tested. If you made a mistake while selecting it, restart the equipment and redo the autotest.*
- *In case of the use of the distal flow sensor, do not use suction system in the exhaust gas outlet.*

To perform the autotest, be sure to perform the following steps:

- The ventilator shall be disconnected from the patient.
- The ventilator should preferably be connected to the electrical network or, if this is not possible, shall have enough battery power.
- The ventilator shall be supplied with both gases (compressed air and O₂) within the recommended pressure range.
- The red LED indicating priority alarms shall remain on.
- It shall be possible to hear the alarm speaker test.
- A RESPIRATORY CIRCUIT SUITABLE FOR THE TYPE OF PATIENT SHALL BE MOUNTED AND THE OUTPUT SHALL BE OBSTRUCTED.



WARNING

- ***Never start the autotest procedure with the ventilator connected to the patient.***
-

If you heard the speaker test, answer YES. Otherwise, answer NO. In this case the equipment shall not be used and therefore, it will be inoperative until a new autotest is carried out, that is, until the equipment is restarted.

To start the autotest, press the confirm button or, if you want to start the ventilation immediately, press the ventilation button.



WARNING

- ***Never respond YES if you could not hear the speaker test under the risk of malfunctioning of the priority alarms during ventilation.***
 - ***Always perform the autotest procedure before connecting the ventilator to a patient.***
-

The autotest checks all important items for proper ventilation:

- Gas regulated pressure
- Proportional Valves
- Flow Sensors
- Expiratory valve
- O₂ cell
- Leakage
- Resistance and compliance of the respiratory circuit
- Flowair system¹

At the end of the leak test, the RESPIRATORY CIRCUIT OUTPUT SHALL BE UNOBSTRUCTED to perform the resistance test and, where appropriate, the proximal flow sensor test.



ATTENTION

- ***Do not forget to unobstructed the circuit output before the resistance test.***
-

¹ Only for Fleximag Max 300 and Fleximag Max 700



ATTENTION

- *If there is a recurring failure in one or more items of the autotest, suspend use of the equipment until the problem is resolved.*
-

The autotest can be canceled at any time by the operator. To perform only some of the tests, you must enter the RESTRICTED menu.



ATTENTION

- *The restricted menu is accessible only with a password.*
-

Once the autotest process is complete, press VENTILATION to go to the ventilation screen.




ATTENTION

- *Position the patient close to the Fleximag Max so that the breathing circuit is not stressed.*
 - *When using the respiratory circuit with water traps, keep them at a lower level to ensure correct operation of the collectors.*
-

6.2 Restricted menu

This menu is password-only accessible and allows you to make advanced settings, calibrations and autotests, and allows you to view system information, battery and status.

6.2.1 Restricted menu



6.2.1.1 Patient

Clicking this button will return you to the patient screen.

6.2.1.2 Autotest

In this menu, you can perform the tests, and cancel and view the results individually. This menu displays the air pressure, O₂ information, date and time of the last test performed.

6.2.1.3 Calibration

In this menu, it is possible to perform the calibration of the distal flow sensor, expiratory valve and oxygen cell. This menu displays the air pressure, O₂ information, date and time of the last test performed.

To perform the calibrations, occlude the breathing circuit, select the items you want to calibrate, and press the calibrate button.

6.2.1.4 Status

In this menu, network connection information is displayed or, depending on the type of access, the engineering screen is displayed.

6.2.1.5 Battery

Percentage and battery status information is displayed in this menu.

6.2.1.6 System

In this menu you can set the date, time, language, brightness, day or night mode and audio volume.



ATTENTION

- *Make sure that the volume adjusted for the audio of the alarms is compatible with the distance that the equipment will be from the clinicians.*

6.2.1.7 Advanced

In this menu you can adjust time and O₂+ concentration, inspiratory and expiratory hold time, keypad lock time, alarm mute time, type of controlled volume adjustment, type of flow sensor and pressure unit.

6.2.1.8 Ventilation

This menu gives access to the main ventilation screen of the patient.


To ventilate, remove the ventilator from standby mode by pressing and holding the respective button located in the quick access function area for 2 (two) seconds.



WARNING

- *Whenever it is restarted, the ventilator will enter **STAND BY** mode and in this condition the patient will not be ventilated.*

6.3 Main screen



6.3.1 Battery status area

Percentage and battery status information is displayed in this menu.

6.3.2 Patient information area

Patient information such as ID, name, initials, last name, sex, height, and volume / weight are set in this menu.

6.3.2.1 Volume x weight patient definition

Once a patient is selected, the ventilator automatically estimates an adequate tidal volume. However, to obtain the best current volume is important to know the ideal weight for each patient.

It is possible to obtain the ideal weight for pediatric and adult patients according to sex and height. To do so, simply select the patient's sex and then adjust the height value, so that the equipment calculates your ideal weight. The parameter setting volume by weight (ml per kilogram) complements the information necessary to better fit the current volume.

For neonatal patient, the ventilator estimates the adequate tidal volume according to the body weight of the patient. To do so, simply select the patient's sex and then adjust the weight value.

Observation

- *To reduce the risk of lung injury, the ventilator uses the ideal weight as a reference for adjusting ventilatory parameters. If you want to change the parameters set by the ventilator, simply select the parameter and change it.*

6.3.3 Historical events area

This menu displays the history of events related to alarms, ventilation, maneuver, settings, battery, calibration, actions and tests.

6.3.4 Alarm setting area

For each of the alarms directly related to the ventilation process, there are one or two limits (high value and low value) to be adjusted. These limits are set directly in the ALARM menu.

To make these adjustments, touch the button corresponding to the alarm to be adjusted and it will be selected (the color will be modified), allowing to change its value using the turns and confirms button. To confirm the set value, press the button corresponding to the alarm again or press the turn and confirm button (ENTER).

In this menu it is also possible to set the maximum allowed apnea time of the patient, which will determine the start of backup ventilation (backup).

There is also the option of automatic adjustment of the alarms, to enable it is necessary that the ventilator is not in STAND BY and that the ventilation is stabilized, aiming for greater patient safety.

6.3.5 Mute alarm button

By pressing this button, the alarms are silenced for the time set in the ADVANCED menu. If a new alarm occurs while the mute button is active, the mute is disabled and the alarm sounds.

6.3.6 Manual cycle button

Manually trigger an inspiratory cycle, depending on the ventilatory mode selected.

6.3.7 Inspiratory hold button

Pressing this button pauses the inspiratory time for the time set in the ADVANCED menu. This function allows the performance of inspiratory suspension maneuvers, used in the case of chest X-rays. During the execution of this maneuver, the apnea alarm is not displayed. Only the message "Insp Pause" is displayed.

6.3.8 Expiratory hold button

Pressing this button pauses the Expiratory time for the time set in the ADVANCED menu. This function allows for extended expiration time maneuvers (prolong expiration time). During the execution of this maneuver, the apnea alarm is not displayed. Only the message "Exp Pause" is displayed.

6.3.9 O₂+ button (O₂ flush)

When this button is pressed, an O₂ flush is performed with the concentration and time set in the ADVANCED menu. During this time, the high FiO₂ alarm will be inhibited. This feature can be used for pre- and post-aspiration of airway secretion procedures and is available in all ventilatory modes. Use with an aspiration catheter may be performed in any mode and does not require a specific adjustment.

6.3.10 Nebulizer / TGI button

By pressing this button, nebulization or TGI is performed according to the setting in the RESOURCES menu.

6.3.11 Graphic freeze button

When you press this button, the graphics are frozen.


6.3.12 Menu access area

Pressing this button displays the following menus: Layout, Resources, Maneuver, Trend, Calibration, Advanced, Status and Engineering mode.


6.3.12.1 Layout

Pressing this button displays the available graphic layouts: 2 waveforms, big numbers, monitor, 3 waveforms, lung and 2 loops.


2 waveforms




Big numbers




Monitor




Lung



3 waveforms



2 loops



Legend:

A – Graphics: SpO₂, CO₂, CO₂(%), Flow, Volume, Pressure

B – Monitored parameters: SpO₂, Pulse, Perf, PVI, EtCO₂, RR, iCO₂, Ppeak, PEEP, Rrate, O₂, Vte, MV, Vti, Pplat, Pmean, I:E, Ti, Te, Vte sp, MV sp, RR sp, Ti/Ttot, RSBi, WOBi, Ri, Re, C dyn, C stat, Leakage (%), Leakage (F), Tc, E, iPEEP, O₂ consumption, Driving Pressure

C – Bargraph

D – Lung

E – Loops: PxV, PxF, VxF, VxCO₂, VxFCO₂

6.3.12.2 Resources

Pressing this button displays the following resources: NIV, External Auxiliary Pressure, Tube Compensation, Sigh, Neb-TGI and Humidifier.

6.3.12.2.1 **Noninvasive Ventilation (NIV)**

Noninvasive ventilation (NIV) refers to the application of ventilatory support without invasive methods of the airways, such as orotracheal intubation or tracheostomy. Nasal or oronasal masks are the most frequently used interfaces for the application of NIV in a hospital environment.

When applying NIV in pressure-controlled modes, the pressure value must not be set to 0 (ZERO) and the pressure drop cycle trigger must be active. Flow trigger remains off.

In NIV the ventilator automatically compensates for higher leakage flows and ignores high minute volume, high current tidal and flow sensor checks.

NIV is available for all ventilatory modes.



WARNING

- *The default values are only an initial reference.*
- *Readjust ventilation parameters as required by the patient.*
- *Use the appropriate mask for each type of patient to avoid excessive leakage.*
- *The exhaled volume of the PATIENT may differ from the volume exhaled measured due to leaks in the mask.*
- *For non-invasive ventilation, use expiratory carbon dioxide concentration measuring devices in accordance with ISO 80601-2-55 (See Chapter 11 - IRMA CO₂ Sensor (optional)).*



Observation

- *Flow triggering is disabled during noninvasive ventilation.*
- *Controlled or support pressure (ΔPS) is a value above PEEP and can be adjusted between + 5 cmH₂O and P_{MAX} .*
- *The continuous flow, which apparently 'leaks' through the expiratory valve is normal and serves to reduce the response time of the patient's ventilation control system.*

6.3.12.2.2 External Auxiliary pressure

On the ventilator's front panel, there is a channel for external auxiliary pressure measurement. To use this channel, connect one end of a suitable tube to the auxiliary pressure nozzle (P.AUXILIARY) and the other end to the pressure channel you wish to measure.

It is possible to use this feature with an esophageal balloon or to obtain the pressure of the carina, among other forms.

Activating this option will draw a new pressure curve on the screen with its instantaneous values. The patient's pressure curve continues to be drawn normally.

6.3.12.2.3 Tube compensation

The main purpose of this feature is to compensate the work imposed on the patient by the endotracheal tube.



WARNING

- ***Wrong selection of the endotracheal tube type or gauge can cause injury to the patient.***

In the Tube Compensation window, select the intubation mode and then adjust the endotracheal diameter and compensation percentage.

After making sure that the adjustment is appropriate to the patient, close the setup window and activate the tube compensation.

By activating this option, a new patient pressure curve will be drawn, with its instantaneous values.

This pressure value is estimated based on algorithms that consider the pipe diameter and the percentage of compensation.

The patient's pressure curve, read by the equipment, remains traced and usually tends to be greater than the estimated pressure.

6.3.12.3 Maneuver

The Maneuver menu is available for all patients, except for neonatal patients. This menu basically provides some processes that aid diagnosis, providing data related to the patient's respiratory mechanics.

6.3.12.3.1 P0.1

The P0.1 index can be considered as the pressure drop, below the baseline pressure, that is generated by the patient's inspiratory effort and measured in the first 100 ms of the onset of the inspiratory phase.

As soon as the maneuver is started, the ventilator will step in spontaneously to identify the inspiratory efforts of the patient.

Whenever a cycle is triggered by patient inspiratory effort, the ventilator will calculate P0.1 and display it in the table on the top panel to the left of the chart.

6.3.12.3.2 Slow Vital Capacity

A slow vital capacity is the capacity of expiration after a maximum inspiration by the patient and serves as parameter for the evaluation of its ventilatory reserve.

To obtain this parameter, it is necessary that the patient is aware, since their collaboration is paramount.

To perform the maneuver, the ventilator will enter CPAP mode, without pressure support (delta PS = 0).

The patient should then perform successive breaths, extending the inspiratory phase to the maximum, and then exhale slowly, to the maximum extent possible.

6.3.12.3.3 P/V Flex

The inflection points of the PV curve (pressure x volume) can be used to obtain the most suitable adjustment values for PEEP and plateau pressure.

Through this feature, the lower and upper inflection points are obtained, the first being the basis for determining the optimal level of PEEP, while the second serves as a parameter for the appropriate maximum pressure and volume level, preventing pulmonary hyperdistention.

For this procedure, the patient should be intubated and without exerting respiratory efforts.

Initially, the user shall adjust the appropriate pressure, volume and flow values to the patient and then press START.

The equipment will provide adjustable constant flow and will pass instantaneous pressure and volume monitoring.

As soon as one of the monitored parameters is reached, the flow is reset, and the ventilator will cycle again normally.

6.3.12.3.4 *Pi Max*

The Pi Max or NIF (negative inspiratory force) is only available for spontaneous modalities and is used to assess inspiratory muscle strength during the weaning process of mechanical ventilation.

Before beginning the maneuver, orient the patient to make the maximum inspiratory effort possible as soon as it is requested.

For this maneuver it is necessary to keep the START button pressed during the entire process.

Once the maneuver is initiated, the PEEP value will be temporarily reset, and the inspiratory branch occlusion will occur.

The negative pressure values relative to the patient's inspiratory effort will be measured to the limit of -60 cmH₂O.

The best value reached will always be displayed in the table to the left of the pressure chart, i.e. the largest pressure drops identified.

6.3.12.3.5 *Gas Trapping Volume*

The trapped volume is only available for assisted-controlled modes.

An unwanted volume of air may eventually become trapped in the lungs in cases of pulmonary hyperinflation or when the interval between breaths is not enough for complete patient expiration to re-establish the balance of the respiratory system. Most notably, when the presence of intrinsic PEEP is detected.

To accomplish this procedure, it is ideal that the patient does not exert respiratory efforts. Therefore, it is recommended that the operator advises the patient if the patient is conscious.

To perform the maneuver, simply press START and then, with each cycle, the ventilator will compare the target volume value (desired) and the total volume value reached.

6.3.12.4 Trend

By pressing this button, the trend settings will be displayed. The trend holds up to 240 hours of ventilation.

6.3.12.5 Calibration

Pressing this button displays the available calibrations: O₂, CO₂ and touch. To perform the calibration, put the equipment in Standby mode and select the item you want to calibrate.

It is recommended to calibrate some components before use in the following situations:

6.3.12.5.1 O₂ cell (galvanic cell only)

- Cell replacement.
- The monitored concentration values (FiO₂) do not seem correct.
- The lower and upper limits do not reach 21 and 100% O₂, respectively.
- Change of patient.

Observations

- *To access the calibration screen, press the CALIBRATION button on the ventilator home screen.*
- *It is not necessary to dispose calibration gas.*

6.3.12.6 Advanced

Pressing this button displays the following settings: O₂+, ins and exp hold, keypad lock, mute alarm, controlled volume adjustment, flow sensor, and pressure unit.

6.3.13 Lock Button

Accidental change protection system. Locks or unlocks the touch screen. When the display commands are locked, press this key to unlock them IMMEDIATELY.

To lock again, simply press this key once or wait for some time without touching the screen and the screen will lock automatically, the time it takes for the screen to lock is set in the general settings.

The time lock can be turned off in the advanced menu.

6.3.14 Stand by button

Pressing this button turns the standby mode on or off. In standby mode, alarms are interrupted, and ventilation is paused. For security reasons, to activate or deactivate the standby mode, press the button for 2 seconds.

6.3.15 Bargraph or monitored parameters area

According to the layout, the pressure bargraph or 6 monitored parameters will be displayed.

6.3.16 Graphics area

In this area the adjusted graphics layout is displayed. The available graphics are:

- Pressure x Time Curve
- Flow x Time Curve
- Volume x Time Curve
- Pressure x Volume Loop
- Volume x Flow Loop
- Pressure x Flow Loop
- Volume x CO₂ Loop
- Volume x FCO₂ Loop
- CO₂ x Time curve ⁽¹⁾
- SpO₂ x Time curve ⁽¹⁾
- Instant pressure bar graph with numerical indicator of peak, plateau or instantaneous pressure

(1) This graph option is only available when an external sensor (oximeter or capnograph) is connected.

The monitored parameter display option is still available, where up to 24 of the 36 monitored parameters are displayed.

6.3.17 Ventilation Mode Setting Area

Pressing this button displays the ventilator modes available.

6.3.17.1 Ventilatory modes available

Table 5 - Ventilatory modes

Mode	Backup	Backup Mode ⁽¹⁾	
		Neo	Ped and Adu
VCV	✓	—	Auto
PCV	✓	Auto	Auto
PRVC	✓	—	Auto
PLV	✓	Auto	—
V-SIMV	✓	—	Auto
P-SIMV	✓	Auto	Auto
CPAP/PS	✓	Adjustable PLV + Auto	Adjustable VCV and PCV + Auto
DualPAP	✓	Adjustable PLV + Auto	Adjustable VCV and PCV + Auto

Mode	Backup	Backup Mode ⁽¹⁾	
		Neo	Ped and Adu
APRV	✓	Adjustable PLV + Auto	Adjustable VCV and PCV + Auto
MMV	✓	—	Adjustable VCV and PCV + Auto
VS	✓	—	Adjustable PRVC + Auto
Nasal CPAP	✓	Adjustable PLV + Auto	—
NIV	✓	Adjustable PLV + Auto	Adjustable VCV and PCV + Auto
VG	✓	Auto	—
O ₂ Therapy	X	—	—

Observations

- For modes in which the backup is set to "Auto", whenever the set apnea time is reached, the ventilator starts a ventilatory cycle, this configuration is based on the settings of the current ventilator mode.

6.3.17.2 Adjusting ventilatory modes

To change or reconfigure a ventilatory mode, simply touch the button that indicates the active mode in the upper right corner of the screen. This button will change color and then the ventilator modes setting screen will be loaded.

Observations

- The ventilatory modes available will be determined according to the ventilator model and the selected patient (see **Table 5**).

To select a ventilatory mode, simply touch the tab with the acronym of the desired mode. All the adjustable parameters required for this ventilatory mode, including those for backup ventilation, will be displayed.

Observations

- *Adjustment of backup ventilation parameters (backup) is only available in spontaneous ventilatory modes. In the other modes, the backup ventilation is automatic and considers the parameters adjusted for the ventilatory mode itself.*
- *When the operator sets pressure or flow trigger, the controlled modes (VCV, PCV, PRVC and PLV) will be monitored. In this case, such information will be displayed on the active mode button.*

After adjusting the parameters, for them to be activated, it is necessary to press the CONFIRM button.

To cancel the settings and to continue with the previous settings, including the ventilatory mode, simply press the CANCEL button. This way the ventilator will ignore the settings made on this screen and will return to the main screen.

ATTENTION

- *If the new value is not confirmed, it will be discarded after the screen lock time.*

6.4 Calibrations

It is recommended to calibrate some components before use in the following situations:

6.4.1 Distal Flow Sensor

- Replacement of the distal flow sensor;
- Active alarm with message "CHECK FLOW SENSOR".

6.4.2 Integrated valve

- Replacement of integrated valve;
- Replacing the diaphragm;
- Incorrect control of PEEP;
- Excessive leakage.

6.4.3 O₂ cell (galvanic cell only)

- Cell replacement;

-
- The monitored concentration values (FiO₂) do not seem correct;
 - The lower and upper limits do not reach 21 and 100% O₂, respectively;
 - Change of patient.

Observation

- *To access the calibration screen, access the menu on the ventilation screen or restricted on the ventilator home screen.*

6.5 Automatic compensation of the respiratory circuit

For volumes and pressures, delivered and monitored, leaks are always considered because of circuit compliance and resistance.

In the event of failure of the compliance and resistance autotest, the efficiency of the compensations is maintained if the circuit used is compatible with the values defined in Table 41 and Table 42.

7 Troubleshooting

In this topic, the main occurrences and the possible solutions are presented.

Occurrence	Probable causes	Proposed solution
Low battery alarm	Depletion of internal battery charge.	<i>Connect the equipment to the electrical outlet or provide other means of ventilatory support.</i>
	Internal battery charging system failure, even with power.	<i>Contact technical assistance for service.</i>
Disconnect alarm	Disconnection in the respiratory circuit.	<i>Locate the disconnect point and firmly connect all parts of the breathing circuit.</i>
	Lack of inspiratory flow.	<i>Check for inspiratory flow and increase if necessary.</i>
	Alteration of the patient's respiratory mechanics.	<i>Establish adequate parameters for ventilatory support.</i>
	Integrated valve diaphragm placed incorrectly or damaged.	<i>Replace or return the diaphragm to the correct position.</i>
	Failure of the electronic control system.	<i>Contact technical assistance for service.</i>
Communication failure alarm	Electronic fault.	<i>Contact technical assistance for service.</i>
High pressure alarm	Alteration of the patient's respiratory mechanics.	<i>Establish adequate parameters for ventilatory support.</i>
	Obstruction in the expiratory branch of the respiratory circuit or in the integrated valve.	<i>Clear the circuit of any obstruction or reposition the integrated valve diaphragm.</i>
	Obstruction of the patient's airway.	<i>Unclog or aspirate the patient's airway.</i>
	Monitored inspiratory pressure is higher than expected.	<i>Check the inspiratory pressure (absolute) setting, which is the sum of the controlled (relative) pressure and PEEP.</i>
	Tube Compensation is on.	<i>Turn off or reconfigure tube compensation.</i>
Low pressure alarm	Alteration of the patient's respiratory mechanics.	<i>Establish adequate parameters for ventilatory support.</i>
	Excessive leakage in the respiratory circuit.	<i>Locate the leak and correct.</i>
Electrical network alarm	Disconnected power cord.	<i>Properly connect the power cord to the equipment or in case of in-hospital transport, check that there is enough battery power.</i>
	Power failure.	<i>Try to reestablish a power source.</i>
Inoperable Alarms	Electronic failure.	<i>Contact technical assistance for service.</i>
Incorrect control of PEEP	Calibration of the integrated valve.	<i>Restart the equipment and calibrate the integrated valve.</i>

Occurrence	Probable causes	Proposed solution
Curves and / or ventilation loops appear with inadequate scaling or tracing speed.	Automatic adjustment of scales or graphic speed off.	<i>Tap on the graph screen and select the automatic adjustment in the corresponding window.</i>
Curves and trend values do not appear or are incorrect.	Failed to set the system clock.	<i>Contact technical assistance for service.</i>
Equipment does not start ventilation	Equipment in standby mode.	<i>Press the STAND BY button for 2 seconds to remove the ventilator from standby mode.</i>
Proportional valve test failure.	Gas network pressure below the minimum limit.	<i>Check and readjust the network pressure to the specified range.</i>
Distal flow sensor test failed.	The respiratory circuit output was not obstructed.	<i>Restart the ventilator and redo the autotest with the respiratory circuit obstructed.</i>
Proximal flow sensor test failure and resistance.	Respiratory circuit output was not open.	<i>Restart the ventilator and redo the autotest remembering to open the respiratory circuit when prompted.</i>
Integrated valve test failure	The respiratory circuit output was not obstructed.	<i>Restart the ventilator and redo the autotest with the respiratory circuit obstructed.</i>
	Integrated valve diaphragm	<i>Reposition the integrated valve diaphragm, restart the equipment, and redo the autotest.</i>
No maneuvers available	Neonatal patient selected.	<i>Maneuvers available only for pediatric and adult patients.</i>
	Equipment in standby mode.	<i>Press the STAND BY button for 2 seconds to remove the ventilator from standby mode.</i>
	Ventilatory mode not compatible with the maneuver.	<i>Set an appropriate ventilator mode according to the desired maneuver.</i>
It is not possible to activate automatic alarm adjustment.	Equipment in standby mode.	<i>Automatic adjustment of alarms can only be calculated with the ventilator cycling. Press the STAND BY button for 2 seconds to remove the ventilator from standby mode, wait for the ventilation to stabilize and then activate the automatic adjustment.</i>
Events cannot be viewed.	There is still no loaded trend curve.	<i>Select a trend range before attempting to query the events.</i>
The setting of the parameters returns to the previous value.	Adjustment was not confirmed.	<i>Confirm the setting by pressing the turn and confirm button or by tapping on the parameter that has been reset.</i>
Inspiratory or expiratory hold do not end as soon as the respective button is released.	The minimum set pause time is longer than desired.	<i>Press the ADVANCED MENU button and reset the minimum pause time.</i>
Alarm sound is off	Alarm silence time activated.	<i>Turn off the alarm silence.</i>
	Electronic failure.	<i>Contact technical assistance for service.</i>
Backup ventilation is not activated.	Apnea time alarm off.	<i>Set an interval for the apnea alarm.</i>



WARNING

- **Never operate the equipment if a problem cannot be solved.**

8 Cleaning, Disinfection and Sterilization

It is important to establish a routine for the cleaning, disinfection or sterilization of the equipment and its components.

The following are the main ways of cleaning, disinfecting or sterilizing according to the characteristics of each component and the equipment.



WARNING

- *Before use, perform the cleaning, disinfection and sterilization procedures as specified in this manual.*

8.1 External Parts

External ventilator surfaces of Fleximag Max shall be cleaned with a clean, soft cloth moistened with the enzymatic detergent.



Observations

- *Make sure that the equipment is turned off to perform cleaning of the display.*
- *Make sure that no residue builds up in the equipment connections.*
- *For cleaning, do not use non-compatible polymer products.*

8.2 Components

Components that come in contact with respiratory gases shall be periodically disassembled for cleaning, disinfection or sterilization, which are: proximal flow sensor, silicone line and integrated valve.

8.2.1 PROXIMAL flow sensor, silicone line and integrated valve

Parts of silicone shall be cleaned according the following steps:

8.2.1.1 Wash

- a) Always use potable water for this procedure;

-
- b) Use neutral and enzymatic detergent. Dilution should be performed as recommended by the manufacturer.
 - c) Immerse the entire body of the flow sensor and the silicone line in the detergent solution, keeping the solution in contact with the accessories for at least 3 minutes;
 - d) The external parts of the parts should be cleaned with a clean, soft cloth moistened with the enzymatic detergent. The internal parts must be cleaned by immersion.

8.2.1.2 Rinse

- a) Always use potable water for rinsing;
- b) Thoroughly rinse the external surface of the accessories with potable water.
- c) Rinse the internal surface by injecting potable water under pressure at least 5 times.



ATTENTION

- *Do not use for cleaning or disinfection phenol (> 5%), ketones, formaldehyde, hypochlorite, chlorinated hydrocarbonates, aromatic hydrocarbonates, inorganic acids and quaternary ammonium compounds.*
- *Never use saline solutions, especially sodium hypochlorite (bleach) and saline, disinfectants, hydrogen peroxide for cleaning or rinsing the accessories.*

8.2.1.3 Drying

Drying of the external parts should be done with a clean, soft and dry cloth and the drying of the internal parts should be done so that the solution drains by gravity.

Composition of the enzymatic detergent

Amylase, Lipase, Protease, Carbohydrase, Isopropyl Alcohol, Control Agent and pH, Preservatives, nonionic surfactant, coadjuvants, stabilizers, additives and water. Minimum amylolytic activity: 0.02 ua.ml⁻¹.min.⁻¹; Minimal proteolytic activity: 0.12 up.ml⁻¹.min.⁻¹; Surfactant content: Ethoxylated alcohol: 3.5%. Pure PH: 7.0.

8.3 Disinfection

8.3.1 External Parts

The external part should be disinfected using a clean cloth moistened with alcohol 70 °.

8.3.2 Exhalation valve, proximal flow sensor and silicone line

After cleaning, the items should be disinfected with alcohol 70°. The external part should be disinfected using a clean cloth moistened with alcohol 70 ° and internal part by immersion.

After disinfection, the external parts should be dried with a clean, soft and dry cloth and the internal parts should be dried so that the solution drains by gravity.

8.3.3 Distal flow sensor (Envitec SpiroQuant A+)

It is recommended to disinfect the flow sensor before use. For disinfection:

1. Disconnect the distal sensor from the exhalation valve and the connector cable and wait 30 minutes.
2. Immerse the sensor in 70% ethanol solution for 1 hour
3. Allow to dry naturally for 30 minutes in the environment before reassembly.



ATTENTION

- *The distal flow sensor must not be sterilized by steam.*
- *For the internal disinfection of DISTAL flow sensor, do not use tools that may generate mechanical forces, such as air or water jet, at the risk of damaging the filament.*
- *Drying of DISTAL flow sensor should occur naturally in ambient air, so avoid the use of compressed air or dryer.*
- *If there is a possibility of the sensor remain infected, replace it promptly.*

8.3.4 IRMA CO₂ sensor

The IRMA CO₂ sensor may be disinfected with a cloth moistened with 70% ethanol.



ATTENTION

- *The airways adapters IRMA CO₂ sensor are unsterile supplies, so the steam sterilization can damage these accessories.*

-
- *The airway adapter should not be reused.*
 - *The reuse of a disposable adapter may cause cross infection.*
 - *Never sterilize or immerse the IRMA CO₂ sensor in liquid.*
-

8.3.5 Oximetry sensor (oximeter)

Use the cleaning instructions as section 8.1.



ATTENTION

- **Caution at cleaning:**
 - *Do not steam-sterilize, pressure-sterilize, or gas-sterilize the Oximeter*
 - *Do not wet up or dive the monitor into any liquid.*
 - *Use cleaning solutions with moderation. Excess solution may run into the monitor and cause internal damages to the components.*
 - *Do not touch, press, or scrub the display panel with abrasive cleaning components, brushes, cleaning instruments, and do not leave it in contact with anything to scratch the panel.*
 - *Do not use solutions derived by petroleum or acetones, or other abrasive solvents to clean the oximeter. These substances attack the device's materials and may result in its failure.*
- **Inaccurate measures may be caused by:**
 - *Incorrect sensor application or use*
 - *Significant levels of dysfunctional hemoglobin. (e.g., carboxyhemoglobin or methemoglobin)*
 - *Intravascular dyes such as indocyanine green or methylene blue.*
 - *Exposure to excessive illumination, such as surgical lamps (especially ones with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, or direct sunlight (exposure to excessive illumination can be corrected by covering the sensor with a dark or opaque material)*
 - *Excessive patient movement.*
 - *Venous pulsations.*
 - *Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.*
- **Loss of pulse signal can occur in any of the following situations:**
 - *The sensor is too tight*
 - *There is excessive lighting from light sources with a surgical lamp, a bilirubin*

lamp, or sunlight.

- *The blood pressure cuff is inflated at the same extremity of where the SpO₂ sensor attached.*
 - *The patient has hypotension, severe vasoconstriction, severe anemia, or hypothermia.*
 - *There is an arterial occlusion next to the sensor.*
 - *The patient is under cardiac arrest or shock.*
-

8.4 Sterilization

- Components that come in contact with respiratory gases shall be dismantled for cleaning and sterilization;
- Do not use abrasive cleaning agents;
- Do not use alcohol to clean plastic parts;
- Do not immerse the equipment in any liquid.

8.4.1 Steam Sterilization

The following accessories shall be sterilized in hot steam (135°C for 5 minutes).

Table 6 – Sterilizable accessories

Description	Sterilization cycles (lifetime)
PROXIMAL flow sensor	50
Silicone Lines	50
Integrated valve	50
Integrated valve diaphragm	50

8.5 Important Advices



WARNING

- *Prior to first use, the equipment and its components shall be cleaned, disinfected and sterilized as specified.*
- *All parts of Magnamed ventilators that come in contact with patient fluids (e.g. breathing circuit) and that are potentially contaminated, are called semi-critical, and*

shall be endure a high-level disinfection or sterilization process before disposal (at the end of their useful lives) or before it is sent in for maintenance.

- *When sending the ventilator for maintenance or repair, observe the process of disinfection RIGOROUSLY.*
 - *Do not immerse the sensor in water, solvents, or cleaning solution (sensors and connectors are not waterproof). Do not sterilize by irradiation, steam, or corrosion sterilization. See the cleaning instructions in the instructions for use for Masimo LNOP®/LNCS® reusable sensors.*
 - *Do not use damaged patient cables.*
-



ATTENTION

- *Removable Magnamed accessories and components that exhibit damage or signs of wear should be replaced to prevent use.*
 - *Packaging of non-sterile devices (breathing circuits, expiratory valves, and connectors) is designed to keep these products at the proper cleaning level to be sterilized prior to use and to minimize microbial contamination.*
 - *Cleaning*
 - *Do not sterilize the oximeter by pressure, gas or steam.*
 - *Do not soak or immerse the monitor in any liquid.*
 - *Use the cleaning solution sparingly. The excess may leak into the monitor and cause internal damage to the components.*
 - *Do not use solutions derived from petroleum, ketones or other harsh solvents to clean the oximeter. These substances attack the material and can cause problems to the device.*
-



Observations

- *Do not use abrasive cleaners.*
 - *Do not use alcohol to clean the plastic parts, except when specified.*
 - *Do not immerse the ventilator in any liquid.*
-

8.6 Processing Methods


Table 7 – Processing methods

Component	Processing method		
	Neutral Enzymatic Detergent	70% Alcohol	Steam Sterilization 135°C for 5 min
Ventilator surface	✓	x	x
Touch screen	✓	✓	x
Silicone pressure line	✓	✓	✓
Integrated valve	✓	✓	✓
Diaphragm	✓	x	✓
Proximal Flow Sensor (Adu, Inf and Neo)	✓	✓	✓
Distal flow sensor (heated filament)	x	✓	x
SpO ₂ sensor	✓	x	x
EtCO ₂ Sensor	✓	x	x

9 Preventive maintenance



WARNING

- The symbol  on the ventilator screen indicates that the equipment has entered the preventive maintenance period. For ICU equipment this period is 5,000 hours or 12 months, whichever occurs first.
- Failure to perform maintenance may affect the safety and performance of the ventilator.
- Maintenance shall be performed as indicated by the manufacturer and only by an authorized technician service. The noncompliance will result in loss of warranty and manufacturer's obligations related to the ventilator.
- Schedule preventive maintenance only with the authorized Magnamed Service.
- Before sending the equipment to the service technician, observe the cleaning and disinfection process **RIGOROUSLY**.

9.1 Verifications

The following verifications shall be made daily and whenever the equipment is to be used:

- AC/DC converter power cable integrity;
- Operation of the alarm system, including audio;
- Air/O₂ Filters installed and unobstructed;
- LCD display;
- Batteries charged;
- Touchscreen;
- Turns and confirm Button;
- Correct installation of the respiratory circuit (including the diaphragm of the integrated valve);
- Mesh filter installed.



WARNING

- The daily check must be performed with the ventilator disconnected from the patient.
- The Fleximag Max, its parts and accessories should not be serviced during use.

9.2 Preventive Maintenance Schedule

Magnamed recommends performing preventive maintenance of ICU ventilators with your authorized network throughout the country.

Below spreadsheet with the schedule of maintenance and exchange of preventive parts. If you need more details, contact the MagnaService (Magnamed Technical Assistance).

Table 8 – Preventive maintenance schedule

ITEM	Period									
	5000 h or 1 year	10000 h or 2 years	15000 h or 3 years	20000 h or 4 years	25000 h or 5 years	30000 h or 6 years	35000 h or 7 years	40000 h or 8 years	45000 h or 9 years	50000 h or 10 years
Assessment according to procedures MAGNAMED	X	X	X	X	X	X	X	X	X	
Evaluate replacement of consumables	X	X	X	X	X	X	X	X	X	
LI-ION Battery Pack		X		X		X		X		
O ₂ cell (galvanic)		X		X		X		X		
Disposal of equipment										X

9.3 Consumable items

Table 9 – Consumable items replacement

ITEM	Period
Silicon line	2 years or after 50 hot steam sterilization cycles
Diaphragm	2 years or after 50 hot steam sterilization cycles
Adult proximal flow sensor	5 years or after 50 hot steam sterilization cycles
Pediatric proximal flow sensor	5 years or after 50 hot steam sterilization cycles
Neonatal proximal flow sensor	5 years or after 50 hot steam sterilization cycles
Integrated valve	5 years or after 50 hot steam sterilization cycles
Mesh filter	500 h of use
Air/O ₂ inlet filter	Replace if it is obstructed

9.4 Internal Batteries

These batteries are responsible for maintaining the operation of the equipment even in the absence of electricity and its duration in normal operation is specified in chapter 10.



WARNING

- *For there to be enough battery power during a power outage, it is important that the equipment remains connected to a mains power supply **WHENEVER POSSIBLE**.*



ATTENTION

- *For the capacity of the batteries in normal operation to be full, they shall be replaced as indicated in the technical specification.*
- *Replacement of the internal batteries shall be performed by trained and qualified personnel.*

9.5 Water Collectors with Coalescing Filter






ATTENTION

- *Do not expose the filter container to non-polycarbonate compatible materials.*
- *Replace the filter when it is clogged so that it does not slow down the inlet flow of the equipment.*

To remove accumulated water, simply press the pin found on the bottom of the collector.

To perform filter replacement, consider the following instructions:

Instruction	Image
Remove the screws with an Allen wrench	


Instruction	Image
Remove the manifold with the O-ring	
Unscrew the filter and thread the new filter	

9.6 Mesh filter

ATTENTION

- *The mesh filter is used to protect the equipment. To protect the patient from bacteria and virus contamination, use the HEPA filter as specified in chapter 10.29.*

The mesh filter is used to protect the equipment from particulate matter in suspension in the environment. To perform filter replacement, consider the following instructions:

Instruction	Image
<p>Turn the filter to the left side to remove.</p> <p>Introduce a new filter and turn to the right side to close.</p>	

9.7 O₂ Concentration Cell

This ventilator is equipped with O₂ monitoring means for measurements of inspiratory oxygen concentration according to ISO 80601-2-55. This equipment has two ways of measuring the concentration of oxygen:

Galvanic Cell - generates an electric signal proportional to the concentration of oxygen in the gas mixture administered to the patient and the intensity of this electrical signal is due to the chemical reaction. This measuring medium is consumable, and the life of the cell as specified by the original manufacturer is 10,000 hours at 100% O₂, i.e. greater than one year of continuous use. However, we recommend switching during preventive maintenance on a 24-month or 10,000-hour schedule (whichever comes first).

Paramagnetic Cell - generates an electric signal proportional to the concentration of oxygen in the gas mixture administered to the patient and the intensity of this electrical signal is due to the torque created in the magnetic arrangement of the cell. This sensor uses the paramagnetic susceptibility of oxygen that distinguishes it from other gases. This measuring medium is not consumable.



ATTENTION

- *The galvanic cell for measuring the oxygen concentration should be replaced as indicated in chapter 10.21.*
- *Replacement should be performed by trained and qualified personnel.*
- *The galvanic O₂ cell undergoes degradation of less than 1% per month in the*

measurement accuracy.

10 Technical specifications

10.1 Classification

10.1.1 Class I Equipment

According to NBR - IEC - 60601, internally energized, type B for continuous operation. IP31 splash proof equipment.

10.1.2 Protection class of applied parts

Patient circuit and flow sensor type B, oximetry sensor and capnography sensor type BF.

10.1.3 CE – Class IIB

According to Annex IX to Directive 93/42 / EEC, Rule 9 - All therapeutic devices that are active because of their characteristics are likely to provide or exchange energy for the body in a potentially dangerous manner, considering the nature, density and location of energy, in which case they belong to class IIB.

10.1.4 ANVISA – Class III

According to RDC 185/01 - Classification rule 11 - All active medical products intended to administer or extract medicinal products, body fluids or other substances from the organism fall within Class II, unless this is done in a potentially dangerous manner, considering the nature of the substances, the part of the body involved and the mode of application, in this case fall under Class III.

10.1.5 FDA – Class II

According to the Code of Federal Regulations (CFR), Title 21 - Food and Drugs, Chapter I - Food and Drug Administration, Department of Health and Human Services, subchapter H - medical devices, Part 868 - Anesthesiology Devices, Subpart F - Therapeutic Devices, Sec. 868.5895 Continuous ventilator - A continuous ventilator (respirator) is a device intended to mechanically control or assist patient breathing by delivering a predetermined percentage of oxygen in the breathing gas. Adult, pediatric, and neonatal ventilators are included in Class II.

10.2 Applicable Standards

- **IEC 60601-1:2005/A1:2012 (EN 60601-1:2006 + A1:2013)** - Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
- **ISO 80601-2-12:2011 / IEC 60601-2-12: 2001 (EN 60601-2-12:2006)** - Medical electrical equipment - Part 2-12: Particular requirements for the safety of lung ventilators - Critical care ventilators

- **ISO 5359:2008/Amd 1:2011 (EN ISO 5359:2008+A1:2011)** - Low-pressure hose assemblies for use with medical gases
- **NBR IEC 60601-1-2 (2017) / IEC 60601-1-2(2014) / (EN IEC 60601-1-2:2015)** - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
- **IEC 62304:2006 +AMD1:2015 (EN 62304:2006/2008)** - Medical device software - Software life-cycle processes
- **IEC 60601-1-8 Ed. 2.0 (2006)/A1:2012 (EN 60601-1-8:2007/A11:2017)** - Medical electrical equipment - Part 1-8: General requirements for basic safety and essential performance - Collateral Standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems
- **IEC 60601-1-4: 1996/A1:1999 (EN 60601-1-4: 1996/A1: 1999)** - Medical electrical equipment - Part 1-4: General requirements for safety - Collateral standard: Programmable electrical medical systems
- **IEC 60601-1-6: 2010 (EN 60601-1-6:2010)** - Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability
- **IEC 62366: 2007 (EN 62366:2008)** - Medical devices - Application of usability engineering to medical devices
- **EN ISO 17665-1:2006** - Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices
- **EN ISO 17664:2004** - Sterilization of medical devices - Information to be provided by the manufacturer for the processing of resterilizable medical devices
- **EN 1041:2008** - Information supplied by the manufacturer of medical devices
- **ISO 15223-1: 2016 (EN ISO 15223-1:2016)** - Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
- **ISO 80601-2-61:2011 (EN ISO 80601-2-61:2011)** – Medical electrical equipment: Particular requirements for basic safety and essential performance of pulse ox equipment
- **ISO 80601-2-55:2011 (EN 80601-2-55: 2011)** - Medical electrical equipment -- Part 2-55: Particular requirements for the basic safety and essential performance of respiratory gas monitors
- **ISO 780:2015** - Packaging -- Distribution packaging -- Graphical symbols for handling and storage of packages

10.3 Physical and Environmental Specifications

Table 10 - Physical and environmental specifications

Parameter	Specification	Tolerance	Unit
Dimensions and weight (basic unit)			
Height	1427	± 5	mm

Parameter	Specification	Tolerance	Unit
Dimensions and weight (basic unit)			
Width	453	± 5	mm
Depth (with handle)	544	± 5	mm
Weight	20	± 0.1	kg
Operation¹			
Temperature	-10 to 50	---	°C
Barometric Pressure	600 to 1100	---	cmH ₂ O
Relative air humidity (non-condensing)	15 to 95	---	%
Storage¹			
Temperature	-20 to 75	---	°C
Barometric Pressure	500 to 1200	---	cmH ₂ O
Relative air humidity (non-condensing)	5 to 95	---	%
Lifetime			
Fleximag Max 300, Fleximag Max 500, Fleximag Max 700	10	---	years

¹ Permissible operating and storage conditions for the whole electromedical system.

10.4 Electrical Specifications

10.4.1 Electrical network

Table 11 – Electrical network

Item	Specification	Tolerance	Unit
Electrical Network (50/60Hz)	100 to 240	± 10%	V _{AC}
Maximum power consumption	80	± 10%	VA

10.4.2 Batteries

Table 12 - Batteries

Item	Specification	Tolerance	Unit
Internal Batteries Li-Ion 11.8 VDC	4000	± 15%	mAh
Battery life (with full charge and normal use) ⁽¹⁾	210 minutes	± 15%	min
Average time to recharge to full (module in operation) ⁽¹⁾	4.0 hours	± 15%	h

(1) The battery charge should be made at the ambient temperature of 5 to 35 ° C

10.4.3 External DC Power Supply

Table 13 – External DC power supply

3.Item	Specification	Tolerance
Power supply ⁽¹⁾	Voltage: 12 to 15 V _{DC}	± 10%
	Current: 11.5 A	

(1) OPTIONAL external power supply



ATTENTION

- ***It is not possible to recharge the internal batteries of the equipment through the external DC power source.***
- ***The sole purpose of this input is to allow the equipment to be temporarily powered by a compatible external power source when there are no other alternatives.***
- ***In ventilation, before disconnecting an external power source, make sure that there is enough charge in the internal batteries or connect the equipment to the electrical network.***

10.4.4 Connectors

Table 14 - Connectors

Item	Specification
External power supply (network)	Connector 3 (three) pins, central pin ground. According to ABNT NBR 14136: 2012
External power supply	12-15V external input. Housing connector 3.96mm - 4 pins 180° female Green color
External sensor: capnograph ⁽¹⁾	Redel connector - 5 pin female receptacle Blue color
External sensor: oximeter ⁽¹⁾	Redel connector - 5 pin female receptacle Blue color
RS-232 standard connector (EIA RS-232C)	DB9 female type (on top) Used for maintenance and data transfer services through ARM (Assistance Remote Magnamed) only by authorized and trained Magnamed personnel.
Data output connector (network)	Standard RJ-45 Ethernet Connector Used to send data to an electronic health recorder. Use a CAT 5E in accordance with ANSI / TIA / EIA-568 or higher with a maximum length of 3 meters to connect the ventilator network port.
Video output connector	Standard connector HDMI
Standard USB connector	Used to transfer print screen, trends, logs and recordings to an external USB flash memory storage device ("flash drive"). It can also be used to software update, only authorized by Magnamed.
Nursing Call	Redel connector - 3 pin female receptacle Color-Red Only use cables that are checked by Magnamed.

(1) Optional



ATTENTION

- **Connected devices shall be approved medical devices that complies with IEC 60601-1.**
- **Use only certified cables on the equipment connectors.**
- **Connecting the ventilator to an IT network may result in risks to the patient, operator, or third party that has not been identified previously. The responsible organization shall identify, analyze, evaluate and control those risks.**
- **Subsequent changes to the IT network may introduce new risks and require further analysis by the responsible organization. Changes to the IT network include: configuration changes, connection of additional items, disconnection of items, updating of equipment connected to the IT network and improvement of equipment**

connected to the data communication port.

- *Failure to implement the communication protocol will result in failure to send data to other equipment.*
-

10.4.4.1 Protocol used for data communication with external devices

The Ethernet port may be used to share ventilator data, such as adjusted parameters, monitored parameters, waveforms, and alarm logging for electronic health recorders. The data has on average a delay of 8 seconds between the instant of data generation and the data output connector.

To send data to other equipment, the IT network must be scalable, with high availability and low data propagation delay.

The required network settings include a Dynamic Host Configuration Protocol (DHCP) enabled network server, so that the SEMP receives a valid Internet Protocol (IP). Communication is done through the TCP protocol in the IT network. For communication with the electronic health recorder, an appropriate communication protocol shall be implemented. To obtain the implementation guide for the communication protocol, contact MagnaService.

The information transitions from the following form: Fleximag Max sends the data to the responding client it received. The customer can ask questions or an operational request for Fleximag Max that immediately responds or confirms the request.



ATTENTION

- *To obtain the implementation guide for the communication protocol, contact MagnaService.*
 - *This implementation shall be performed in a network with the characteristics described in 10.4.4.1 by an IT specialist.*
 - *IT network failures to provide the required features may lead to delays in data communication or incorrect, incomplete or corrupted data transmission, resulting in incorrect information for the user.*
-

10.4.5 Gas Inlet Connector

Table 15 - Gas Inlet Connector

Item	Specification
Connections	According to ABNT NBR 11906: 2011
Hoses and extensions	According to ISO 5359: 2008
Compressed air and oxygen inlet pressure	200 to 600 kPa (29 to 87 psi)
Respiratory circuit	According to ISO 5367: 2014



ATTENTION

- *Input pressures exceeding the specified limit can damage the equipment.*
- *For inlet pressures less than 250 kPa, the maximum flow shall be 120 L / min.*




Observations

- *All materials that make up the product are compatible with oxygen gas, ambient air and medical compressed air.*

10.5 Pneumatic Specifications

10.5.1 Pneumatic layout¹



Legend

1	Filter	20	Proximal Flow Sensor
3	150 psi Pressure Sensor	21	Linear actuator
5	45 psi pressure regulated	22	Exp pressure measurement point
6	200 LPM Proportional	23	Distal flow sensor
7	Flow Air System ²	24	HEPA filter
8	One-way valve (anti-asphyxiation)	25	x-1 / 6 psi Solenoid valve
9	Solenoid Valve On/Off (NIF)	26	Measurement point (proximal pressure differential)
10	x-5 100 psi Solenoid Valve (Nebulizer)	27	Respiratory circuit
15	O ₂ concentration cell	28	HME filter
16	Internal Flow Sensor	29	Integrated valve with diaphragm
17	One-way valve – low pressure (insp output)	30	Humidifier (optional)
18	Insp pressure measurement point	31	Nebulizer (optional)
19	On/Off Overpressure Valve	32	Capnograph (optional)

Figure 1 – Pneumatic layout

¹ Pneumatic layout for Fleximag Max 700 model

² FlowAir system only for Fleximag 300 and Fleximag 700 models

10.6 Internal Flow Transducer Specifications

Table 16 - Internal Flow Transducer - General Specifications

General specifications	
<p>The internal flow transducer contains two sensors, one flow sensor and the other to measure the temperature.</p> <p>Each sensor has an independent nonlinear voltage output. To determine the mass flow of gas passing through the transducer, the output voltage of each of the sensors is measured.</p> <p>A microprocessor processes the results and calculates the flow using a specific algorithm.</p> <p>The circuit that measures the flow is usually known as a thermal sensor or hot wire anemometer.</p> <p>This flow transducer uses a heated wire sensor and maintained at a temperature of 150 ° C.</p> <p>The velocity of the gas passing through the sensor determines the rate of heat transfer between the sensor and the gas.</p> <p>This heat transfer rate is translated into a voltage necessary to maintain the temperature at 150 ° C.</p> <p>Consequently, this voltage is a function of the mass flow of gas passing through the sensor.</p> <p>The heat transfer rate is also influenced by the gas temperature.</p> <p>A thermistor circuit is used to measure the temperature of the gas and a correction is also made through specific algorithm.</p>	
Reading range	Air: 0 to 300 SLPM O ₂ : 0 to 300 SLPM
Specified Tolerance	Air: 2.0% or 0.05 SLPM (whichever is greater) O ₂ : 2.0% or 0.05 SLPM (whichever is greater)
Resistance	< 2.5mbar
Gas temperature range	5 to 46°C
Humidity Range	Dry gas (< 10% UR)
Operating Pressure	Atmospheric pressure
Power supply	5V ± 10% sensor and 2.7V - 5.5V Eeprom
Response Time	< 2.5ms
Burst pressure	Above 100 psi
Weight	21g

10.7 Distal Flow Sensor Specifications

Table 17 - Distal Flow Sensor - General Specifications

General specifications	
Intended use	Measure the flow expired by the patient.
Working principle	Thermal sensor or hot wire anemometer, thermistor.
Reading range	0 to 160 SLPM
Tolerance	$\pm 8\%$
Resistance	< 2.5mbar
Pressure range	$\pm 100\text{mbar}$
Lifetime	While calibration is successful
Material	MABS

10.8 Proximal Flow Sensor Specifications

Table 18 - Proximal Flow Sensor - General Specifications ADULT

General specifications	
Intended use	Measure the inspired and expired flow of the patient.
Working principle	Pressure differential
Reading range	-180 to 180 SLPM
Tolerance	$\pm 10\%$
Material	PSU

10.9 Pressure Sensor Specifications

Table 19 - Pressure sensor - General specifications

General specifications	
Intended use	Measure the inspired pressure of the patient
Working principle	Pressure differential
Reading range	-60 to 120 cmH ₂ O
Tolerance	$\pm 5\%$ (0 to 85°C)

Sensitivity	90mV/kPa
Response Time	< 1ms

10.10 Galvanic O₂ cell specifications

Table 20 – Galvanic O₂ cell specifications

General specifications	
Intended use	Measure the O ₂ concentration delivered from the equipment to the patient.
Measuring Range	0 to 100%
Exit signal	9 – 13 mV
Response time 90%	13 s
Accuracy	± 2%
Linearity	± 2%
Recommended flow rate	0.1 – 10 lpm
Data sampling rate	7 Hz
Method for calculating the gas level reading	Simple moving average (SMA) of 64 positions acquired every 140ms
Respiratory rate	The Respiratory Rate is displayed every 3 breaths and the mean value is updated with each breath.
Effects of gas and vapor interference	
Gas or steam	Gas level
Response to 80% NO	< 5%
Response to 7,5% Halothane	< 5%
Response to 7,5% Isoflurane	< 5%
Response to 7,5% Enflurane	< 5%
Response to 9% Sevoflurane	< 5%
Response to 20% Desflurane	< 5%
Response to 10% CO ₂	< 5%

10.11 Paramagnetic O₂ cell Specification


Table 21 – Paramagnetic O₂ cell specification

General specifications	
Intended use	Measure the O ₂ concentration delivered from the equipment to the patient.
Measuring Range	-15% to +200%
Exit signal	0 - 50 mV
Response time 90%	< 8 s
Accuracy	< ± 0.2%
Linearity	< ± 0.2%
Data sampling rate	7 Hz
Method for calculating the gas level reading	Simple moving average (SMA) of 64 positions acquired every 140ms
Respiratory rate	The Respiratory Rate is displayed every 3 breaths and the mean value is updated with each breath.
Effects of gas and vapor interference	
Gas or steam	Gas level
N ₂ O	- 0,20%
CO ₂	- 0,26%
H ₂ O	- 0,03%
Methane	- 0,16%
CO	0,06
Helium	0,29
NO	42,56
NO ₂	5,00

10.12 Ventilation Modes Specifications

10.12.1 VCV

VCV - Mandatory volume-controlled ventilation	
<p>Description:</p> <p>In this mode, the respiratory rate, tidal volume, and inspiratory flow (or inspiratory ratio or time) are set.</p> <p>The onset of inspiration (triggering) occurs according to the pre-established respiratory rate. triggering occurs exclusively for time, if the sensitivity setting is off.</p> <p>The transition between inspiration and expiration (cycling) occurs after release of the pre-set tidal volume at given velocity through the flow (or ratio or inspiratory time).</p>	<p>Adjusted Parameters:</p> <ul style="list-style-type: none"> • VOLUME • RATE • FLOW or RATIO or INS TIME • PEEP • O₂ CONCENTRATION • LIMIT PRESSURE • PAUSE (% or s) • PRESSURE TRIGGER • FLOW TRIGGER • FLOW WAVE <p>Note: Automatic backup ⁽¹⁾ 1 - Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.</p>



Once all ventilation parameters are received by the ventilator, it calculates the T_{INS}, T_{EXP}, T_{PAUSE} and I: E Ratio, as a function of the adjusted Flow, Pause, Wave and Frequency, thus obtaining all ventilation control times.

-
1. Ventilation without Inspiratory Pause, after the T_{INS} the ventilator cycles to expiration. The inspiratory pressure reached is a consequence of the volume delivered and the resistance and compliance of the patient's respiratory circuit.
 2. Ventilation with Inspiratory Pause, after delivery of the set volume the ventilator keeps the expiration interrupted until T_{INS} is completed, after which the ventilator cycles to expiration, the characteristic is the formation of pressure plateau (the difference between the peak and the plateau depends of airway resistance).
 3. If pressure or flow triggers are activated, then the ventilator attempts to synchronize the onset of the next inspiration with the patient's effort, according to the established levels. The information on which trigger type triggered the inspiratory cycle is reported in the status and message area. The patient's inspiratory effort detection for synchronization occurs at any point in the expiratory time.

Observations

- *If the patient makes inspiratory efforts and the triggers are properly adjusted, the ventilatory mode becomes assisted-controlled. In this situation, the monitored respiratory rate may be significantly greater than that adjusted.*

-
4. ASCENDENT (or accelerated) flow waveform.
 5. DECENDENT (or decelerated) flow waveform.
 6. SINUSOIDAL flow waveform.
 7. Representation of Pressure Limitation. In this situation the ventilator limits the pressure to the adjusted value because of factors such as patient compliance and imposed pressure limit, the adjusted volume is NOT DELIVERED, and this condition is informed in the alarms area (LIMITED PRESSURE message) on the screen.

WARNING

- ***Upon reaching the pressure limit set at the maximum pressure setting (LIMITED PRESSURE alarm), the set volume is NOT DELIVERED.***
- ***The default values are initial reference only.***
- ***Readjust ventilation parameters as required by the patient.***

10.12.2 PCV

PCV - Mandatory pressure-controlled ventilation

Description:

In this ventilatory mode, the respiratory rate, the inspiratory time, and the inspiratory pressure limit are set. The onset of inspiration (triggering) occurs according to the pre-established respiratory rate. The trigger, if the sensitivity adjustment is deactivated, is determined exclusively by the respiratory rate, and cycling (end of inspiration and beginning of expiration) is determined by the inspiratory time.


The tidal volume depends on the preset inspiratory pressure, the respiratory system impedance conditions and the inspiratory time selected by the operator.

Adjusted Parameters:

- CONTROL PRESSURE
- RATE
- INSPIRATORY TIME
- PEEP
- O₂ CONCENTRATION
- RISE TIME
- PRESSURE TRIGGER
- FLOW TRIGGER
- FLOW (\dot{V} - NEONATAL)

Note: Automatic backup ⁽¹⁾

1 - Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.



Once all ventilation parameters are set on the ventilator, it calculates the period, T_{EXP} and I: E ratio in function of T_{INS} and Rate, thus obtaining all ventilation control times.

1. Pressure Controlled Ventilation - The ventilator seeks to achieve adjusted inspiratory pressure in the shortest possible time, and this is accomplished by controlling the inspiratory flow.
2. Volume delivered to the patient is a consequence of the resistance and complacency of the patient's breathing circuit. The ventilator remains at the inspiratory pressure level set during T_{INS} after which the ventilator cycles to expiration, keeping the PEEP pressure set.
3. If pressure or flow triggering is activated, then the ventilator attempts to synchronize the onset of the next inspiration with the patient's effort, according to the established levels. The information

on which trigger type triggered the inspiratory cycle is reported in the status and message area of the screen. The patient's inspiratory effort detection for synchronization occurs at any point in the expiratory time.

Observation

- *If the patient makes inspiratory efforts and the sensitivities are properly adjusted, the ventilatory mode becomes assisted-controlled. In this situation, the monitored respiratory rate may be significantly greater than that adjusted.*

-
4. The rise time of the pressure can be adjusted by T_{RISE} (RISE TIME), the initial peak flow is generally lower than that in which $T_{\text{RISE}} = 0$ (depends on the resistance and compliance of the patient's respiratory circuit).

WARNING

- *The default values are initial reference only.*
- *Readjust ventilation parameters as required by the patient.*

10.12.3 PLV

PLV - Pressure limited ventilation

Description:

In this bias-flow ventilatory mode, the respiratory rate, inspiratory time, and inspiratory pressure limit are set. The onset of inspiration (triggering) occurs according to the pre-established respiratory rate. The trigger, if the sensitivity adjustment is disabled, is determined exclusively by the respiratory rate, but cycling (end of inspiration and beginning of expiration) occurs according to the inspiratory time.

The tidal volume depends on the preset inspiratory pressure, the respiratory system impedance conditions and the inspiratory time selected by the operator.


Usually by observing the flow curve, you see a peak of flow that decreases as time goes by.

Adjusted Parameters:

- INSPIRATORY PRESSURE
- RATE
- INSPIRATORY TIME
- PEEP
- O₂ CONCENTRATION
- FLOW (\dot{V})
- PRESSURE TRIGGER
- FLOW TRIGGER

Note: Automatic backup ⁽¹⁾

1 - Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.



Once all ventilation parameters are set in the ventilator, it calculates T_{EXP} as a function of Rate and T_{INS}, thus obtaining all ventilation control times.

1. Pressure Limited Ventilation - The ventilator seeks to achieve adjusted inspiratory pressure, and this is accomplished by occluding the integrated valve. It is important to note that the rise time of the pressure is dependent on the adjusted bias flow.

-
2. Volume delivered to the patient is a consequence of the resistance and complacency of the patient's breathing circuit. The ventilator remains at the inspiratory pressure level set during T_{INS} after which the ventilator cycles to expiration, keeping the PEEP pressure set.
 3. If pressure or flow triggering is activated, then the ventilator attempts to synchronize the onset of the next inspiration with the patient's effort, according to the established levels. The information on which trigger type triggered the inspiratory cycle is reported in the status and message area of the screen. The patient's inspiratory effort detection, for synchronization, occurs at any time during the expiratory time.



WARNING

- *The default values are initial reference only.*
- *Readjust ventilation parameters as required by the patient.*
- *The operator shall consider the inspiratory time and respiratory mechanics of the patient to define the bias flow adjusted. If the flow is not enough, the airway pressure may not reach the adjusted value.*



Observation

- *If the patient makes inspiratory efforts and the trigger are properly adjusted, the ventilatory mode becomes assisted-controlled. In this situation, the monitored respiratory rate may be significantly greater than that adjusted.*

10.12.4 PRVC

PRVC - Pressure regulated volume controlled

Description:

Time-limited, pressure-limited cycling mode that uses tidal volume as feedback to continuously adjust the pressure limit.

The first three respiratory cycles are in volume-controlled mode, allowing the ventilator to calculate respiratory mechanics. In the next cycles the ventilation is distributed with pressure limit and cycled by time to reach 60% of the adjusted volume. The onset of inspiration (triggering) occurs according to the pre-established respiratory rate. The end of the inspiration and beginning of the expiration (cycling) happens according to the adjusted inspiratory time.

At each cycle the ventilator adjusts the pressure limit (5 cmH₂O up) according to the tidal volume distributed in the previous cycle, until the tidal volume indicated by the operator is reached.

The maximum pressure limit is 5 cmH₂O below the pressure limit indicated by the operator.

Adjusted Parameters:

- VOLUME
- LIMIT PRESSURE
- RATE
- INSPIRATORY TIME
- PEEP
- O₂ CONCENTRATION
- RISE TIME
- PRESSURE TRIGGER
- FLOW TRIGGER

Note: Automatic backup ⁽¹⁾

1 - Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.



Once all ventilation parameters are set in the ventilator, it calculates the T_{EXP} as a function of Rate and T_{INS} , thus obtaining all ventilation control times.

1,2,3. Evaluation phase of lung compliance. After obtaining the compliance value the ventilator automatically sets a pressure value to reach 60% of the set volume, and then the ventilator adjusts the pressure automatically every three cycles of PCV.

4,5. Start of automatic pressure control to reach defined volume.

6. If the pressure trigger and flow trigger are active, then the ventilator seeks to synchronize the onset of the next inspiration to the patient's effort, according to the configured trigger. Detection of the patient's "window" of effort for synchronization begins in the last quarter of the period of controlled ventilation.

7. Volume reached.



WARNING

- *Upon reaching the pressure limit set at the maximum pressure setting (LIMITED PRESSURE alarm), the set volume is NOT DELIVERED.*
- *The default values are initial reference only.*
- *Readjust ventilation parameters as required by the patient.*




Observations

- *If the patient makes inspiratory efforts and the triggers are properly adjusted, the ventilatory mode becomes assisted-controlled. In this situation, the monitored respiratory rate may be significantly greater than that adjusted.*
- *Automatic pressure control occurs with $PEEP + 5\text{cmH}_2\text{O} \leq P_{Limit}$*

10.12.5 V-SIMV

V-SIMV - Volume controlled synchronized intermittent mandatory ventilation	
<p>Description:</p> <p>In V-SIMV, the respiratory rate, tidal volume and inspiratory flow or inspiratory time or ratio are set, as well as the trigger criterion for the occurrence of ventilator triggering by the patient.</p> <p>This mode allows the ventilator to apply predetermined preset cycles in sync with the patient's inspiratory effort.</p> <p>Mandatory cycles occur in the pre-determined window of time (the onset of inspiration occurs according to the pre-established respiratory rate), but synchronized with the patient's trigger.</p> <p>If there is apnea, the next cycle will be triggered by time until the patient's inspiratory incursions return.</p> <p>The transition between inspiration and expiration in the mandatory cycles occurs after release of the pre-set tidal volume at given velocity through the flow (or ratio or inspiratory time). In spontaneous cycles, cycling (transition between inhalation and expiration) occurs by flow as it falls to a certain percentage of its maximum reached value.</p> <p>In order to obtain the IMV in this mode, simply deactivate the support pressure by setting the support pressure equal to zero ($\Delta PS = 0$) or the flow and pressure sensitivities, equal to zero.</p>	<p>Adjusted Parameters:</p> <ul style="list-style-type: none">• VOLUME• RATE• FLOW or RATIO or INS TIME• PEEP• O₂ CONCENTRATION• FLOW WAVE FORM• PAUSE (%)• PRESSURE TRIGGER• FLOW TRIGGER• ΔPS (Support Pressure)• RISE TIME• FLOW CYCLING (% FLOW)• LIMIT PRESSURE <p>Note: Automatic backup ⁽¹⁾ 1 - Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.</p>



Once all ventilation parameters are set in the ventilator, it calculates the T_{INSP} and T_{EXP} in function of Flow, Pause, Wave and Rate, thus obtaining all ventilation control times.

1. It represents a cycle of VCV (volume controlled) with inspiratory pause;
2. Represents a spontaneous breathing cycle WITHOUT SUPPORT PRESSURE;
3. It represents a VCV (volume controlled) cycle after the SIMV period;
4. Represents a spontaneous breathing cycle WITH SUPPORT PRESSURE, whose cycling occurs by flow, as soon as it falls to a certain percentage of its maximum reached value.
5. The percentage of peak flow at which inspiratory phase cycles to the expiratory phase is programmable. Rise time (T_{RISE}) is also applied to support pressure (see PCV).
6. If the patient performs inspiratory effort, at the end of the SIMV (T_{SIMV}) period there is a window for synchronization of the controlled ventilation cycle, which is 'open' from $0.75 \times T_{SIMV}$, that is, in the last quarter of the period SIMV opens a synchronization window of the mandatory ventilation cycle. The information on which trigger type triggered the inspiratory cycle is reported in the message and status area of the screen.

WARNING

- **Adjusted support pressure (ΔPS) is a value above PEEP. Therefore, the inspiratory pressure of support will be the sum of PEEP and ΔPS .**
- **The default values are initial reference only.**


-
- *Readjust ventilation parameters as required by the patient.*
-

Observation

- *The monitored respiratory rate may be higher than the adjusted respiratory rate, as the patient may breathe spontaneously during the mandatory ventilation cycles.*
-

10.12.6 P-SIMV

P-SIMV – Pressure controlled synchronized intermittent mandatory ventilation	
<p>Description:</p> <p>In P-SIMV, the respiratory rate, inspiratory pressure and inspiratory time are set, as well as the sensitivity criterion for the occurrence of ventilator triggering by the patient.</p> <p>This mode allows the ventilator to apply predetermined preset cycles in sync with the patient's inspiratory effort.</p> <p>Mandatory cycles occur in the pre-determined window of time (the onset of inspiration occurs according to the pre-established respiratory rate) but synchronized with the patient's trigger.</p> <p>If there is an apnea, the next cycle will be triggered by time until the patient's inspiratory incursions return.</p> <p>The transition between inspiration and expiration in the mandatory cycles occurs after release of the pre-set tidal volume at given velocity through the flow (or ratio or inspiratory time). In spontaneous cycles, cycling (transition between inhalation and expiration) occurs by flow as it falls to a certain percentage of its maximum reached value.</p> <p>In order to obtain the IMV in this mode, simply deactivate the support pressure by setting the support pressure equal to zero ($\Delta PS = 0$) or the flow and pressure sensitivities, equal to zero.</p>	<p>Adjusted Parameters:</p> <ul style="list-style-type: none">• CONTROL PRESSURE• RATE• INSPIRATORY TIME• PEEP• O₂ CONCENTRATION• RISE TIME• ΔPS (Support Pressure)• FLOW CYCLING (% FLOW)• PRESSURE TRIGGER• FLOW TRIGGER• Flow (\dot{V} - NEONATAL) <p>Note: Automatic Backup ⁽¹⁾ 1 - Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.</p>



Once all the ventilation parameters are set in the ventilator, it calculates the T_{EXP} according to T_{INS} and Rate, thus obtaining all ventilation control times.

1. Represents a PCV (pressure controlled) cycle during T_{INS}
2. Represents a spontaneous breathing cycle WITHOUT SUPPORT PRESSURE;
3. Represents a PCV (pressure controlled) cycle after the SIMV period;
4. Represents a spontaneous breathing cycle WITH SUPPORT PRESSURE, whose cycling occurs by flow, as soon as it falls to a certain percentage of its maximum reached value.
5. The percentage of peak flow at which inspiratory phase cycles to the expiratory phase is programmable. Rise time (T_{RISE}) is also applied to support pressure (see PCV).
6. If the patient performs inspiratory effort, at the end of the SIMV (T_{SIMV}) period there is a window for synchronization of the controlled ventilation cycle, which is 'open' from $0.75 \times T_{SIMV}$, that is, in the last quarter of the period SIMV opens a synchronization window of the mandatory ventilation cycle. The information on which trigger type triggered the inspiratory cycle is reported in the message area and status of the screen.

WARNING

- **Adjusted support pressure (ΔPS) is a value above PEEP. Therefore, the inspiratory pressure of support will be the sum of PEEP and ΔPS .**

-
- *The default values are initial reference only.*
 - *Readjust ventilation parameters as required by the patient.*
-


Observation

- *The monitored respiratory rate may be higher than the adjusted respiratory rate, as the patient may breathe spontaneously during the mandatory ventilation cycles.*
-

10.12.7 CPAP/PS

CPAP / PS - Spontaneous ventilation with positive airway pressure and pressure support	
<p>Description:</p> <p>In CPAP / PS, the ventilator allows the patient to breathe spontaneously, but provides continuous pressurization in both inspiration and expiration, and assist ventilation during inspiration by maintaining a supporting pressure until the patient's inspiratory flow is reduced to a critical (adjustable) level of inspiratory peak flow reached.</p> <p>The onset of inspiration occurs when the patient makes an effort that is recognized by the ventilator according to the adjusted trigger. If the backup is active and the patient goes into apnea, the onset of inspiration will occur according to the set apnea time.</p> <p>This allows the patient to control the respiratory rate and inspiratory time, and thus the volume of inspiratory air.</p> <p>Cycle (transition between inspiration and expiration) occurs by flow as it falls to a certain percentage of its maximum reached value.</p> <p>If the support pressure value (ΔPS) is set to 0 (ZERO) and the means for cycle triggering are both turned off, ventilation with pure CPAP mode will be activated, which is a mode of spontaneous ventilation not assisted by the ventilator.</p> <p>The tidal volume depends on the respiratory effort of the patient and the conditions of respiratory mechanics of the lung and chest wall.</p>	<p>Adjusted Parameters:</p> <ul style="list-style-type: none"> • PEEP / CPAP • O₂ CONCENTRATION • ΔPS (Support Pressure) • CYCLE (% FLOW) • RISE TIME • PRESSURE TRIGGER • FLOW TRIGGER • FLOW (NEONATAL) • BACKUP <ul style="list-style-type: none"> ○ VCV Backup <ul style="list-style-type: none"> ▪ VOLUME ▪ RATE ▪ FLOW ▪ LIMIT PRESSURE ○ PCV Backup <ul style="list-style-type: none"> ▪ CONTROL PRESSURE ▪ RATE ▪ INSP TIME ○ PLV-NEONATAL Backup <ul style="list-style-type: none"> ▪ INSP PRESSURE ▪ RATE ▪ INSP TIME ○ Automatic Backup⁽¹⁾ 1 - Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.

CPAP / PS - Spontaneous ventilation with positive airway pressure and pressure support



1 and 2. represent spontaneous cycles with the support pressure in ZERO.

3, 4 and 5. They represent spontaneous breathing cycles of the patient with support pressure greater than zero. T_{RISE} (Rise Time) of the support pressure can be adjusted so that the initial flow is smoothed.


The cycling occurs by flow as soon as it falls to a certain percentage of its maximum reached value.

6. If the patient goes into apnea, after T_{APNEIA} (s) the ventilator will present this condition through an alarm in its message area on the screen and will start the selected backup ventilation, according to the settings and programmed parameters.



WARNING

- **Adjusted support pressure (ΔPS) is a value above PEEP. Therefore, the inspiratory pressure of support will be the sum of PEEP and ΔPS .**
- **The default values are initial reference only.**
- **Readjust the ventilation parameters according to the patient's need**





1 and 2. They represent spontaneous cycles.

1. If the patient does not breathe after the time for apnea, the ventilator enters the backup and activates the apnea alarm.

10.12.8 DualPAP

DualPAP - Two-stage positive pressure ventilation			
<p>Description:</p> <p>In DualPAP, the ventilator operates at two pressure levels set by the operator, Superior Pr and Inferior Pr:</p> <p>The change to the lower pressure level (end of inspiration) occurs at the end of Superior T (time set for the superior pressure level). Likewise, restoring the superior pressure level (beginning of inspiration) occurs as soon as the Lower T (the time for the lower pressure level) is exhausted.</p> <p>Consequently, the respiratory rate and the I: E ratio are directly related to this alternation between levels.</p> <p>DualPAP allows spontaneous cycles at both pressure levels and counts with the possibility of synchronization with the patient's inspiratory effort. Cycle (transition between inspiration and expiration) occurs by flow as it falls to a certain percentage of its maximum reached value.</p> <p>Without spontaneous breathing, DualPAP is similar to the controlled pressure mode, differing from this by adjusting the times (Superior T and Inferior T) rather than respiratory rate.</p> <p>The onset of inspiration occurs when the patient makes an effort that is recognized by the ventilator according to the adjusted trigger. If the backup is active and the patient goes into apnea, the onset of inspiration will occur according to the set apnea time.</p>	<p>Adjusted Parameters:</p> <ul style="list-style-type: none"> • SUPERIOR PR. • SUPERIOR T. • INFERIOR PR. • INFERIOR T. • O₂ CONCENTRATION • ΔPS (Support Pressure) • LIMIT PRESSURE • CYCLE (% FLOW) • PRESSURE TRIGGER • FLOW TRIGGER • RISE TIME • FLOW (\dot{V} - NEONATAL); • BACKUP <table style="width: 100%; border: none;"> <tr> <td style="vertical-align: top;"> <ul style="list-style-type: none"> ○ VCV Backup <ul style="list-style-type: none"> ▪ VOLUME ▪ RATE ▪ FLOW ▪ PEEP ▪ FLOW WAVE ▪ PAUSE ○ PCV Backup <ul style="list-style-type: none"> ▪ CONTROL PRESSURE ▪ RATE ▪ INSPIRATORY TIME ▪ PEEP </td> <td style="vertical-align: top;"> <ul style="list-style-type: none"> ○ PLV-NEONATAL Backup <ul style="list-style-type: none"> ▪ INSP PRESSURE ▪ RATE ▪ INSPIRATORY TIME ○ Automatic Backup⁽¹⁾ <p>1 – Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.</p> </td> </tr> </table>	<ul style="list-style-type: none"> ○ VCV Backup <ul style="list-style-type: none"> ▪ VOLUME ▪ RATE ▪ FLOW ▪ PEEP ▪ FLOW WAVE ▪ PAUSE ○ PCV Backup <ul style="list-style-type: none"> ▪ CONTROL PRESSURE ▪ RATE ▪ INSPIRATORY TIME ▪ PEEP 	<ul style="list-style-type: none"> ○ PLV-NEONATAL Backup <ul style="list-style-type: none"> ▪ INSP PRESSURE ▪ RATE ▪ INSPIRATORY TIME ○ Automatic Backup⁽¹⁾ <p>1 – Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.</p>
<ul style="list-style-type: none"> ○ VCV Backup <ul style="list-style-type: none"> ▪ VOLUME ▪ RATE ▪ FLOW ▪ PEEP ▪ FLOW WAVE ▪ PAUSE ○ PCV Backup <ul style="list-style-type: none"> ▪ CONTROL PRESSURE ▪ RATE ▪ INSPIRATORY TIME ▪ PEEP 	<ul style="list-style-type: none"> ○ PLV-NEONATAL Backup <ul style="list-style-type: none"> ▪ INSP PRESSURE ▪ RATE ▪ INSPIRATORY TIME ○ Automatic Backup⁽¹⁾ <p>1 – Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.</p> 		

DualPAP - Two-stage positive pressure ventilation





WARNING

- *The support pressure (ΔPS) is a value above the superior or inferior pressure. Therefore, the maximum support pressure will be the sum of that reference pressure and ΔPS .*
- *The default values are initial reference only.*
- *Readjust ventilation parameters as required by the patient.*
- *In the absence of spontaneous cycles when in DualPAP, adjust the superior and inferior pressures so that the minute volume delivered to the patient is enough.*



Observation


- *Changes in pressure levels are synchronized.*

10.12.9 APRV

APRV - Continuous positive pressure ventilation and airway pressure relief

<p>Description:</p> <p>In APRV, the ventilator operates at two pressure levels set by the operator, Superior Pr and Inferior Pr.</p> <p>The transient relief for the lower pressure level (end of inspiration) occurs at the end of the superior T (time set for the superior pressure level). Likewise, restoring the superior pressure level (beginning of inspiration) occurs as soon as the inferior T (pressure relief time) is exhausted.</p> <p>Consequently, the respiratory rate and the resulting I: E ratio are directly related to this alternation between levels.</p> <p>APRV has the characteristic of inversion of the I: E ratio, where the time of the inferior pressure level is usually less than that of the superior pressure level, functioning only as a temporary relief. Cycle (transition between inspiration and expiration) occurs by flow as it falls to a certain percentage of its maximum reached value.</p> <p>Without spontaneous breathing, APRV is similar to the controlled pressure mode, differing from this by adjusting the times (superior and inferior) rather than the respiratory rate.</p> <p>The onset of inspiration occurs when the patient makes an effort that is recognized by the ventilator according to the adjusted trigger. If the backup is active and the patient goes into apnea, the onset of inspiration will occur according to the set apnea time.</p>	<p>Adjusted Parameters:</p> <ul style="list-style-type: none"> • SUPERIOR PR. • SUPERIOR T. • INFERIOR PR. • INFERIOR T. • O₂ CONCENTRATION • ΔPS (Support Pressure) • LIMIT PRESSURE • CYCLE (% FLOW) • PRESSURE TRIGGER • FLOW TRIGGER • RISE TIME • FLOW (\dot{V} - NEONATAL); • BACKUP • ○ VCV Backup <ul style="list-style-type: none"> ▪ VOLUME ▪ RATE ▪ FLOW ▪ PEEP ▪ FLOW WAVE ▪ PAUSE ○ PCV Backup <ul style="list-style-type: none"> ▪ CONTROL PRESSURE ▪ RATE ▪ INSPIRATORY TIME ▪ PEEP ○ PLV-NEONATAL Backup <ul style="list-style-type: none"> ▪ INSP PRESSURE ▪ RATE ▪ INSPIRATORY TIME ○ Automatic Backup⁽¹⁾ <p>1 – Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.</p>
--	---

APRV - Continuous positive pressure ventilation and airway pressure relief



WARNING

- *The support pressure (ΔPS) is a value above the superior or inferior pressure. Therefore, the maximum support pressure will be the sum of that reference pressure and ΔPS .*
- *The default values are initial reference only.*
- *Readjust ventilation parameters as required by the patient.*

10.12.10 MMV

MMV – Spontaneous ventilation with mandatory minute volume	
<p>Description:</p> <p>Initially, the ventilator allows spontaneous test cycles with a support pressure of 5 or 10 cmH₂O above the adjusted PEEP.</p> <p>Then the minute volume is measured, and approximate compliance is calculated.</p> <p>For each subsequent cycle, the ventilator recalculates the compliance of the previous cycle and adjusts the pressure level to the next cycles, reaching the adjusted minute volume.</p> <p>The increase in pressure between the cycles never exceeds 3 cmH₂O and the maximum level reached does not exceed the value of the set limit pressure.</p> <p>If this value is reached, and the adjusted minute volume is not reached, the limited pressure alarm will be displayed.</p> <p>On the other hand, if the set apnea time is reached, the backup ventilation will be activated.</p> <p>The onset of inspiration occurs when the patient makes an effort that is recognized by the ventilator according to the adjusted trigger. If the backup is active and the patient goes into apnea, the onset of inspiration will occur according to the set apnea time. Cycle (transition between inspiration and expiration) occurs by flow as it falls to a certain percentage of its maximum reached value.</p>	<p>Adjusted Parameters:</p> <ul style="list-style-type: none">• MINUTE VOLUME• PEEP• CYCLE (% FLOW)• PRESSURE TRIGGER• FLOW TRIGGER• O₂ CONCENTRATION• RISE TIME• LIMIT PRESSURE• BACKUP <ul style="list-style-type: none">○ VCV Backup<ul style="list-style-type: none">▪ VOLUME▪ RATE▪ FLOW▪ FLOW WAVE▪ PAUSE○ PCV Backup<ul style="list-style-type: none">▪ CONTROL PRESSURE▪ RATE▪ INSPIRATORY TIME○ Automatic Backup⁽¹⁾ 1 – Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.

10.12.11 VS

VS - Spontaneous ventilation with assured volume	
<p>Description:</p> <p>Initially, the ventilator allows spontaneous test cycles with a support pressure of 5 or 10 cmH₂O above the adjusted PEEP.</p> <p>Then the delivered volume is measured, and the approximate compliance is calculated.</p> <p>For each subsequent cycle, the ventilator recalculates the compliance of the previous cycle and adjusts the pressure level for the next cycles, reaching the adjusted volume.</p> <p>The increase in pressure between the cycles never exceeds 3 cmH₂O and the maximum level reached does not exceed the value of the set limit pressure.</p> <p>If this value is reached, without the adjusted volume being reached, the limited pressure alarm will be displayed.</p> <p>On the other hand, if the set apnea time is reached, the backup ventilation will be activated.</p> <p>The onset of inspiration occurs when the patient makes an effort that is recognized by the ventilator according to the adjusted trigger. If the backup is active and the patient goes into apnea, the onset of inspiration will occur according to the set apnea time. Cycle (transition between inspiration and expiration) occurs by flow as it falls to a certain percentage of its maximum reached value.</p>	<p>Adjusted Parameters:</p> <ul style="list-style-type: none">• VOLUME• PEEP• CYCLE (% FLOW)• PRESSURE TRIGGER• FLOW TRIGGER• O₂ CONCENTRATION• RISE TIME• LIMIT PRESSURE• BACKUP <ul style="list-style-type: none">○ PRVC Backup<ul style="list-style-type: none">▪ RATE▪ INSPIRATORY TIME○ Automatic Backup⁽¹⁾ 1 – Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.



WARNING

- ***The apnea time alarm can be turned off. In this condition, BACKUP VENTILATION WILL NEVER BE ACTIVATED.***
- ***If backup ventilation is turned off, the operator must be sure that this adjustment is necessary and aware of the clinical implications involved.***

10.12.12 Nasal CPAP

Nasal CPAP - Spontaneous ventilation with positive nasal airway pressure	
<p>Description:</p> <p>As in CPAP mode, in Nasal CPAP the ventilator allows the patient to breathe spontaneously, with the difference that the ventilator compensates for leaks automatically and ignores high minute volume, high current tidal and flow sensor verification.</p> <p>This mode is only available for neonatal patients and requires a compatible nasal interface.</p> <p>The onset of inspiration occurs when the patient makes an effort that is recognized by the ventilator according to the adjusted trigger. If the backup is active and the patient goes into apnea, the onset of inspiration will occur according to the set apnea time. Cycle (transition between inspiration and expiration) occurs by flow as it falls to a certain percentage of its maximum reached value.</p>	<p>Adjusted Parameters:</p> <ul style="list-style-type: none"> • PEEP / CPAP • O₂ CONCENTRATION • ΔPS (Support Pressure) • CYCLE (% FLOW) • RISE TIME • PRESSURE TRIGGER • FLOW • BACKUP <ul style="list-style-type: none"> ○ PLV-NEONATAL Backup <ul style="list-style-type: none"> ▪ INSP PRESSURE ▪ RATE ▪ INS TIME ○ Automatic Backup⁽¹⁾ <p>1 - Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.</p>

10.12.13 NIV

NIV - Non-Invasive Ventilation	
<p>Description:</p> <p>As in CPAP mode, in NIV the ventilator allows the patient to breathe spontaneously, with the difference that the ventilator compensates for leaks automatically</p>	<p>Adjusted Parameters:</p> <ul style="list-style-type: none"> • PEEP / CPAP • O₂ CONCENTRATION • ΔPS (Support Pressure) • CYCLE (% FLOW) • RISE TIME • PRESSURE TRIGGER • FLOW (Neonatal) • BACKUP

NIV - Non-Invasive Ventilation

<p>and ignores the high minute volume, high current tidal and flow sensor verification.</p> <p>This mode is available to all patients.</p> <p>The onset of inspiration occurs when the patient makes an effort that is recognized by the ventilator according to the adjusted trigger. If the backup is active and the patient goes into apnea, the onset of inspiration will occur according to the set apnea time. Cycle (transition between inspiration and expiration) occurs by flow as it falls to a certain percentage of its maximum reached value.</p>	<ul style="list-style-type: none"> ○ VCV Backup <ul style="list-style-type: none"> ▪ VOLUME ▪ RATE ▪ FLOW ▪ LIMIT PRESSURE ▪ FLOW WAVE ▪ PAUSE ○ PCV Backup <ul style="list-style-type: none"> ▪ CONTROL PRESSURE ▪ RATE ▪ INSP TIME ○ PLV-NEONATAL Backup <ul style="list-style-type: none"> ▪ INSP PRESSURE ▪ RATE ▪ INS TIME ○ Automatic Backup⁽¹⁾ <p>1 - Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.</p>
---	---

10.12.14 VG

VG - Guaranteed Volume Ventilation

<p>Description:</p> <p>Available only for neonatal patient, this ventilatory mode is similar to the PLV, except that a tidal volume value can be set to be delivered by the ventilator.</p> <p>The onset of inspiration (triggering) occurs according to the pre-established respiratory rate. Cycle (transition between inspiration and expiration) occurs by flow as it falls to a certain percentage of its maximum reached value.</p>	<p>Adjusted Parameters:</p> <ul style="list-style-type: none"> • INSPIRATORY PRESSURE • RATE • INSPIRATIONAL TIME • PEEP • O₂ CONCENTRATION • FLOW (\dot{V}) • PRESSURE TRIGGER • FLOOD TRIGGER • VOLUME <p>Obs.: Automatic Backup⁽¹⁾ 1 - Whenever the set apnea time is reached, the ventilator triggers a ventilator cycle, the setting of which is based on the current mode settings.</p>
--	---



WARNING

- *This mode is only available in neonatal patient with proximal flow sensor.*
- *It is mandatory to perform auto test to enable this mode.*

10.12.15 O₂ THERAPY

O ₂ THERAPY	
<p>Description:</p> <p>Available for all patient types, this special mode offers continuous O₂ flush according the O₂ concentration and flow set. This mode requires a compatible interface.</p>	<p>Adjusted Parameters:</p> <ul style="list-style-type: none">• FLOW• O₂ CONCENTRATION

10.13 Oxygen adjustment response time

The table below indicates the time required for oxygen concentration in the delivered volume to change from a fraction of 21% to 90% of the maximum oxygen concentration adjustable using the maximum internal volume respiratory system configuration.

Table 22 – Oxygen adjustment response time

Delivered volume (mL)	Maximum response time (s)
500	98
150	83
30	360

10.14 Accuracy of settings

The following table shows the maximum error between the adjusted value and the value applied by the ventilator.

Table 23 – Accuracy of settings

Item	Parameter	Accuracy ^{1,2}
1	Tidal volume (All range) ¹	$\pm (4 \text{ mL} + 15\% \text{ of set volume})$
2	Inspiratory pressure	$\pm (2 \text{ cmH}_2\text{O} + 4\% \text{ of set pressure})$
3	PEEP	$\pm (2 \text{ cmH}_2\text{O} + 4\% \text{ of set PEEP})$
4	FiO ₂	$\pm (5\% + 2,5\% \text{ of set FiO}_2)$

1 The volume and pressure accuracy are preserved for respiratory circuits with resistance up to 1.9 cmH₂O and compliance up to 5 mL/cmH₂O

2 Performance accuracies were determined in the worst-case configuration, with respiratory circuit with water trap, heated humidifier, proximal flow sensor and CO₂ sensor, because they have greater compliance and a greater number of potential leakage points.

The performance accuracies were determined through a test system with the measurement uncertainties described in the table below:

Table 24 – Uncertainty of settings

Item	Parameter	Uncertainty
1	Tidal volume	± 2,5%
2	Inspiratory pressure	± 2,0%
3	PEEP	± 2,0%
4	FiO ₂	± 2,0%

10.15 Adjustable Parameter Specifications

Table 25 - Adjustable parameters

Item	Parameter	Specification	Resolution		Unit
1	Tidal volume	2.0 to 3000	Neonatal ⁽¹⁾	2.0 to 10.0: 0.1	mL
				10 to 99: 1	
			Pediatric	10 to 100: 5	
				100 to 300: 10	
			Adult	100 to 1000: 10	
				1000 to 3000: 50	
2	Respiratory rate ^{(2) (3)}	0 to 200	Neonatal	0 to 200: 1	rpm
			Pediatric	0 to 200: 1	
			Adult	0 to 100: 1	
3	Rise time	0 to 2.0	0.1		s
4	Pause (plateau)	0 to 70	10		%
		0 to 2	0.1		s
5	Inspiratory pressure and limit	0 to 120	1		cmH ₂ O
6	ΔPS	0 to 120	1		cmH ₂ O
7	PEEP	0 to 50	1		cmH ₂ O
8	Pressure trigger	0.0 to -20	0.0 to -2.0: - 0.2		cmH ₂ O
			-2 to -10: - 1		
9	Flow trigger	0.0 to 30	0.0 to 1.0: 0.1		L/min

Item	Parameter	Specification	Resolution	Unit
			1.0 to 30.0: 0.5	
10	Cycle	5 to 80 (maximum 3 s)	5	%
11	O ₂ Concentration	21 to 100	1	% vol
12	Inspiratory time	0.05 to 30	0.05 to 0.70: 0.01	s
			0.70 to 1.00: 0.05	
			1.0 to 30.0: 0.1	
13	Flow waveform	Square, Descending or Decelerating, Ascending or Accelerated, Sinusoidal	---	---
14	CPAP	0 to 50	1	cmH ₂ O
15	Superior Pressure (DualPAP/APRV)	5 to 90	1	cmH ₂ O
16	Inferior Pressure (DualPAP/APRV)	0 to 45	1	cmH ₂ O
17	Superior Time (DualPAP/APRV)	0.10 to 59.8	0.20 to 0.70: 0.01	s
			0.70 to 1.00: 0.05	
			1.00 to 59.80: 0.10	
18	Inferior Time (DualPAP/APRV)	0.20 to 59.9	0.20 to 0.70: 0.01	s
			0.70 to 1.00: 0.05	
			1.00 to 59.90: 0.10	
19	Ratio I: E	1:599 to 299:1 ⁽³⁾	1:0.1	---
20	Backup ⁽⁴⁾	OFF, PLV, PCV, VCV and PRVC	---	---
21	Inspiratory flow	1 to 180	1	L/min
22	Patient Height	64 to 132 (Pediatric)	1	cm
		133 to 250 (Adult)	1	cm
23	Patient weight	0,1 to 5,9 (Neonatal)	0,1	kg
24	Nebulizer flow - 100% oxygen ⁽⁵⁾	5 to 8 (without direct adjustment)	---	L/min
25	TGI (Tracheal Gas Insufflation) Flow - 100% Oxygen ⁽⁵⁾	5 to 8 (without direct adjustment)	---	L/min
26	Sigh ⁽⁶⁾	1 to 3	1	gasp

Item	Parameter	Specification	Resolution	Unit
27	Volume of the sigh ⁽⁶⁾	10 to 100	10	% Vt
28	Rate of the sigh ⁽⁶⁾	20 to 100	10	cycles
29	Tube compensation	Endotracheal Tracheostomy	---	---
30	Diameter of the tube	2.5 to 12.0	2.5 to 10 : 0.5	mm
			10 to 12 : 1	
31	% Tube Compensation	10 to 100	10	%
32	Minimum inspiratory hold ⁽⁷⁾	0,1 to 30	0,1 to 1 : 0.1	s
	Minimum expiratory hold		1 to 30 : 1	
33	Alarm muting time	OFF, 10 to 120	10	s
34	Screen Lock Time	OFF, 1 to 30	1 to 5: 1	min
			10 to 30: 5	
35	Flow (only in O ₂ Therapy mode)	0 to 60	0 to 60: 1	L/min

(1) Neonatal patient volume only with NEONATAL proximal flow sensor and in VG mode.

(2) Respiratory rate 0 (zero) will only be reached in spontaneous modes, with triggers and backup off.

(3) The minimum and maximum values of rate and I: E ratio depend on the adjusted ventilatory mode.

(4) Adjustable backup options for spontaneous modes, for other modes, backup is automatic.

(5) Nebulizer and TGI flows cannot be activated simultaneously.

(6) Sigh adjustment available only in VCV and V-SIMV modes

(7) Hold time when immediately pressing and releasing button.



ATTENTION

- ***The FlexiMag Max ventilator cares for any patient, from premature to morbidly obese, however, the patient height adjustment used for optimal weight calculation is limited.***
- ***For patients who exceed this limit, the parameters can be adjusted directly by the operator.***

10.16 Monitored Parameter Specifications

The monitored parameters are calculated using the ordinary least squares (OLS) technique, which calculates the mean through the most recent values in a data series. Thus, for each value included in the averaging calculation, the oldest value is excluded.

Table 26 - Monitored ventilatory parameters

Parameter	Range	Resolution	Measurement Accuracy ⁽¹⁾
Instantaneous pressure	0 to 120 cmH ₂ O	1	± (2 cmH ₂ O + 4% of reading)
Peak pressure (P _{peak})	0 to 120 cmH ₂ O	1	± (2 cmH ₂ O + 4% of reading)
Mean pressure (P _{mean})	0 to 120 cmH ₂ O	1	± (2 cmH ₂ O + 4% of reading)
Plateau Pressure (P _{plat})	0 to 120 cmH ₂ O	1	± (2 cmH ₂ O + 4% of reading)
PEEP	0 to 120 cmH ₂ O	1	± (2 cmH ₂ O + 4% of reading)
Intrinsic PEEP (iPEEP)	0 to 99.9 cmH ₂ O	1	± (2 cmH ₂ O + 4% of reading)
Measured Flow ⁽⁴⁾	-180 to 180 L/min	1	± (50mL/min + 10 % of reading)
Tidal Volume in Volume Guaranteed (V _{te}) ^{(2) (3) (4)}	0.0 to 10.0 mL	0.1	±(4,0 mL + 15% of reading)
	10 to 100 mL	1	
Tidal volume (V _{te}) ^{(3) (4)}	0 to 999 mL	1	±(4,0 mL + 15% of reading)
	1,00 to 3,00 L	0.01	
Minute volume ⁽⁴⁾ (MV)	0.001 to 0.999 L	0.001	±(4,0 mL + 15 of reading)
	1.00 to 9.99 L	0.01	
	10.0 to 99.9 L	0.1	
Inspiratory time	0.05 to 9.99 s	0.01	± (0,10 s + 10% of reading)
	10.0 to 60.0 s	0.1	± (0,10 s + 10% of reading)
Expiratory time	0.05 to 9.99 s	0.01	± (0,10 s + 10% of reading)
	10.0 to 60.0 s	0.1	± (0,10 s + 10 of reading)
I:E Ratio	1:599 to 599:1	1:0.1	± (0,1 + 10 % of reading)
Total respiratory rate	0 to 200 bpm	1	± (1bpm + 10% of reading)

Parameter	Range	Resolution	Measurement Accuracy ⁽¹⁾
Spontaneous respiratory rate	0 to 200 bpm	1	± (1bpm + 10% of reading)
O ₂ Concentration (O ₂)	12,0 to 99.9 %	0.1	± (2,5% + 2,5% of reading)
	100 to 110 %	1	
Airway Resistance ⁽³⁾ (Ri and Re)	0 to 99.9 cmH ₂ O/L/s	0.1	± (5cmH ₂ O/L/s +20% of reading)
	100 to 200 cmH ₂ O/L/s	1	
Dynamic Compliance	0 to 99.9 mL/cmH ₂ O	0.1	± (1mL/cmH ₂ O + 10% of reading)
	100 to 200 mL/cmH ₂ O	1	
Static Compliance	0 to 99.9 mL/cmH ₂ O	0.1	± (1mL/cmH ₂ O + 10% of reading)
	100 to 200 mL/cmH ₂ O	1	
Auxiliary pressure	0 to 120 cmH ₂ O	1	± (2 cmH ₂ O + 4% of reading)
Estimated tracheal pressure	0 to 120 cmH ₂ O	1	± (2 cmH ₂ O + 4% of reading)
Elastance	0 to 100 cmH ₂ O/L	1	± (1cmH ₂ O/L + 10% of reading)
Leakage flow ⁽⁴⁾	0.0 to 19.9 L/min	0.1	± (50mL/min + 10% of reading)
	20 to 180 L/min	1	
Percentage leakage ⁽⁴⁾	0 to 100 L/min	1	± (10% + 10% of reading)
Time constant (TC)	Calculated (s)	0.1	± (0,1s + 10% of reading)
Ti / Ttotal	Calculated (s)	0.1	± (0,1s + 10% of reading)
RSBi - Surface Respiration Index (IRRS, Tobin Index)	Calculated (cycles/min/L)	1	± (1 cycle/min/L + 10% of reading)
WOBi (Respiratory work imposed)	Calculated (J/min)	0.01	± (0,1 J/min + 10% of reading)
WOBi (Respiratory work imposed)	Calculated (J/L)	0.01	± (0,1 J/L + 10% of reading)
Pi Max	-60 to 120 cmH ₂ O	1	± (2 cmH ₂ O + 4% of reading)
Driving Pressure	Calculated (Pplateau – PEEP)	1	± (2 cmH ₂ O + 4% of reading)
O ₂ consumption	0 to 180 L/min	0.1	± (50mL/min + 10% of reading)

-
- (1) When two tolerances are indicated, consider the one with the highest value. The accuracy of volume and pressure is preserved for respiratory circuits with resistance up to 1.9 cmH₂O and compliance up to 5 mL / cmH₂O.
 - (2) Only in VG mode.
 - (3) For airway resistance greater than 150 cmH₂O / L / s, the expired volume will have its tolerance outside the declared. In this condition, the measured inspiratory volume remains unchanged.
 - (4) Volumes, flows and leakage specifications associated with the ventilator respiratory system are expressed in BTPS.



WARNING

- ***Ventilation with cyclic pressure up to 100 cmH₂O can add tolerance error of up to 2%.***
- ***The accuracy of the parameters can be affected under the following conditions:***
 - ***Reuse of single-use accessories;***
 - ***Incorrect sensor uses for patient type;***
 - ***Secretion in the circuit and flow sensor;***
 - ***Leakage in the circuit;***
 - ***Oxygen cell and flow sensor not calibrated;***
 - ***Condensation in the circuit, flow sensor and gas inlet;***
 - ***Use of nebulizer with the capnography sensor;***
 - ***Proximal flow sensor with tubes facing down.***
- ***To maintain the accuracy, keep the flow sensor, breathing circuit, and gas inlet dry, clean, and non-condensing.***



Observations

- *In practice, the units of measure of pressure are equivalent, being able to adopt that 1 mbar = 1 hPa ≈ 1 cmH₂O.*
- *Oxygen monitoring is obtained accurately in up to 20 seconds after initialization.*

10.17 Safety System Specifications

- Anti-asphyxia valve for protection against gas supply failure.
- Safety relief valve of 100 cmH₂O, according to basic ventilator norm in order to avoid overpressure in the respiratory circuit.

- Active overpressure valve that, when detecting obstructions, is activated to reduce the pressure in the respiratory circuit.
- High Pressure Relief Valve is activated when the mains pressure is above 800 kPa (8.0 kg / cm²) by sending the excess gas into the atmosphere. This will limit the ventilator supply pressure.

10.18 Alarm System Specification

The alarm priority is determined by the device's risk management process.

Table 27 - Alarm Priority

Potential result of failure to respond to alarm cause	Potential Damage Onset ⁽¹⁾		
	Immediate ⁽²⁾	Promptly ⁽³⁾	Late ⁽⁴⁾
Death or Irreparable injury	HIGH PRIORITY	HIGH PRIORITY	MEDIUM PRIORITY
Repairable injury	HIGH PRIORITY	MEDIUM PRIORITY	-
Mild injury or discomfort	MEDIUM PRIORITY	-	-

- (1) Initiation of potential damage refers to the occurrence of the lesion and not to its manifestation
 (2) There is potential for the event to develop over a period usually not enough for manual corrective action.
 (3) There is potential for the event to develop over a period usually enough for manual corrective action.
 (4) There is potential for the event to unfold in an unspecified period greater than that provided at the "prompt".

In this alarm system there is no change in the priority of the alarm condition and in the occurrence of more than one alarm simultaneously:

- High priority alarm messages are displayed alternately.
- In the absence of high priority alarms, the medium priority alarm messages will be displayed alternately

Alarm messages are displayed as soon as the alarm condition is detected, so there is no delay in displaying the messages.

Table 28 - Alarms Characteristics

Alarm	Characteristic	High priority	Medium Priority
Visual	Color	Red	Yellow
	Intermittent frequency	1,6 hz	0,7 hz

Sound	Number of saved pulses	10 pulses	3 pulses
	Interval between saved	5s	5s
	Sound pressure range	61dBA	60dBA
	Pulse frequency	686 hz	686 hz

Observations

- To identify the occurrence of an alarm, the operator should preferably be 1 meter from the front of the equipment at an angle of 30 ° with the horizontal axis at the center of the monitor's viewing plane.

10.18.1 Adjustable alarm specifications

Table 29 - Adjustable alarms

Alarm	Adjustment	Limit	Default Values ⁽¹⁾			Unit
			Neonatal	Pediatric	Adult	
Maximum Pressure	OFF, 0 to 120	High	30	30	40	cmH2O
		Low	OFF	OFF	OFF	
PEEP	OFF, 0 to 80	High	10	15	20	cmH2O
		Low	OFF	OFF	OFF	
Total Volume	OFF, 0 to 3000	High	50 mL	500 mL	1.0 L	L or mL
		Low	OFF	OFF	OFF	
Minute volume	OFF, 0.0 to 99	High	1.0	10	20	L
		Low	0.5	2	3.6	
Respiratory rate	OFF, 0 to 200	High	80	60	60	rpm
		Low	5	5	5	
O ₂	OFF, 18 to 100	High	80	80	80	%
		Low	OFF	OFF	OFF	
EtCO ₂ ⁽²⁾	OFF, 0 to 80	High	45	45	45	mmHg
		Low	OFF	OFF	OFF	
CO ₂ Ins ⁽²⁾	OFF, 0 to 80	High	3	3	3	mmHg
Heart rate ⁽²⁾	OFF, 0 to 240	High	150	120	100	bpm
		Low	OFF	OFF	OFF	

Alarm	Adjustment	Limit	Default Values ⁽¹⁾			Unit
			Neonatal	Pediatric	Adult	
SpO ₂ ⁽²⁾	OFF, 0 to 100	Low	85	85	85	%
Apnea	OFF, 0 to 60	High	15	15	15	s
Automatic ⁽³⁾	OFF, 10, 20 and 30		OFF			%

(1) Whenever the equipment is started, or the patient type is changed, or the battery power is drained without the ventilator being connected to the power supply, the alarms will reset to default values.

(2) Alarms available only with optional external sensors.

(3) The alarm limits will be adjusted according to the monitored values. Valid only for basic ventilation alarms (maximum pressure, PEEP, volume, minute volume, rate and FiO₂).



WARNING

- **Alarms will return to default values whenever the device restarts or the patient changes.**
- **The apnea time can be turned off and, in this condition, there will be no backup ventilation.**
- **THE OPERATOR SHALL BE AWARE OF THE RISKS OF KEEPING THE APNEIA ALARM OFF.**
- **Automatic alarm adjustment is based on the monitored values, so it can only be used when the ventilator is NOT in the STAND BY mode and preferably when the parameters are stable.**

10.18.2 Ventilator alarm messages

In the event of one or more ventilator-related alarms, the following messages may be displayed, depending on their respective priorities:

Table 30 - High priority alarms

High priority alarm	Delay	Description
INOPERATIVE EQUIPMENT	< 1 second	Indicates that there was a technical failure in the equipment that needs to be replaced.
LOW BATTERY	< 1 second	When the internal battery charged is at its end. Appropriate means of ventilatory support of the patient should be provided.
APNEA	< 1 second	Means that the time elapsed since the last inspiration is greater than the alarm value set as the maximum apnea time.
LOW O₂ SUPPLY PRESSURE	< 1 second	The oxygen network pressure is below the specified range. This alarm will not be triggered if the O ₂ % parameter is at 21% (air) and the airline is operating within the required specifications.

High priority alarm	Delay	Description
LOW AIR SUPPLY PRESSURE	< 1 second	The air network pressure is below the specified range. This alarm will not be triggered if the O ₂ % parameter is at 100% and the Oxygen network is operating within the required specifications.
COMMUNICATION FAILURE	< 1 second	Indicates that there was a technical failure in the equipment that needs to be replaced.
OBSTRUCTION	< 2 cycles	There is some obstruction in the respiratory circuit that prevents the patient's complete or adequate expiration.
DISCONNECTED	< 5 cycles	There was disconnection of the respiratory circuit or flow sensor lines (if any), which prevents adequate ventilation of the patient.
HIGH PRESSURE	< 2 cycles	The pressure reached exceeded the alarm value set as the upper pressure limit.
LOW PRESSURE	< 2 cycles	The pressure did not reach the alarm value set as the lower pressure limit.
CHECK INTERNAL SENSOR	< 1 second	Indicates that there was a technical failure in the equipment that needs to be replaced.
CHECK BATTERY	< 1 second	Indicates that there may be an issue with the battery.
HIGH VOLUME	< 3 cycles	The patient volume exceeded the alarm value set as its upper limit.
LOW VOLUME	< 3 cycles	The patient volume is below the alarm value set as its lower limit.
FLOWAIR: VERY HIGH TEMPERATURE ¹	< 1 second	Indicates that the temperature of the FlowAir system is too high.
FLOWAIR: FAILURE ¹	< 1 second	Indicates a failure in the FlowAir system
FIO₂ BELOW 18%	< 3 cycles	The inspired fraction of O ₂ is less than 18%.

Table 31 - Medium Priority Alarms

Medium priority alarm	Delay	Description
LIMITED PRESSURE	< 1 second	When the monitored pressure reaches the set maximum pressure. The volume delivered by the ventilator module does not reach the adjusted volume due to the pressure limitation.
CHECK FLOW SENSOR	< 3 cycles	indicates that the flow sensor is disconnected. In these conditions all the monitoring that depends on this sensor (VT, MV, Rate, Vins, Tinsp, I: E, T exp, C stat, C dyn, Res, τ, iT, Volume Leak, VxTime Chart) will NOT be displayed. In volume-controlled ventilation modes, delivered volumes of the equipment will vary by up to ± 10%.
HIGH PEEP	< 3 cycles	The positive end expiratory pressure (PEEP) exceeded the alarm value set as its upper limit.
LOW PEEP	< 3 cycles	The positive end expiratory pressure (PEEP) did not reach the alarm value set as its lower limit.
HIGH MINUTE VOL	< 3 cycles	The minute volume delivered to the patient exceeded the alarm value set as its upper limit.
LOW MINUTE VOL	< 3 cycles	The minute volume delivered to the patient is below the alarm value set as its lower limit.

¹ Only for Fleximag Max 300 and Fleximag Max 700 models

Medium priority alarm	Delay	Description
HIGH RATE	< 3 cycles	The respiratory rate of the patient exceeded the alarm value set as their upper limit.
LOW RATE	< 3 cycles	The respiratory rate of the patient did not reach the alarm value set as their lower limit.
HIGH LEAKAGE	< 2 cycles	The measured leakage flow has exceeded the maximum compensation limit.
HIGH O ₂ SUPPLY PRESSURE	< 1 second	Indicates that the network pressure is above the specified value.
HIGH AIR SUPPLY PRESSURE	< 1 second	Indicates that the network pressure is above the specified value.
HIGH DRIVING PRESSURE	< 2 cycles	The driving pressure exceeded the alarm value set as the upper limit of the driving pressure.
LOW DRIVING PRESSURE	< 2 cycles	The expired driving pressure is below the alarm value set as the lower limit of the driving pressure.
HIGH FIO ₂	< 3 cycles	The inspired fraction of O ₂ has exceeded the alarm value set as its upper limit.
LOW FIO ₂	< 3 cycles	The inspired fraction of O ₂ did not reach the alarm value set as its lower limit.
LOW BATTERY	< 1 second	The battery is less than half of its total capacity.
INVERTED I:E RATIO	< 1 second	It indicates that the ratio I: E is inverse, i.e., the time of the inspiratory phase is greater than the time of the expiratory phase.
ASSURED VOLUME NOT REACHED	< 1 second	The adjusted volume could not be reached.
FLOWAIR: HIGH TEMPERATURE ¹	< 1 second	Indicates that the temperature of the FlowAir system is increasing.

Table 32 - Low priority alarms

Low priority alarm	Delay	Description
APNEA WITH BACKUP	< 1 second	Indicates that apnea occurred and triggered a backup ventilation because there was one available or enabled.
NO AC POWER	< 1 second	Indicates that the ventilator is disconnected from the electrical network.
LOW AIR SUPPLY PRESSURE	< 1 second	The air network pressure is below the specified range and the flowair system was activated.



WARNING

- **When receiving alarm information, promptly provide service to resolve the problem.**
- **Once the situation that requires the total silence of the audible alarm is terminated, the alarms shall be reactivated for patient safety.**

¹ Only for Fleximag Max 300 and Fleximag Max 700 models



ATTENTION

- *To silence the audible alarm, press the MUTE button. Audible alarms will be disabled for the set period or until a new alarm occurs.*
- *It may be dangerous if different alarm pre-configurations are used for the same or similar equipment in the same area, such as an intensive care unit or cardiac surgery room.*
- *The equipment will always start with the audio volume set to the maximum level (6) regardless of the level set when it was turned off.*
- *If an alarm occurs and the audio volume is set to a value below the maximum level (6), the volume will increase gradually every 15 seconds until it reaches its maximum limit.*

10.18.3 Ventilator Alert Messages

In the event of one or more ventilator related alerts, the following messages may be displayed:

Table 33 - Alert messages

Message	Delay	Description
ASSIST TRIGGER: FLOW	< 1 second	Indicates the occurrence of an assisted trigger, generated by increased inspiratory flow.
ASSIST TRIGGER: PRESSURE	< 1 second	Indicates the occurrence of an assisted trigger, generated by a pressure drop.
MANUAL TRIGGER	< 1 second	Indicates the occurrence of an assisted trigger, generated manually by the operator.
SPONT TRIGGER: FLOW	< 1 second	Indicates the occurrence of a spontaneous trigger, caused by increased inspiratory flow.
SPONT TRIGGER: PRESSURE	< 1 second	Indicates the occurrence of a spontaneous trigger, generated by a pressure drop.
SPONT TRIGGER: MANUAL	< 1 second	Indicates the occurrence of a spontaneous trigger, generated manually by the operator.

Message	Delay	Description
TO VENTILATE, HOLD STANDBY	< 1 second	Indicates that the unit is in standby and the button must be used to resume ventilation.
CO₂: CONNECTED	< 1 second	CO ₂ monitoring is active.
VENTILATION LOCKED	< 1 second	Indicates that the settings are locked to prevent accidental changes. If any adjustments need to be made the lock button must be pressed.
AUTOTEST LEAK FAILED	< 1 second	Indicates that the leak test has not been performed or the leak is above the tolerated limit.
ACTIVE SIGH	< 1 second	Sigh feature is enabled.

10.18.4 IRMA CO₂ sensor alarm messages

In the event of one or more alarms related to the IRMA CO₂ sensor, the following messages may be displayed, according to their respective priorities:

Table 34 - High priority alarms

High priority alarm	Delay	Description
HIGH EtCO₂	< 3 seconds	The expired CO ₂ rate exceeded the alarm value set as the upper limit of EtCO ₂ .
LOW EtCO₂	< 3 seconds	The expired CO ₂ rate is below the alarm value set as the lower limit of EtCO ₂ .
HIGH ICO₂	< 3 seconds	The inspired CO ₂ rate exceeded the alarm value set as the upper limit of CO _{2i} .
CO₂: REBOOT IRMA	< 3 seconds	Indicates that the CO ₂ sensor must be disconnected and reconnected.
CO₂: HARDWARE ERROR	< 3 seconds	Indicates that the CO ₂ sensor must be replaced.
CO₂: REPLACE ADAPTOR	< 3 seconds	Indicates that the airway adaptor should be replaced.
CO₂: NO ADAPTOR	< 3 seconds	Indicates that the airway adaptor must be properly connected.
CO₂: CALIBRATION LOST	< 3 seconds	Indicates that the sensor has lost its original factory calibration.
CO₂: SPEED ERROR	< 3 seconds	Indicates that the sensor motor is out of range.

Table 35 - Medium Priority Alarms

Medium Priority Alarm	Delay	Description
CO ₂ : ZERO REQUIRED	< 3 seconds	Indicates the need to reset the IRMA CO ₂ sensor.
CO ₂ : OUT OF RANGE	< 3 seconds	Indicates that the CO ₂ reading is incorrect or inaccurate.
CO ₂ : PARAMETER OUT OF RANGE	< 3 seconds	Indicates that some parameter is outside the specified range and prevents the correct reading of CO ₂ .
CO ₂ : SPEED ERROR	< 3 seconds	Problem with the CO ₂ motor sent by Phase In.
CO ₂ : CALIBRATION LOST	< 3 seconds	Problem with factory calibration of CO ₂ sent by Phase In.
CO ₂ : NO ADAPTOR	< 3 seconds	Indicates that the airway adaptor must be properly connected.
CO ₂ : CHANGE ADAPTER	< 3 seconds	Indicates that the IRMA CO ₂ sensor must be replaced.
CO ₂ : REPLACE IRMA	< 3 seconds	Indicates that there is a problem with Hardware and should be replaced.
CO ₂ : RESTART IRMA	< 3 seconds	Indicates that the IRMA CO ₂ sensor must be disconnected and reconnected.

Table 36 - Low priority alarms

Low priority alarm	Delay	Description
CO ₂ : DISCONNECTED	< 3 seconds	CO ₂ monitoring was interrupted during use.



WARNING

- *When receiving alarm information, promptly provide service to resolve the problem.*
- *Once the situation that requires the mute of the audible alarm is terminated, the alarms shall be reactivated for patient safety.*



ATTENTION

- *It may be dangerous if different alarm pre-configurations are used for the same or similar equipment in the same area, such as an intensive care unit or cardiac surgery room.*

10.18.5 Oximeter Alarm Messages

In the event of one or more alarms related to the oximeter, the following messages may be displayed, according to their respective priorities:

Table 37 - High priority alarms

High priority alarm	Delay	Description
HIGH PULSE	< 3 seconds	The patient's heart rate exceeded the alarm value set as the upper limit.
LOW PULSE	< 3 seconds	The patient's heart rate is below the alarm value set as the lower limit.
LOW SpO₂	< 3 seconds	The saturation rate of O ₂ is below the alarm value set as the lower limit of SpO ₂ .

Table 38 - Medium Priority Alarms

Medium priority alarm	Delay	Description
SPO₂: CABLE	< 3 seconds	Problem with the SpO ₂ cable sent by Masimo.
SPO₂: FALIURE	< 3 seconds	Problem with the SpO ₂ sensor sent by Masimo.
SPO₂: NO PATIENT	< 3 seconds	Indicates that the sensor is not on the patient's finger.
SPO₂: NO SENSOR	< 3 seconds	Indicates that the sensor is not connected.
SPO₂: LOW PERFUSION	< 3 seconds	Indicates signal quality too low.
SPO₂: NO PULSE	< 3 seconds	Indicates that it is not detecting the beat signal.
SPO₂: UNKNOWN SENSOR	< 3 seconds	Indicates that the connected sensor is not responding as expected.
SPO₂: LIGHT NOISE	< 3 seconds	Indicates that light interference is disrupting the reading.
SPO₂: EXCEDED OF LIGHT	< 3 seconds	Indicates that it is not possible to make the measurement due to excessive light.
SPO₂: NO ADHESIVE SENSOR	< 3 seconds	Indicates that the adhesive sensor is not connected.
SPO₂: LOW SIGNAL QUALITY	< 3 seconds	Indicates that signal quality is not good enough for monitoring.
SPO₂: DETECTING PULSE	< 3 seconds	Indicates that SpO ₂ is initiating calculations.
SPO₂: BOARD FAILURE	< 3 seconds	Indicates that the SpO ₂ board has stopped working.
SPO₂: CONNECTED	< 3 seconds	Indicates that the sensor is connected.

Medium priority alarm	Delay	Description
SPO ₂ : DISCONNECTED	< 3 seconds	Indicates that SpO ₂ has been disconnected during monitoring.

Table 39 - Low priority alarms

Low priority alarm	Delay	Description
SPO ₂ : DEMO	< 3 seconds	Indicates that SpO ₂ is in demo mode.



WARNING

- *When receiving alarm information, promptly provide service to resolve the problem.*
- *Once the situation that requires the mute of the audible alarm is terminated, the alarms shall be reactivated for patient safety.*



ATTENTION

- *It may be dangerous if different alarm pre-configurations are used for the same or similar equipment in the same area, such as an intensive care unit or cardiac surgery room.*

10.18.6 Alarm Testing

To perform the alarm tests, a breathing circuit and a simulator balloon are necessary.



ATTENTION

- *Never perform the alarm test with the patient connected to the equipment.*

10.18.6.1 FiO₂ alarm

To test the high O₂ concentration alarm, set the alarm's upper limit concentration below the concentrating adjusted in the mode parameters. To test the low O₂ concentration alarm, set the alarm's lower limit concentration above the concentration adjusted in the mode parameters.

10.18.6.2 Airway pressure alarm

To test the high pressure alarm, put the ventilator in PCV mode, set PEEP to zero, Control Pr to 5 and set the Ppeak alarm's upper limit to 5. Assemble the respiratory circuit, put the equipment to cycle and press the balloon so that the monitored pressure is greater than the set pressure. To test the low pressure alarm, set the Ppeak alarm's lower limit above inspiratory pressure monitored in the equipment.

10.18.6.3 Volume alarm

To test the high volume alarm, set the alarm's upper limit volume below the volume adjusted in the mode parameters. To test the low volume alarm, set the alarm's lower limit volume above the volume adjusted in the mode parameters.

10.18.6.4 No AC power alarm

Remove the power supply.

10.18.6.5 Low battery alarm

Leave the equipment on, unplugged and cycling for 3 hours.

10.18.6.6 Gas supply failure alarm

Disconnect the air and O₂ supply.

10.18.6.7 PEEP alarm

To test the high PEEP alarm, set the alarm's upper limit PEEP below the adjusted in the mode parameters. To test the low PEEP alarm, set the alarm's lower limit PEEP above the PEEP adjusted in the mode parameters.

10.18.6.8 Disconnection alarm

Disconnect the simulator from the respiratory circuit and leave the equipment cycling.

10.18.6.9 Obstruction alarm

The obstruction alarm is triggered when some form of obstruction occurs in the respiratory circuit that prevents complete expiration of the patient. In PEDIATRIC and ADULT mode, the criteria for triggering this alarm is based on the ratio of average values obtained from PEEP and the limiting pressure (P_{max}). When the pressure value is above the average of the reference parameters (PEEP and P_{max}), the alarm is triggered.

In NEONATAL mode, the obstruction alarm is activated when the airway pressure is above the SET PRESSURE + 5cmH₂O. When an occlusion occurs in the respiratory circuit, the ventilator drives an overpressure valve system that relieves pressure in the circuit in order to preserve the integrity of the patient's lungs.

To test the obstruction alarm operation, place the ventilator in normal operation and press the expiratory branches of the patient circuit to simulate the obstruction and observe the activation of the alarm.

10.18.6.10 Minute volume alarm

To test the high minute volume alarm, set the alarm's upper limit minute volume below the monitored in the equipment. To test the low minute volume alarm, set the alarm's lower limit minute volume above the minute volume monitored in the equipment.

10.18.6.11 Respiratory rate alarm

To test the high rate alarm, set the alarm's upper limit rate below the monitored in the equipment. To test the low rate alarm, set the alarm's lower limit rate above the rate monitored in the equipment.

10.18.6.12 Apnea alarm

To test the apnea alarm, set the apnea time to 3 seconds and adjust a low rate in the ventilatory mode. Put the ventilator in normal operation and wait the set time to occur the alarm.

10.18.6.13 EtCO₂ alarm

To test the EtCO₂ alarm, connect a capnograph to the ventilator. Assemble the respiratory circuit to use with the capnograph and blow the airway adapter so that the ventilator monitors the capnograph parameters. To test the high EtCO₂ alarm, set the alarm's upper limit EtCO₂ below the monitored in the equipment. To test the low EtCO₂ alarm, set the alarm's lower limit EtCO₂ above the EtCO₂ monitored in the equipment.



WARNING

- *Airway adapters are single use. The re-use may result in cross contamination.*
-

10.18.6.14 iCO₂ alarm

To test the iCO₂ alarm, connect a capnograph to the ventilator. Assemble the respiratory circuit to use with the capnograph and blow the airway adapter so that the ventilator monitors the capnograph parameters. To test the high iCO₂ alarm, set the alarm's limit iCO₂ below the monitored in the equipment.



WARNING

- *Airway adapters are single use. The re-use may result in cross contamination.*
-

10.18.6.15 Heart rate alarm

To test the heart rate alarm, connect an oximeter to the ventilator. Assemble the respiratory circuit and connect the oximeter in the finger to monitor the oximeter parameters. To test the high heart rate alarm, set the alarm's upper limit heart rate below the monitored in the equipment. To test the low heart rate alarm, set the alarm's lower limit heart rate above the heart rate monitored in the equipment.

10.18.6.16 SpO₂ alarm

To test the SpO₂ alarm, connect an oximeter to the ventilator. Assemble the respiratory circuit and connect the oximeter in the finger to monitor the oximeter parameters. To test the SpO₂ alarm, set the alarm's limit SpO₂ below the monitored in the equipment.

10.18.7 Battery test

To check the battery status, enter the Battery menu, where the voltage and battery status information will be displayed. To test the operation of the battery, disconnect the equipment from the electrical power and check that the remaining battery power is displayed in percentage. Wait a few minutes, reconnect the machine to the electrical power and check that the battery is recharged.

10.19 Performance Specifications

Table 40 - Performance specifications

Parameter	Specification		Unit	Tolerance
Valve Response Time $T_{0.90}$	10		ms	± 20%
Maximum Flow in Supporting Pressure and Spontaneous Breathing	Adult	180	L/min	± 10%
	Neonatal	60	L/min	± 10%

Observation

- Pressure-controlled ventilation is recommended for leakage flows greater than the limit specified above.
- In this case the maximum compensated flow may be greater than 100 L / min.

10.20 Respiratory Circuit Specifications

Table 41 - Expiratory and inspiratory branch resistance specifications

Respiratory circuit	Flow (L/min)	Resistance (hPa or cmH ₂ O) ¹			
		Circuit	Circuit + Flow Sensor	Circuit + Flow Sensor + HME Filter	Circuit + Flow Sensor + CO ₂ Sensor + HME Filter
Neonatal	2.5	0.3	0.85 ⁽¹⁾		
Pediatric	15.0	0.2	1.7	1.8	1.9 ⁽¹⁾
Adult	30.0	0.4	0.7	1.55	1.75 ⁽¹⁾

¹ Maximum resistance to assure the accuracy.

Table 42 - Respiratory circuit compliance specification

Respiratory circuit	Pressure (cmH ₂ O)	Compliance default ⁽¹⁾ (mL/cmH ₂ O)	Maximum Compliance ⁽²⁾ (mL/cmH ₂ O)
Neonatal	60 ± 3	0,5	1.5
Pediatric	60 ± 3	1	4
Adult	60 ± 3	2	5

¹ Default compliance will be used if the auto test is not performed or if a failure occurs.

² Maximum compliance to assure the accuracy

Table 43 - Operation Specifications

Description	Specification	Tolerance
Recommended Maximum Operating Temperature	37 °C	± 3 °C
Recommended maximum operating pressure	120 cmH2O	± 2 cmH2O

10.21 Maintenance and Calibration Specifications

Table 44 - Maintenance and calibration specifications

Description	Specification	Tolerance
Inspection and replacement of integrated valve diaphragm	Under inspection or 5,000 hours or 12 months (whichever comes first)	± 500 h / ± 1 month
Inspection and replacement of the O ₂ galvanic cell	Recommended replacement if there are problems with calibration or 10,000 hours or 24 months (whichever comes first)	
inspection and Replacement of Internal Battery	10,000 hours or 24 months (whichever comes first)	
Equipment Inspection	5,000 hours or 12 months (whichever comes first)	
Calibration of equipment	5,000 hours or 12 months (whichever comes first)	

10.22 IRMA CO₂ Sensor Specifications

Table 45 - IRMA CO₂ sensor - General specifications

General specifications	
Description	Mainstream monitoring sensor with infrared technology.
Dimensions (L x P X A)	38 x 37 x 34mm (1.49" x 1.45" x 1.34")
Cable length	2.50m (± 0.02m)
Weight	<25g (without cable) <38g (with cable)
Atmospheric Operating Pressure	525 to 1200cmH2O (525cmH2O corresponds to an altitude of 4572m or 15,000 feet).
Storage and Transport Atmospheric Pressure	500 to 1200cmH ₂ O.
Mechanical resistance	Supports 1m repeated drops on a hard surface.
Power Supply	4.5 to 5.5 VDC, Max 1.0W (power measured with 5V and less than 350mA for 200ms).
Surface Temperature (Ambient Temperature of 23 ° C)	Max: 41°C / 106°F.
Airway Adapter	<u>Adult / Pediatric (Disposable):</u> Adds less than 6ml of dead space;

General specifications	
	Loss of pressure less than 0.3cmH2O at 30L / min.
	<u>Neonatal (Disposable):</u>
	Adds less than 1ml of dead space;
	Loss of pressure less than 1.3cmH2O at 10L / min.

Table 46 - IRMA CO₂ Sensor - Outputs

Outputs	
Breath Detection	Adaptive threshold, minimum 1 vol% change in CO ₂ concentration.
Respiratory rate	0–150 ±1 bpm. The respiration rate is displayed after three breaths and the average value is updated every breath.
Fi and ET	Fi and ET are displayed after one breath and have a continually updated breath average.
Wave forms	CO ₂ .
Diagnostic Parameters	Atmospheric pressure, software and hardware revision, serial number.
Information	New Breath Detection, Apnea, Check Adapter, Unspecified Accuracy, and Sensor Error.
Method to calculate the gas level reading	The highest concentration of CO ₂ during one breathing cycle with a weight function applied to favor values closer to the end of the cycle.

Table 47 - IRMA CO₂ sensor - Gas analyzer

CO ₂ Gas Analyzer	
Sensor	Gas Analyzer with 2 to 9 NDIR Channels (Non-Dispersive Infrared) that measures in the range of 4 to 10µm. Makes pressure, temperature and interference correction throughout the spectral range.
Data acquisition rate	10 kHz
Data sampling rate	20 Hz
Calibration	Zeroing recommended when changing Airway adapter. No span calibration required for the IR bench.
Warm-up Time	Information on the concentration is analyzed and sent every 10 seconds. Total accuracy in measurements: 1 minute.
Rise time (at 10 L/min)	CO ₂ ≤ 90ms.
Total System Response Time	< 1s.

Observations

- Accurate CO₂ monitoring is achieved 1 minute after initialization.

Table 48 - IRMA CO₂ sensor – Accuracy I

Exactitude / accuracy of measurements (under standard conditions)		
Type of Gas	Range (AX+)	Accuracy
CO ₂	0 to 15	±(0,2 vol% + 2% of reading)
	15 to 25	Not specified
<i>Note: Concentration of gases expressed in units of percentage volume.</i>		

Table 49 - IRMA CO₂ sensor - Accuracy II

Exactitude / accuracy of measurements (under all conditions)	
Type of Gas	Accuracy
CO ₂	± (0.3 vol% + 4% of reading)
<i>Note: The accuracy specification is valid for any specified environmental condition, except as described in the table below with "Effects of Gas and Steam Interference".</i>	

Table 50 - IRMA CO₂ sensor - Effects of gas and vapor interference

Effects of gas and vapor interference		
Gases or Steam	Gas Level	CO ₂
N ₂ O	60 vol%	--- (1 and 2)
HAL	4 vol%	--- (1)
ENF, ISO, SEV	5 vol%	+8% of the read measurement. ⁽³⁾
DES	15 vol%	+12% of the read measurement. ⁽³⁾
Xe (Xenon)	80 vol%	-10% of the read measurement. ⁽³⁾
He (Helium)	50 vol%	-6% of the read measurement. ⁽³⁾
Measured dose inhaler propellant	Not designed for use with a measured-dose propellant inhaler.	
C ₂ H ₅ OH (Ethanol)	0,3 vol%	--- (1)
C ₃ H ₇ OH (Isopropanol)	0,5 vol%	--- (1)
CH ₃ COCH ₃ (Acetone)	1 vol%	--- (1)
CH ₄ (Methane)	3 vol%	--- (1)
CO (Carbon Monoxide)	1 vol%	--- (1)
NO (Nitrogen Monoxide)	0,02 vol%	--- (1)
O ₂	100 vol%	--- (1 and 2)

(1) Despicable interference. Effects of interference do not change the values in the "Accuracy of Measurements (under all conditions)" table above.

(2) For sensors that are not measuring N₂O and / or O₂, concentrations must be entered manually by the user.

- (3) Interference at the indicated gas level. For example, 50 vol% Helium typically decreases the values read in CO₂ by 6%. This means that if the mixture contains 5.0 vol% CO₂ and 50 vol% Helium, the CO₂ concentration measurement will normally be calculated as follows:

$$(1 - 0.06) * 5.0 \text{ vol\%} = 4.7 \text{ vol\% CO}_2.$$

According to ISO 80601-2-55.

Table 51 – Quantitative effects of humidity and condensation

Temp [C]	RH [%]	P [hPa]	H ₂ O part.pres [hPa]	Err _{rel} [%]	Err _{rel} ATPD [%]	Err _{rel} [%] BTPS
10	20	1013	2	0	-0,2	+6,0
20	20	1013	5	0	-0,5	+5,7
25	0	1013	0 (ATPD)	0	0	+6,2
25	23	1013	7,3	0	-0,7	+5,5
25	50	1013	16	0	-1,6	+4,6
30	80	1013	42	0	-4,1	+2,0
37	100	1013	63 (BTPS)	0	-6,2	0
37	100	700	63	0	-9,0	-2,8

10.23 Oximeter specifications

Table 52 – Performance specifications

Specification Criteria	Functional SpO ₂ (%)	Pulse Rate (bpm)	Perfusion index (%)	PVI (%)
Display Range	0.0 - 100.0 %	25 – 240 bpm	0.02 – 20.0 %	0 – 100 %
Calibration Range	70 – 100 %	25 – 240 bpm	0.10 – 20.0 %	-
Calibration Standard	Invasive Co-oximeter	ECG and Patient Simulator	Patient Simulator	-
No motion accuracy (ms)	≤ 2.0 %	≤ 3.0 bpm	-	-
Motion Accuracy (rms)	≤ 3.0 %	≤ 5.0 bpm	-	-
Resolution	≤ 0.1 %	≤ 1 bpm	≤ 0.01 %	≤ 1 %
Time to Display	≤ 8, ≤ 12 s	≤ 8, ≤ 12 s	≤ 8, ≤ 12 s	-
Asystole Detection Time	≤ 8 s	≤ 8 s	≤ 8 s	-
Delay	≤ 10 s	≤ 10 s	≤ 10 s	-
Response Time	≤ 20 s	≤ 20 s	≤ 20 s	-
Display Update Frequency	≥ 1 Hz	≥ 1 Hz	≥ 1 Hz	≥ 1 Hz
Averaging Time (s)	2-4, 4-6, 8, 10, 12, 14, 16	-	-	-

Table 53 – Environmental specifications

Characteristic	Specification
Operating Conditions	
Incandescent light intensity	100 k Lux (Sunlight)
Fluorescent Light Intensity	10 k Lux
Fluorescent Light Frequency	50, 60 Hz ± 1.0 Hz
Temperature	5 to 40 °C
Humidity	15 to 95 %, non-condensing
Pressure	500 to 1060 mbar

Storage Conditions	
Temperature	-40 to 70 °C
Humidity	15 to 95 %, non-condensing

10.24 Electromagnetic Compatibility

Changes or modifications to this equipment not expressly approved by MAGNAMED may cause EMC problems with this equipment. Contact MAGNAMED for technical assistance.

This equipment has been designed and tested to comply with the applicable EMC standards as described below:

- Immunity: IEC 60601-1-2
- Emission: CISPR11 (Group 1 – Class A)
- Approvals: OS / IEC 60601-1

This equipment has been designed and tested to meet the following essential requirements: deliver volume within alarm limits or generate alarm condition; FiO₂ monitoring; generate PEEP alarm conditions when PEEP is above or below the alarm limit; monitor expired volume; alarm condition to indicate when the electrical voltage has fallen below what is necessary to maintain normal operation; alarm condition to indicate when the battery charge is near the end; technical alarm condition to indicate failure of the air and oxygen network; limiting the reverse flow between the gas port to values below 100 mL/h; ensure the accuracy of the oxygen level within the specified range.



WARNING

- ***Use of mobile phones or other radio-frequency (RF) emitting equipment near the system may cause unexpected or adverse results. Monitor operation if there are sources of radiofrequency emission in the vicinity.***
- ***Use of other electrical equipment in or near the system may cause interference. Before use on the patient, you should check that the equipment functions normally in the set configuration.***
- ***The use of this adjacent equipment or other equipment shall be avoided as it may result in improper operation. If this is necessary, this and other equipment shall be observed to verify that they are operating normally.***
- ***The use of accessories, transducers, and cables other than those specified or provided by Magnamed may result in high electromagnetic emissions or reduced***

electromagnetic immunity from this equipment and result in improper operation.

- *Portable RF communication equipment (including peripherals such as antenna cables and external antennas) shall not be used within 30 cm of any part of the Fleximag Max, including cables specified by Magnamed. Failure to do so may cause performance degradation.*
- *If essential performance is requested or degraded due to electromagnetic disturbances, the ventilator may stop ventilating. In this case, the operator shall provide manual ventilation.*

10.24.1 Guidelines and manufacturer's declaration - Electromagnetic emissions

The Fleximag Max is intended for use in the electromagnetic environment specified below. It is recommended that the customer or user of Fleximag Max ensure that it is used in such an environment.

Emissions Testing	Compliance	Electromagnetic environment - guidelines
RF emissions CISPR 11	Group 1	Fleximag Max uses RF energy only for its internal functions. Therefore, RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The Fleximag Max is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings for domestic purposes, as long as the following warning is respected:
Harmonic emissions IEC 61000-3-2	Class A	Warning: This equipment / system is intended for use by healthcare professionals only. This equipment / system may cause radio interference or may impair the operation of nearby equipment. It may be necessary to take measures such as reorienting or relocating the Fleximag Max or shielding the site.
EF voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

10.24.2 Guidelines and manufacturer's declaration - Electromagnetic immunity

The Fleximag Max is intended for use in the electromagnetic environment specified below. It is recommended that the purchaser or user of Fleximag Max ensure that it is used in such an environment.

Immunity Test	IEC-60601 Test Level	Level of Compliance	Electromagnetic environment - guidelines
IEC 61000-4-2 - Electrostatic discharge (ESD)	± 8 kV per contact ± 15 kV through the air	± 8 kV per contact ± 15 kV through the air	Floors should be made of wood, concrete or ceramic. If the floors are covered with synthetic material, the relative humidity should be at least 30%.

Immunity Test	IEC-60601 Test Level	Level of Compliance	Electromagnetic environment - guidelines
IEC 61000-4-4 - Electrical Fast Transit / Pulse Train ("Burst")	± 2 kV at the input power input c.a.	± 2 kV at the input power input c.a.	Mains power quality should be that of a typical commercial or hospital environment.
IEC 61000-4-5 - Surge	± 1 kV line (s) to line (s) ± 2 kV line (s) to ground	± 1 kV line (s) to line (s) ± 2 kV line (s) to ground	Mains power quality should be that of a typical commercial or hospital environment.
IEC 61000-4-11 - Voltage dips	0% UT; 0.5 cycle at 0 °, 45 °, 90 °, 135 °, 180 °, 225 °, 270 ° and 315 ° 0% UT; 1 cycle (single phase: at 0 °) 70% UT; 25/30 cycles (single phase: at 0 °) <5% UT (> 95% UT voltage drop) for 5 seconds.	0% UT; 0.5 cycle at 0 °, 45 °, 90 °, 135 °, 180 °, 225 °, 270 ° and 315 ° 0% UT; 1 cycle (single phase: at 0 °) 70% UT; 25/30 cycles (single phase: at 0 °) <5% UT (> 95% UT voltage drop) for 5 seconds.	Quality of energy supply should be that of a typical hospital or commercial environment.
IEC 61000-4-11 - Voltage interruptions	0% UT; 250/300 cycles	0% UT; 250/300 cycles	Quality of energy supply should be that of a typical hospital or commercial environment.
IEC 61000-4-8 - Magnetic field of the frequency of mains (50/60 Hz)	30 A/m	30 A/m	Magnetic fields in the frequency of feeding should be at levels characteristic of a typical location in a typical hospital or commercial setting

Note: UT is the a.c. mains voltage prior to application of the test level.

The **Fleximag Max** is intended for use in the electromagnetic environment specified below. It is recommended that the purchaser or user of **Fleximag Max** ensure that it is used in such an environment.

Immunity Test	IEC-60601 Test Level	Level of Compliance	Electromagnetic environment - guidelines
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz outside the ISM bands (a)	3 V	Portable and mobile RF communications equipment should be used no closer to any part of Fleximag Max, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2\sqrt{P}$
	10 Vrms 150 kHz to 80 MHz in the ISM bands (a)	10 V	
radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	10 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz

where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).^b
Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^c, should be less than the compliance level in each frequency range^d
Interference may occur in the vicinity of equipment marked with the symbol:



NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people..

^a The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.

^b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.

^c Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Fleximag Max is used exceeds the applicable RF compliance level above, the Fleximag Max should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Fleximag Max.

^d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V / m.

Recommended separation distances between portable and mobile RF communications equipment and the Fleximag Max

The **Fleximag Max** is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the **Fleximag Max** can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **Fleximag Max** as recommended below, according to the maximum output power of the communications equipment.

Separation distance according to frequency of transmitter (m)

Rated maximum output power of transmitter (W)	150 kHz to 80 MHz outside ISM bands	150 kHz to 80 MHz in ISM bands	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.12	0.23
0.1	0.38	0.38	0.38	0.73
1	1.2	1.2	1.2	2.3
10	3.8	3.8	3.8	7.3
100	12	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.

NOTE 3 An additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.

NOTE 4 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Fields in the vicinity of RF wireless communication equipment

The Fleximag Max enclosure interface was tested as specified in the table below using the test methods specified in IEC 61000-4-3

Band [MHz]	Freq. test [MHz]	Modulation	Test level [V/m]
380 to 390	385	Pulse, 18 Hz	27
430 to 470	450	FM, 1 kHz, Deviation from \pm 5kHz	28
704 to 787	710 745 780	Pulse, 217 Hz	9
800 to 960	810 870 930	Pulse, 18 Hz	28
1.700 to 1.990	1.720 1.845 1.970	Pulse, 217 Hz	28
2.400 to 2.570	2.450	Pulse, 217 Hz	28
5.100 to 5.800	5.240 5.500 5.785	Pulse, 217 Hz	9

Observation

- ***The EMISSION characteristics of this equipment make it suitable for use in industrial areas and hospitals (ABNT NBR IEC / CISPR 11 class A). If used in a residential environment (for which ABNT NBR IEC / CISPR 11 class B is normally required), this equipment may not provide adequate protection for radio frequency communication services. The user may need to take mitigation measures, such as relocating or reorienting the equipment.***

10.24.3 Electrical safety

The following are the precautions that shall be observed when combining these items (non-medical equipment) with the system.



WARNING

- *Items that do not comply the requirements of IEC 60601-1 cannot be placed within 1.5 m of the patient.*
- *Do not connect non-medical electrical equipment directly to the wall outlet. Use AC power supply with own transformer. Otherwise, current leakage will increase above the levels accepted by IEC 60601-1 under normal conditions and single fault conditions. This may cause dangerous electrical shock to the patient or the operator.*
- *After connecting any equipment to these outlets, subject the system to a complete leakage current test (in accordance with IEC 60601-1).*
- *The electromedical system operator shall not touch non-medical electrical equipment and the patient simultaneously. This may cause dangerous electrical shock to the patient or the operator.*
- *Do not connect a multiple socket or extension cable to the electromedical system.*

10.25 Audible acoustic energy

Table 54 – Audible acoustic energy specification

Configuration	Sound Pressure Level	Sound Power Level
Volume \geq 300 mL	47 dB \pm 2 dB	55 dB \pm 2 dB
300 mL \geq Volume \geq 50 mL	46 dB \pm 2 dB	54 dB \pm 2 dB
Volume \leq 50 mL	45 dB \pm 2 dB	53 dB \pm 2 dB

10.26 Mask to non-invasive ventilation

Table 55 – Mask to NIV ventilation specification

Specification	
Adult/ pediatric connection	22 mm
Neonatal connection	15 mm



WARNING

- Use only MASKS specified by MAGNAMED with local registration and CE mark.
- Use MASK suitable for patient type to avoid excessive leakage.

10.27 Respiratory Circuit

Table 56 – Respiratory circuit specification

Specification	
Adult/ pediatric connection	22 mm
Neonatal connection	15 mm
Resistance	≤ 0.3 mbar/L.s-1



WARNING

- Use only respiratory circuit specified by MAGNAMED with local registration and CE mark.
- Use respiratory circuit suitable for patient type to avoid excessive leakage.

10.28 HME Filter



ATTENTION

- It is recommended the use of HME filter following the specifications given below with local registration and CE mark.

Table 57 – HME filter specification

Specification	
Compliance	ISO 23328-1, ISO 23328-2, ISO 9360-1 and ISO 9360-2

Connection to conical connectors	22 mm (female)	
Connection to breathing tubes	22 mm (male) / 15 mm (female)	
Bacterial and viral filtration efficiency	99,99%	
Resistance	30 L/min	2,02 cmH ₂ O
	60 L/min	5,19 cmH ₂ O
	90 L/min	9,37 cmH ₂ O

10.29 HEPA Filter



ATTENTION

- *In order to avoid cross contamination, a HEPA filter following the specifications given below, with local registration and CE mark shall be used.*

Table 58 – HEPA filter specification

Specification		
Compliance	ISO 23328-1 and ISO 23328-2	
Connection to conical connectors	22 mm (female)	
Connection to breathing tubes	22 mm (male) / 15 mm (female)	
Bacterial and viral filtration efficiency	99,99%	
Resistance	30 L/min	1 cmH ₂ O

11 IRMA CO₂ Sensor (optional)

11.1 Intended Use

The IRMA™ CO₂ mainstream sensor has been developed to monitor respiratory gases in adult, pediatric and neonatal patients, in therapy units, surgical centers and emergency rooms.

It consists of a single-track sensor with up to 9 channels of non-dispersive infrared ("NDIR") for gas identification, a barometric pressure sensor, an electric voltage regulator and a microprocessor. The unit weighs less than 25g.

Carbon dioxide (CO₂) concentrations are monitored along with other parameters such as respiratory rate, gas waveform, and concentration of each gas during inspiration and expiration.

The airway adapter fits snugly into the CO₂ sensor. This equipment uses XTP™ window technology. The airway adapter shall be positioned between the endotracheal tube and the breathing circuit allowing XTP windows positioned on the sides of the sensor to measure gas concentrations.

Running on a standard low-voltage continuous electric voltage, the sensor is designed to meet portability requirements and low power consumption, typically below 1 Watt. It is designed to be extremely easy to integrate with any monitoring device, allowing the visualization of gas information in real time.

11.2 Instructions for use

The IRMA CO₂ sensor is designed to be used connected to Magnamed ventilators and any other compatible monitoring device. Its function is to monitor the signal and the concentration value of the gas in real time.

It shall not be used as the sole means of monitoring the patient. It shall always be used in conjunction with other vital signs monitoring equipment and this process should be accompanied by a specialist capable of analyzing the patient's conditions.

The IRMA CO₂ sensor is designed for use by trained and authorized healthcare professionals only.



WARNING

- *The sensor shall not have direct contact with the patient during its use.*
 - *Perform the leak test in the patient circuit with the airway adapter connected to the circuit.*
 - *Always check the gas reading and waveform in the ventilator before to connect the airway adapter in the patient circuit.*
 - *Ventilation with cyclic pressure up to 100 cmH₂O does not interfere with gas reading.*
-

11.2.1 Reset procedure

To ensure high accuracy in the values measured by the IRMA sensor, the following zeroing recommendations shall be followed.



WARNING

- *Incorrect sensor reset will result in incorrect reading of measured values.*
-



Observations

- *The sensor zeroing option is available in the CALIBRATION window as soon as the sensor is identified and ready for use.*
 - *It may take a few seconds until the sensor is ready for the zeroing process.*
-

Zeroing shall be done by connecting an airway adapter to the sensor without connecting it to the respiratory circuit. When the gas monitoring signals are stable, press the button to start the zero reset.

Exceptional care should be taken to avoid breathing near the sensor before or during the reset. The presence of ambient air (21% O₂ and 0% CO₂) in the airway adapter is of crucial importance for successful zeroing. If the calibration error message appears immediately after the end of the zeroing procedure, it shall be repeated.

Zeroing shall be performed every time the airway adapter is replaced. It shall also be executed whenever there is a displacement of the baseline (offset) in any of the gas measurements or when any of the

alarm messages are displayed: "CO₂: PARAM OUT OF RANGE ", "CO₂: OUT OF RANGE" or "CO₂: ZERO REQUIRED".

After connecting the sensor or replacing the airway adapter, wait at least one minute before starting the zeroing procedure to allow the sensor to warm up.

The green LED on the sensor will blink for 5 seconds while the reset process is in progress.

11.2.2 LED status information

Table 59 - IRMA CO₂ Status LEDs

Color (state)	Significance
Green (constantly on)	System OK
Green (blinking)	Zeroing in progress
Blue (constantly on)	Anesthetic agent present
Red (constantly on)	Sensor error
Red (blinking)	Check the adapter

11.3 Preventive maintenance

The calibration of the gas shall be checked at regular intervals by reference instrument.

11.4 Important Notices



WARNING

- *The IRMA CO₂ sensor shall be operated exclusively by trained and authorized medical personnel.*
- *The sensor shall not be used with flammable anesthetic agents.*
- *Airway adapters shall not be reused. Reusing a disposable adapter may cause cross-infection.*
- *Do not use the adult / pediatric airway adapter in neonates as the adapter adds a 6mL dead space to the patient's breathing circuit.*
- *Do not use the neonate airway adapter in adult patients as this adapter may add excessive resistance.*
- *Measurements may be affected by radio-frequency communication equipment or cellular devices.*
- *The user must make sure that the sensor is used in environments that comply with the electromagnetic environment specifications stated in this manual.*
- *Do not use the airway adapter with measured-dose inhalers or nebulized medications as they may affect the transmission of light into the sensor windows.*

-
- *The IRMA CO₂ sensor is designed to be an adjunct device in patient monitoring, so its information should be analyzed along with other measurements and symptoms.*
 - *Incorrect zeroing may result in erroneous measurements.*
 - *Replace the airway adapter if there is condensation inside the adapter.*
 - *Use only airway adapters and sensors produced by Masimo.*
 - *The sensor shall not come in direct contact with the patient during use.*
 - *Do not connect the airway adapter between the endotracheal tube and the respiratory circuit elbow as this can cause patient secretions to block the adapter windows, causing incorrect sensor operation.*



Figure 2 - Incorrect and correct positioning of the airway adapter



ATTENTION

- *Do not apply electrical voltage to the sensor cable.*
 - *Do not use the sensor in environments with specifications outside the limits stated in your technical specification.*
-

12 Pulse Oximeter (optional)

12.1 Intended Use

The Masimo MS-2040 pulse oximeter is a self-contained solution that allows the safe measurement of SpO₂, pulse rate, perfusion index and PVI, even with movement or low perfusion.

12.2 Principle of Operation

The MS board of the Masimo SET pulse oximeter is based on three principles:

1. Differential absorption of oxyhemoglobin and deoxyhemoglobin from red and infrared lights (spectrophotometry).
2. The volume of arterial blood in the tissue and the light absorbed in the blood changes (plethysmography).
3. Arteriovenous shunting is highly variable and its fluctuation of absorbance through the venous blood is a major component of noise during the pulse.

The MS board of the Masimo SET pulse oximeter as well as the traditional pulse dosimetry determines SpO₂ by the passage of red and infrared light in a capillary bed and changes the measurement during the pulsatile cycle. Red and infrared (LED) emitting diodes on the oximetry sensors serve as light source, the photodiode serves as a photodetector.

Traditionally pulse oximetry assumes that all pulsations in the light absorption signal are caused by oscillations in arterial blood volume. Assuming that the blood flow in the sensor region passes entirely through the capillary bed rather than some arteriovenous shunt. Traditional pulse oximetry calculates the ratio of pulsatile absorption (AC) to mean absorption (DC) at each of the two wavelengths, 660nm and 905nm:

$$S(660) = AC(660) / DC(660)$$

$$S(905) = AC(905) / DC(905)$$

oximeter then calculates the ratio between these two pulse signals arterial absorption:

$$R = S(660) / S(905)$$

This value of R is used to find the saturation (SpO₂) in a verification table made by the oximeter software. The values of this table were obtained based on studies on human blood, performed with healthy adult volunteers, in situations of induced hypoxia.

The MS board of the Masimo SET pulse oximeter assumes that the arteriovenous shunt is highly variable in the floating absorption because the venous blood is a component of noise during the pulse. The

MS board decomposes S (660) and S (905) into an arterial signal, plus a noise component and calculates the arterial signal ratio without noise:

$$S(660) = S1+N1$$

$$S(905) = S2+N2$$

$$R = S1 / S2$$


Again, R is the ratio of two pulsed arterial absorption signals and its value is used to find SpO₂ saturation in an empirical derived equation in the oximeter software. The values in the empirical derived equation were obtained based on studies on human blood, performed with healthy adult volunteers, in situations of induced hypoxia.

The above equations are combined and a reference noise (N') is determined:

$$N' = S(660) - S(905) \times R$$

If there is no noise N'= 0: then S (660) = S (905) x R which is the same ratio as the traditional pulse oximeter.

The equation for reference noise is based on the value of R, at the value sought to determine SpO₂. The MS board software scans all possible R values that correspond to the SpO₂ values between 1% and 100% and generates a value N 'for each of these R values. Signals S (660) and S (905) are (ACC) that produces a power output versus a possible SpO₂ value as shown in the following figure where R corresponds to SpO₂ = 97%:



The DST curve has two peaks: a peak corresponding to the highest saturation is selected as the SpO₂ value. The entire sequence is repeated every two seconds in the most recent four seconds of the received data. Concluding the SpO₂ of the MS board, corresponds to the running evaluation of arterial saturation of the hemoglobin updated every two seconds.

12.3 Important Notices



WARNING

- *Risk of explosion. Do not use the MS pulse oximeter in the presence of flammable anesthetics or other flammable substances in contamination such as air, oxygen enriched environments or nitrous oxide.*
- *The pulse oximeter shall not be used as an apnea sensor.*
- *The heart rate is based on the optical detection of the peripheral pulse flow and in this way may not detect certain arrhythmias. Therefore, the pulse oximeter should not be used as a replacement of ECG-based arrhythmia analysis*
- *A pulse oximeter can be considered a warning device. As an indicator of the patient's deoxygenation tendency, blood samples can be analyzed by co-oximetry laboratory to complete understanding of the patient's condition.*
- *The MS board of the pulse oximeter must be operated only by a qualified person.*
- *The manual, instructions for use and all precautions and specifications shall be read before use.*
- *There is danger of electric shock. Do not remove the monitor cover except for battery replacement.*
- *The operator may perform the maintenance procedures described in the product manual.*
- *Contact Magnamed service for oximeter repairs.*
- *As with all medical equipment, position the cable on the patient in order to reduce the possibility of entanglement or strangulation.*
- *Interfering substances:*
 - *Carboxyhemoglobin may erroneously increase readings. The level of increase is approximately equal to the amount of carboxyhemoglobin present.*
 - *Dyes or any substance containing dye, which change usual pigmentation of the artery may cause reading errors.*
- *Do not use the sensor during magnetic resonance (MRI) scanning:*
 - *Induced current can potentially cause burns.*
 - *The oximeter may affect the resonance image and the resonance unit may affect the oximetry measurements.*
- *If the fidelity of any measure does not seem reasonable, first check the patient's vital signs for alternative measures and make sure the oximeter is working properly.*
- *Before using, read the instructions for using the LNOP / LNCS sensors carefully.*
- *Use only Masimo oximetry sensors to measure SpO₂.*
- *Damage to the tissue may be caused by improper application or misuse of LNOP / LNCS sensors.*
- *Inspect the sensor location as directed in the instructions for use of the product to*

ensure skin integrity and correct positioning and adhesion of the sensor.

- *Do not use damaged LNOP / LNCS sensors.*
- *Do not use LNOP / LNCS sensors with exposed optical components.*
- *Do not immerse the sensor in water, solvents, or cleaning solution (sensors and connectors are not waterproof).*
- *Do not sterilize by irradiation, steam or sterilization by oxides.*
- *See the cleaning instructions in the instructions for use for reusable Masimo LNOP / LNCS sensors.*
- *Do not use damaged cables.*
- *Do not immerse the patient cable in water, solvents or cleaning solution (patient cables are not waterproof).*
- *Do not sterilize by irradiation, steam or sterilization by oxides.*
- *See the cleaning instructions in the instructions for use for reusable Masimo LNOP / LNCS patient cables.*
- *Do not use the adult / pediatric sensor in neonate patients. It may cause incorrect measurements of physiological parameters.*
- *Do not use the neonate sensor in adult / pediatric patients. It may cause incorrect measurements of physiological parameters.*

13 Technical Service



WARNING

- *FlexiMag Max is a life-support device and therefore, if any repairs or maintenance on this equipment is needed, only ask for an authorized Magnamed service.*
 - *Failure to perform the preventive maintenance may affect the safety and performance of the ventilator.*
 - *DO NOT OPERATE the equipment if it is not operating according to the specifications contained in this operating manual.*
 - *Before sending the equipment to the service technician, observe the cleaning and disinfection process RIGOROUSLY.*
 - *Magnamed will make available upon request the circuit diagrams, component list, descriptions, calibration instructions and other information required to perform an authorized service.*
-

14 Discard

The Fleximag Max ventilator shall be discarded as electrical and electronic equipment. Accessories and consumables must be disposed of as described in the instruction manual. Follow local government recommendations for proper disposal.


















WARNING










- *Dispose of the removable parts of the equipment in accordance with your institution's parts and pieces arrangement protocol.*
- *Follow local government recommendations for environmental protection, especially in cases of disposal of electronic waste or electronic parts.*
- *All parts of Magnamed ventilators that come in contact with patient fluids (e.g. breathing circuit) and that are potentially contaminated, are called semi-critical, and must endure a high-level disinfection or sterilization process before disposal (at the end of their useful lives) or before it is sent in for maintenance.*
- *When disposing of parts of Magnamed ventilators that come in contact with fluids from a patient, indicate as potentially infectious hospital waste.*
- *Disposal of batteries shall be in accordance with local regulations.*
- *Disposal of the galvanic cells shall follow local regulations.*
- *Airway adapters shall be disposed of in accordance with local regulations for medical disposal.*
- *Do not disassemble the equipment. All service or maintenance on the ventilator can only be carried out by qualified technician, trained and duly authorized by MAGNAMED.*

15 Symbology

15.1 Symbols used in the equipment









Table 60 - Symbols used in the equipment









Symbol	Description
	Preventive maintenance period
	Patient
	Equipment with applied part type B
	Equipment with applied part type BF
IP31	Degree of protection IP31 for protection against solid objects 2.5 mm in diameter or more and protected against falling drops
	EC Conformity: Indicates that the equipment complies with the European Council directive 93/42 / EEC for medical equipment
	Authorized Representative in the European Community
	Attention! Consult accompanying documents
	See operating instructions
	Separate collection of electrical / electronic equipment (do not dispose as ordinary waste)
	Manufacture Date
	Manufacturer's identification
	Gas inlet (AR / O ₂)
	On/Off
	Power connection
	Serial connection

Symbol	Description
	USB connection
	Network connection
	Alternating current
	Direct current
	DC Input - External Power Supply
	Fuse
	Identify or advise cleaning or filter change
	US federal law restricts this device to sale by or on the order of a physician.
	Consult documentation that came with the product.

15.2 Symbols used in packaging and labeling

Table 61 - Symbols used in packaging and labeling

Symbol	Description
	Fragile
	Direction of upper box face
	Keep protected from sunlight
	Keep protected from moisture
	Maximum Stacking Quantity
	Temperature limits
	EC Conformity: Indicates that the equipment complies with the European Council directive 93/42 / EEC for medical equipment
	INMETRO

Symbol	Description
	Single use part or accessory. Discard after use. Reprocessing is prohibited.
	Part or accessory not supplied sterile. Sterilize before use.
	Expiration date
	Do not use if packaging is damaged
	Manufacturer's Catalog Number
	Manufacturer serial number
	Manufacturer batch code
	Consult documentation that came with the product.

16 Abbreviations and Terms Used

Table 62 - Abbreviations and terms used

Abbreviation	Meaning
Δ PS	Support pressure delta (pressure above PEEP)
O ₂	O ₂ Concentration
PEEP	Positive end-expiratory pressure
Control Pr	Controlled pressure (pressure above PEEP)
Inferior Pr	Inferior pressure in APRV / DualPAP
Insp Pr	Inspiratory pressure (absolute pressure in neonatal mode)
Limit Pr	Limit pressure
Superior Pr	Superior Pressure in APRV / DualPAP
FI Trigger	Flow trigger (for triggering)
Pr Trigger	Pressure sensitivity (for triggering)
T Inferior	Inferior Time in APRV / DualPAP
T Rise	Rise time (rise ramp)
T Superior	Superior Time in APRV / DualPAP
Ins Time	Inspiratory time
Minute Vol	Minute volume
Vol/weight	Volume by patient weight
NIV or VNI	Noninvasive Ventilation
O ₂ +	Concentration of 50 to 100% O ₂ for a certain time
O ₂ 100%	100% O ₂ concentration for a certain time
MANUAL CYCLE or INSP MANUAL	Manual Cycle trigger
INSP HOLD	Inspiratory hold
EXP HOLD	Expiratory hold
Leakage	Percentage or leakage flow
Cdyn	Dynamic Compliance
Cstat	Static Compliance
E	Elastance
R Rate or f	Respiratory rate
Rate sp or fspont	Spontaneous respiratory rate
I:E	I: E Ratio
Pmean	Pressure Mean
Ppeak	Pressure Peak
Pplat	Plateau Pressure
RE	Expiratory resistance
RI	Inspiratory resistance
RSBi	Rapid Shallow Breathing Index (IRRS, Tobin Index)
TC	Expiratory Time Constant
Te	Expiratory time
Ti	Inspiratory time
Ti/Ttot	Reason for total inspiratory time
Vte or VTE	Total expired volume
Vte sp or VTE spont	Spontaneous total expired volume
Vti or VTI	Total Inspired Volume
MV or VM	Minute volume
MV sp or VM spont	Spontaneous minute volume
WOBi	Imposed Work of Breathing



WARNING

- *The controlled pressure setting (Pr Control) in pediatric or adult patients, refers to a relative pressure, that is, it adjusts the pressure value ABOVE PEEP.*
 - *The resulting inspiratory pressure will be the sum of the pressure controlled and PEEP.*
-

17 Biocompatibility Statement

We declare under our sole responsibility that all materials used in parts applied to Fleximag Max, such as silicone and polysulfide, have been widely used in the medical field over time, with no effects related to toxicity or tissue effects, ensuring their biocompatibility.



WARNING

- *Common accessories purchased from third parties SHALL be registered with local government.*

18 Warranty

The products manufactured and marketed by **MAGNAMED TECNOLOGIA MÉDICA S/A** are warranted against in material and manufacturer defects throughout Brazil, as described below.

The warranty period of the equipment is 12 months. For batteries and accessories, the period is 3 months, as long as their original characteristics are maintained, these terms begin from the date of acquisition by the first purchaser of the product, as indicated in the purchase receipt issued by MAGNAMED TECNOLOGIA MÉDICA S / A.

Warranty liability is limited to exchange, repair and labor for any parts that are defective or that do not meet the specifications contained in the Product Operation Manual.

The warranty is limited to the product which is used under normal conditions and for the intended purpose and whose preventive maintenance and replacement of parts and repairs are carried out according to the instructions in the Product Operation Manual by personnel authorized by the manufacturer.

The warranty does not cover defects caused by improper use or installation, accident, improper sterilization, service, installation, operation or alteration performed by personnel not authorized by the manufacturer.

Breakage or absence of seals or warranty seals by unauthorized personnel results in loss of product warranty.

Parts subject to wear or deterioration by normal use, adverse conditions of use, misuse or accidents are not covered under warranty.

Any expenses and risks with the transportation of the product are not covered by the guarantee.




For equipment sold with extended warranty, it will only be valid if the preventive maintenance suggested by Magnamed is carried out in accordance with chapter 9 of this manual.

There is no express or implied warranty other than as set forth above.

19 Training

To request training, contact Magnamed Product Specialist who will direct you to the nearest authorized representative. The training is performed in-loco and is expected to be executed in 1 hour.

Annex I – Differences between models

Specification	Fleximag Max 700	Fleximag Max 500	Fleximag Max 300
BACK VIEW			
FLOWAIR	YES	NO	YES
AIR PROPORTIONAL VALVE	YES	YES	NO
Alarm - FLOWAIR: Very High Temperature	YES	NO	YES
Alarm - FLOWAIR: HIGH TEMPERATURE	YES	NO	YES
Alarm - FLOWAIR: Failure	YES	NO	YES
Alarm - LOW AIR SUPPLY PRESSURE	YES	YES	NO
Alarm - HIGH AIR SUPPLY PRESSURE	YES	YES	NO
LOW PRESSURE O ₂ INLET	YES	NO	YES


This product is only intended for use in pulmonary ventilation and shall be operated only by qualified professionals.

MAGNAMED

Manufacturer
Technical assistance
Customer Service



Rua Santa Mônica, 801, 831
06715-865 – Parque Industrial San José - Cotia – SP Brasil
Tel/Fax: +55 11 4615-8500
E-mail: magnamed@magnamed.com.br
Website: www.magnamed.com.br
CNPJ: 01.298.443/0002-54
State registration: 149.579.528.111



European Representative
CMC Medical Devices & Drugs S.L.
Horacio Lengo Nº 18, CP 29.006
Málaga – Spain
Tel.: +34 951214054
Fax: +34 952330100
E-mail: info@cmcmmedicaldevices.com
Website: www.cmcmmedicaldevices.com

Technical Responsible: Toru Miyagi Kinjo
Inscription CREA-SP 5061555031
Legal Responsible: Wataru Ueda